GE Healthcare

IG Option Installation Manual

OPERATING DOCUMENTATION

5161904

Rev 4

© 2005-2009, General Electric Company.
All Rights Reserved.

This page intentionally left blank

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

Rev. June 13, 2006

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS

CAUTION
Risk of
Explosion

Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

Revision History

Revision	Date	Reason for change
1	11/20/05	Initial release.
2	9/30/06	Add Jarrell IG Installation.
3	6/18/07	Updated Jarrell IG Option Kit Materiel List
4	2/03/09	Added IG installation flowchart (Chapter 1) and procedures on All-In-One Console (Chapter 3). Added 2 IG option installation procedure (Chapter 4)

Table of Contents

Chapter 1	
Introduction.....	15
Section 1.0	
Flowchart	16
1.1 Flowchart (For LCGOC/ 2IG Option for GOC4).....	16
1.2 Flowchart (For All-In-One Console).....	17
Section 2.0	
Material List	18
Chapter 2	
Installation Procedure (For LCGOC).....	19
Section 1.0 Overview	19
Section 2.0 Procedure.....	19
2.1 Power-Off (Shut Down) the Console	19
2.2 Install the IG Node.....	19
2.3 Reconfiguration.....	21
2.4 Verify DARC Node is Up.....	21
2.5 Verify IG Node	22
2.6 Start Applications.....	23
2.7 Retro Recon Test.....	23
Section 3.0 Finalization.....	23
Chapter 3	
Installation Procedure (For All-In-One Console)	25
Section 1.0 Overview	25
Section 2.0 Procedure.....	25
2.1 Power-Off (Shut Down) the Console	25
2.2 Install the IG Node.....	25
2.3 Reconfiguration.....	27
2.4 Verify IG Node	27
2.5 Start Applications.....	28
2.6 Retro Recon Test.....	28
Section 3.0 Finalization.....	28
Chapter 4	
Installation Procedure (For 2 IG Option for GOC4)	29
Section 1.0 Overview	29
Section 2.0 Procedure.....	29
2.1 Power-Off (Shut Down) the Console	29

2.2	Install the IG Node	29
2.3	Reconfiguration	32
2.4	Verify DARC Node is Up	33
2.5	Verify IG Node.....	33
2.6	Start Applications	34
2.7	Retro Recon Test	34
	Section 3.0 Finalization	34
Chapter 5		
Attaching Rating Plate		35

Chapter 1

Introduction

By installing this Option, Image Generator (IG) becomes available for LCGOC and All-In-One Console. By installing 2 IG Option, GOC4 Console will be upgrade from 1IG to 3 IGs. The image reconstruction speed will be increased. This manual describes the installation procedure.

IG Option: BrightSpeed and BrightSpeed Select with LCGOC and AIO Console.


2 IG Option: LightSpeed RT16 and BrightSpeed with GOC4 Console (Consist of two IG option kit).

Section 1.0

Flowchart


1.1 Flowchart (For LCGOC/ 2IG Option for GOC4)

Below Illustration shows main flow of the installation procedures of the option on LCGOC and GOC4.

1.2 Flowchart (For All-In-One Console)

Below illustration shows main flow of the installation procedures of the option on All-In-One console.

Section 2.0 Material List

Verify that the contents of this option kit are correct according to the Material List below.

IG Option Kit (5144667) - Westville

ITEM	PART NO.	DESCRIPTION	QTY.
1	2362872	IMAGE GENERATION COMPUTER - GRE (IG)	1
2	2374965-2	CABLE - ETHERNET CAT5E X 610MM LG	1
3	5161904	IG Option Installation Manual (This Manual)	1
4	5146271	Rating Plate of IG	1

IG Option Kit (5192881) - Jarrell Darc/All-In-One Console/GOC4

ITEM	PART NO.	DESCRIPTION	QTY.
1	5147443-2	IMAGE GENERATION COMPUTER - GRE (IG)	1
2	5127222-2	Shield Ethernet Cable	1
3	5161904	IG Option Installation Manual (This Manual)	1

Chapter 2

Installation Procedure (For LCGOC)

Section 1.0 Overview

This procedure describes the installation of the IG Node. After the IG Node has been installed, several tests are required:

- Communication from the DARC Node to the IG Node
- Testing of the RAC (Reconstruction Acceleration Card)

Section 2.0 Procedure

NOTICE Prevent permanent damage to the static-sensitive boards. Attach an anti-static wrist-strap to your wrist and to a bare metal grounding point.

2.1 Power-Off (Shut Down) the Console

- 1.) Select one of the following methods to Power OFF the Operator Console:
- 2.) If Applications are up, click on the [Shut Down] button and select [Shutdown].
- 3.) If Applications are down, open a Unix Shell at the Toolchest. Enter (type) the following commands, at the indicated prompts:

```
ctuser@hostname: su -  
Password: #bigguy  
root@hostname: halt
```
- 4.) Power OFF the console at the front panel switch when 'System Halted' message displays.

2.2 Install the IG Node

Note: Jarrell DARC and IG Nodes require software release 06MW29.x.

- 1.) Remove the front cover and rear cover of the console.
- 2.) Remove and retain the two (2) front mounting screws from the console chassis.
- 3.) Slide the IG Node into the Operator Console.
- 4.) Replace the two (2) front mounting screws on the front of the IG Node.
- 5.) Verify the power switch (located at the rear of the IG Node, if present) is in the ON position. (See [Figure 2-1, Rear of "Westville" IG](#))
- 6.) Connect the power cable to the rear of IG. The Power cable is already in the chassis.
- 7.) Connect the Ethernet cable between DARC Ethernet (IG1 (eth1)) and IG Ethernet (eth0). (See [Figure 2-1](#) to [Figure 2-4](#))

Figure 2-1 Rear of "Westville" IG

Figure 2-2 Rear of "Jarrell" IG

Figure 2-3 Rear of "Westville" DARC2


Figure 2-4 Rear of "Jarrell" DARC2

2.3 Reconfiguration

- 1.) Power on the console.
- 2.) When a pink window displays: "You have 5 seconds to cancel the application startup."
Click [Cancel].
- 3.) Open a Shell.
Type: `su -`
Password: `#bigguy`
Type: `reconfig`
- 4.) Click [Config].
- 5.) Select [Hardware]
- 6.) Select [1] for IG number.

- 7.) Click [Accept].
- 8.) Click [OK] to reboot.
- 9.) After approximately 4 minutes a pink pop up window appears with the message:
`CT Software Auto-Start Disabled.`
Click [OK].

2.4 Verify DARC Node is Up

The DARC Node must be up to check the IG Node. Perform `ping` and `rsh` commands on the DARC Node, as follows:

- 1.) Open a Unix Shell.
- 2.) `ctuser@hostname: ping darc`
Verify ping is successful.
Press `<Ctrl+C>` to stop ping process.
- 3.) `ctuser@hostname: rsh darc`
Verify `ctuser@darc` prompt is present.
- 4.) Exit from the DARC2 Node Remote Shell:
`ctuser@darc exit`

Note: This Unix Shell will be used again, so do not exit or close, at this time.

2.5 Verify IG Node

Note: In the case of an IG Node replacement, or troubleshooting with an IG Node swap, the affected IG node must be manually power-cycled (powered down and then up again), at the front panel. The manual power cycle is required for the IG Node to obtain an IP address and an IG#.

- 1.) Power OFF the IG Node at the front panel power switch.

Figure 2-5 Front Panel of "Westville" IG Node


Figure 2-6 Front Panel of "Jarrell" IG Node

- 2.) Power ON the IG Node at the (top) front panel switch.
- 3.) Wait 2 to 4 minutes.
- 4.) In the open Unix shell (window), verify the IG Node:
 - a.) `ctuser@darc: rsh ig1`
 - b.) If `rsh` to `ig1` is successful, then type: `exit`, to get out of `ig` remote shell
 - c.) Type: `exit`, to get out of DARC Node remote shell.
 - d.) Type: `exit`, to close the Unix Shell.

2.6 Start Applications

- 1.) In the open Unix Shell, start Application Software:
`ctuser@hostname: st`
- 2.) A pink pop-up Attention Message appears:
`OC initializing. Please wait...`
- 3.) Select [OK] for other pink message boxes that appear throughout the reboot process.
- 4.) If any Recon Selftest Failures are encountered, review the error log.

2.7 Retro Recon Test

Reconstruct images to verify the IG functions well.

- 1.) Click on the [RETRO RECON] selection on the left Monitor.
- 2.) Select [RAT GOLD] series.
- 3.) Click [SELECT SERIES] .
- 4.) Click on [CONFIRM] .
- 5.) Verify the image(s) reconstruct and are displayed on the right Monitor.
- 6.) Click on [QUIT] .

Section 3.0 Finalization

- 1.) Take some scans (e.g., scout, axial & helical).
- 2.) Verify the image(s) reconstruct and are displayed.
- 3.) Install the front and rear Operator Console covers.

Chapter 3

Installation Procedure (For All-In-One Console)

Section 1.0 Overview

This procedure describes the installation of the IG Node. After the IG Node has been installed, several tests are required:

- Communication from the Host Computer to the IG Node
- Testing of the RAC (Reconstruction Acceleration Card)

Section 2.0 Procedure

NOTICE Prevent permanent damage to the static-sensitive boards. Attach an anti-static wrist-strap to your wrist and to a bare metal grounding point.

2.1 Power-Off (Shut Down) the Console

- 1.) Select one of the following methods to Power OFF the Operator Console:
- 2.) If Applications are up, click on the [Shut Down] button and select [Shutdown].
- 3.) If Applications are down, open a Unix Shell at the Toolchest. Enter (type) the following commands, at the indicated prompts:

```
ctuser@hostname: su -  
Password: #bigguy  
root@hostname: halt
```
- 4.) Power OFF the console at the front panel switch when 'System Halted' message displays.

2.2 Install the IG Node

Note: AIO and IG Nodes require software release 08BW20.x or later.

- 1.) Remove the front cover and rear cover of the console.
- 2.) Remove and retain the two (2) front mounting screws from the console chassis.
- 3.) Slide the IG Node into the Operator Console.
- 4.) Replace the two (2) front mounting screws on the front of the IG Node.
- 5.) Connect the power cable to the rear of IG. The Power cable is already in the chassis.
- 6.) Connect the Ethernet cable between AIO Ethernet (IG1 (eth0)) and IG Ethernet (eth0). (See [Figure 3-1](#) and [Figure 3-2](#))

Figure 3-1 Rear of "Jarrell/All-In-One" IG


Figure 3-2 Rear of Host Computer in "AIO" Console

2.3 Reconfiguration

- 1.) Power on the console.
- 2.) When a pink window displays: "You have 5 seconds to cancel the application startup."
Click [Cancel].
- 3.) Open a Shell.
Type: `su -`
Password: `#bigguy`
Type: `reconfig`
- 4.) Click [Config].
- 5.) Select [Hardware]
- 6.) Select [1] for IG number.

- 7.) Click [Accept].
- 8.) Click [OK] to reboot.
- 9.) After approximately 4 minutes a pink pop up window appears with the message:
`CT Software Auto-Start Disabled.`
Click [OK].

2.4 Verify IG Node

Note: In the case of an IG Node replacement, or troubleshooting with an IG Node swap, the affected IG node must be manually power-cycled (powered down and then up again), at the front panel. The manual power cycle is required for the IG Node to obtain an IP address and an IG#.

- 1.) Power OFF the IG Node at the front panel power switch.


Figure 3-3 Front Panel of "Jarrell" IG Node

- 2.) Power ON the IG Node at the (top) front panel switch.
- 3.) Wait 2 to 4 minutes.
- 4.) In the open Unix shell (window), verify the IG Node:
 - e.) `ctuser@hostname: rsh ig1`
 - f.) If `rsh` to `ig1` is successful, then type: `exit`, to get out of `ig` remote shell
 - g.) Type: `exit`, to close the Unix Shell.

2.5 Start Applications

- 1.) In the open Unix Shell, start Application Software:
`ctuser@hostname: st`
- 2.) A pink pop-up Attention Message appears:
`OC initializing. Please wait...`
- 3.) Select [OK] for other pink message boxes that appear throughout the reboot process.
- 4.) If any Recon Selftest Failures are encountered, review the error log.

2.6 Retro Recon Test

Reconstruct images to verify the IG functions well.

- 1.) Click on the [RETRO RECON] selection on the left Monitor.
- 2.) Select [RAT GOLD] series.
- 3.) Click [SELECT SERIES] .
- 4.) Click on [CONFIRM] .
- 5.) Verify the image(s) reconstruct and are displayed on the right Monitor.
- 6.) Click on [QUIT] .

Section 3.0 Finalization

- 1.) Take some scans (e.g., scout, axial & helical).
- 2.) Verify the image(s) reconstruct and are displayed.
- 3.) Install the front and rear Operator Console covers.

Chapter 4

Installation Procedure (For 2 IG Option for GOC4)

Section 1.0 Overview

This procedure describes the installation of the 2 IG Option for GOC4 will use two IG option kits (5192881). After the IG Nodes have been installed, several tests are required:

- Communication from the Host Computer to the IG Node
- Testing of the RAC (Reconstruction Acceleration Card)

Section 2.0 Procedure

NOTICE Prevent permanent damage to the static-sensitive boards. Attach an anti-static wrist-strap to your wrist and to a bare metal grounding point.

2.1 Power-Off (Shut Down) the Console

- 1.) Select one of the following methods to Power OFF the Operator Console:
- 2.) If Applications are up, click on the [Shut Down] button and select [Shutdown].
- 3.) If Applications are down, open a Unix Shell at the Toolchest. Enter (type) the following commands, at the indicated prompts:

```
ctuser@hostname: su -  
Password: #bigguy  
root@hostname: halt
```
- 4.) Power OFF the console at the front panel switch when 'System Halted' message displays.

2.2 Install the IG Node

- 1.) Remove the front cover and rear cover of the console.
- 2.) Remove two faceplates covering the slots for Image Generator (IG) # 2 and # 3. Keep screws to secure the two IGs that shipped in the upgrade kit.
- 3.) Insert IG # 2 into the "IG # 2" slot, sliding it on the rails until flush with the GRE chassis.
- 4.) Secure IG # 2 using two screws from Step 2.
- 5.) Insert IG # 3 into the "IG # 3" slot, sliding it on the rails until flush with the GRE chassis.
- 6.) Secure IG # 3 using two screws from Step 2.

Figure 4-1 Slots for Adding IG Nodes

- 7.) Move to the rear of the console to continue the install of the Ethernet cables and power cords to the IGs.
- 8.) Connect the Ethernet cables for IG # 2 and IG # 3 to the DARC as shown [Figure 4-2](#) to [Figure 4-4](#). (IG # 2 ETH0 <-> DARC IG2, IG # 3 ETH0 <->DARC IG3)

Figure 4-2 Rear of "Jarrell" IG

Figure 4-3 Rear of "Westville" DARC2



Figure 4-4 Rear of "Jarrell" DARC2

- 9.) Connect power cable (pre-wired in console chassis) to the added IG's a.c. input receptacle.

2.3 Reconfiguration

- 1.) Power on the console.
- 2.) When a pink window displays: "You have 5 seconds to cancel the application startup."
Click [Cancel].
- 3.) Open a Shell.
Type: `su -`
Password: `#bigguy`
Type: `reconfig`
- 4.) Click [Config].
- 5.) Select [Hardware]
- 6.) Select [3] for IG number.

- 7.) Click [Accept].
- 8.) Click [OK] to reboot.
- 9.) After approximately 4 minutes a pink pop up window appears with the message:
CT Software Auto-Start Disabled.
Click [OK].

2.4 Verify DARC Node is Up

The DARC Node must be up to check the IG Node. Perform `ping` and `rsh` commands on the DARC Node, as follows:

- 1.) Open a Unix Shell.
- 2.) `ctuser@hostname: ping darc`
Verify ping is successful.
Press `<Ctrl+C>` to stop ping process.
- 3.) `ctuser@hostname: rsh darc`
Verify `ctuser@darc` prompt is present.
- 4.) Exit from the DARC2 Node Remote Shell:
`ctuser@darc exit`

Note: This Unix Shell will be used again, so do not exit or close, at this time.

2.5 Verify IG Node

Note: In the case of an IG Node replacement, or troubleshooting with an IG Node swap, the affected IG node must be manually power-cycled (powered down and then up again), at the front panel. The manual power cycle is required for the IG Node to obtain an IP address and an IG#.

- 1.) Power OFF the IG Node at the front panel power switch.


Figure 4-5 Front Panel of "Jarrell" IG Node

- 2.) Power ON the IG Node at the (top) front panel switch.
- 3.) Wait 2 to 4 minutes for the IG Node to boot up.
- 4.) In the Unix Shell, verify the `cat /usr/g/config/INFO` file for number of IGs (NUMIG) actively configured in the Operator Console.
- 5.) In the Unix Shell, verify the following IG Node(s) as applicable:

Note: For an Ethernet cabling connection issue, use the command `'ethtool -p eth#'` (where # represents the Ethernet address) as root on the DARC2 Node or IG Node to determine if the Ethernet cable is correctly connected from the DARC2 Node to the IG Node. The port will blink or flash when activated and `<Ctrl+C>` will stop the test. This command can be used on the Host, DARC2 Node and the IG Node(s).

Example for IG Node # 1 Ethernet cable:

```
[root@ig1] ethtool -p eth0 (<Ctrl+C> to stop the flashing Ethernet port LED)
```

```
{ctuser@hostname} rsh darc
{ctuser@darc} su -
password: #<password>
{root@darc} rsh ig1
{root@ig1} exit
{root@darc} rsh ig2
{root@ig2} exit
{root@darc} rsh ig3
{root@ig3} exit
{root@darc} exit
{ctuser@darc} exit
```

2.6 Start Applications

- 1.) In the open Unix Shell, start Application Software:
ctuser@hostname: **st**
- 2.) A pink pop-up Attention Message appears:
OC initializing. Please wait...
- 3.) Select [OK] for other pink message boxes that appear throughout the reboot process.
- 4.) If any Recon Selftest Failures are encountered, review the error log.

2.7 Retro Recon Test

Reconstruct images to verify the IG functions well.

- 1.) Click on the [RETRO RECON] selection on the left Monitor.
- 2.) Select [RAT GOLD] series.
- 3.) Click [SELECT SERIES] .
- 4.) Click on [CONFIRM] .
- 5.) Verify the image(s) reconstruct and are displayed on the right Monitor.
- 6.) Click on [QUIT] .

Section 3.0 Finalization


- 1.) Take some scans (e.g., scout, axial & helical).
- 2.) Verify the image(s) reconstruct and are displayed.
- 3.) Install the front and rear Operator Console covers.

Chapter 5

Attaching Rating Plate

Attach the provided rating plate to the following location (in the label area left to the rating plate of the OC).

OC Rear View

www.gehealthcare.com