GE Healthcare

**BrightSpeed Elite**  
**16 FPS Option Installation**  
Service Documentation

5322976-1EN

Rev 2

© 2008, General Electric Company,  
All Rights Reserved.

---

# **LEGAL NOTES**

---

## **TRADEMARKS**

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

## **COPYRIGHTS**

All Material, Copyright © 2008-2009 by General Electric Company, Inc. All rights reserved.

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## **DAMAGE IN TRANSPORTATION**

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

Rev. June 13, 2006

## **CERTIFIED ELECTRICAL CONTRACTOR STATEMENT**

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

## **IMPORTANT...X-RAY PROTECTION**

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

---

## LITHIUM BATTERY CAUTIONARY STATEMENTS

**CAUTION**  
Risk of  
Explosion

Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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### OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

# Revision History

Revision	Date	Reason for change
2	03/23/09	Updated Material List on Chapter 1, Section 2.0
1	10/28/08	Initial release.

# Table of Contents

<b>Chapter 1</b>	
<b>Introduction</b> .....	<b>17</b>
<b>Section 1.0</b>	
<b>Flowchart</b> .....	<b>17</b>
<b>Section 2.0</b>	
<b>Material List</b> .....	<b>19</b>
<b>Chapter 2</b>	
<b>Installation Procedures</b> .....	<b>21</b>
<b>Section 1.0 Remove the original Host Computer</b> .....	<b>21</b>
<b>Section 2.0 Install new Host Computer</b> .....	<b>23</b>
<b>Section 3.0 Install DARC Computer</b> .....	<b>24</b>
<b>Section 4.0 Install IG Computers</b> .....	<b>26</b>
<b>Section 5.0 Console Power on and verification</b> .....	<b>27</b>
5.1 Power on Console .....	27
5.2 DARC Computer verification.....	27
5.3 IG Computers verification .....	28
<b>Section 6.0 Check the hardware setting in reconfig</b> .....	<b>29</b>
<b>Section 7.0 Restore system state</b> .....	<b>30</b>
<b>Section 8.0 Install Software Options</b> .....	<b>31</b>
<b>Section 9.0 Load Other Service Software</b> .....	<b>34</b>
9.1 Load Class C .....	34
9.2 Load Class M.....	34
<b>Section 10.0 Retro Recon Test</b> .....	<b>35</b>
<b>Section 11.0 Verify Recon Data Path Diagnostic</b> .....	<b>35</b>
<b>Section 12.0 Re-save Reconfig Info File on System State</b> .....	<b>35</b>
<b>Section 13.0 Finalization</b> .....	<b>37</b>
<b>Chapter 3</b>	
<b>Final Check</b> .....	<b>39</b>
<b>Section 1.0</b>	
<b>Attaching The Rating Plate</b> .....	<b>39</b>

# Chapter 1

## Introduction

---

This option is only applied to BrightSpeed Elite system With All-In-One console.

By installing this option, DARC2 computer and 3 Image Generator (IG) becomes available for BrightSpeed Elite system. The image reconstruction speed will be increased. This manual describes the installation procedure.

### Section 1.0

#### Flowchart

This procedure describes the procedure in below sequence:

- Remove original Host Computer in the console
- Install new Host Computer, shipped in 16 FPS option package
- Install DARC Computer
- Install 3 IG Computer

After all the computer have been installed, several tests are required:

- Communication from Host Computer to DARC, and DARC to 3 IG Computer
- Testing of retro recon and online scan


**NOTICE** Prevent permanent damage to the static-sensitive boards. Attach an anti-static wrist-strap to your wrist and to a bare metal grounding point.

**NOTICE** The Operating software and application software has been loaded on Host Computer and DARC Computer by manufacturer.

Below illustration shows the main flow of the installation procedures of the option.

## Section 2.0 Material List

Verify that the contents of this option kit (5322976) are correct according to the Material List below.

ITEM	PART NO.	DESCRIPTION	QTY.
1	5183547-38	GEHC MI and CT Integrated Buy Level PC Workstation. HP8400 with 4GB Memory	1
2	5147442	Computer - Data Acquisition Recon Control, 2nd Generation	1
3	5147443-2	Jarrell GRE IG with RAC rev.D or Higher	3
4	5127222-2	Cable, Ethernet, 2ft Lg- CAT5e	4
5	2337414	DVD RAM MEDIA, 9.4 GB	1
6	5322976-1EN	16 FPS Option Installation Manual (This Manual)	1
7	46-303268P1	PRODUCT LOCATOR CARD FORM FOR RATING PLATES	1
8	5328866	BSD Delight and RT16 Xtra Vision Software Collector	1
9	5330693	Service pack 1.0-3 for BrightSpeed Delight	1

---

# **Chapter 2**

## **Installation Procedures**

---

### **Section 1.0 Remove the original Host Computer**

The original Host computer in console should be removed before 16 FPS option installation.

- 1.) If Applications are up, click on the [Shut Down] button.
- 2.) Power OFF the Operator Console front switch, when instructed to “Power Down”.
- 3.) Remove the rear Operator Console cover.
- 4.) Visually verify proper labeling of each cable and disconnect all HP Host Computer rear cables. Refer to below schematic for the cable connections.

- 5.) Remove the front Operator Console cover.
- 6.) Remove and retain the four (4) front mounting screws.

**CAUTION**

**Potential for injury.**

**The HP Computer module weighs approximately 20.5kg (45 lbs).  
Use care, when removing the HP module.**

- 7.) Slide out the HP Host Computer module.

## Section 2.0 Install new Host Computer

The new Host Computer, which is packed in 16 FPS option kit, should be installed in console.

**NOTICE:** If the option have already been loaded on original Host Computer, please perform NIC card exchange to avoid option key re-generation. Detailed procedure please refer to console/ replacement/ service manual.

- 1.) If present, remove the top and bottom covers from the new HP Host computer.  
The computer module will not install properly, with these covers attached.

**CAUTION**

**Potential for injury.**

**The HP Computer module weighs approximately 20.5kg (45 lbs).  
Use care, when removing the HP module.**

- 2.) Slide in the new HP Host Computer module.
- 3.) Replace the four (4) front mounting screws.
- 4.) Reconnect all HP Host Computer rear cables, which is removed in section 1.0.

**NOTICE:** Below three cables should NOT be connected to Host Computer at this time. They will be connected to DARC Computer in Section 4.0

- a.) Fiber cable
- b.) DIP X-ray abort cable (From ICOM)
- c.) RHART cable (From ICOM)

**NOTICE:** All cable connector screw-locks must be finger-tightened, only. Do NOT use any tool to tighten cable screw-locks.

- 5.) At the rear of the HP Host Computer, verify the power switch is in the ON position.

## Section 3.0 Install DARC Computer

- 1.) Remove all the metal plate which are installed at upper position of the Host Computer, and retain the front mounting screws.

**Figure 2-1 Metal plates installed at front of console**

- 2.) Slide the DARC Computer into the Operator Console from the front. The DARC Computer shall be at upper position next to the Host Computer.
- 3.) Install the four (4) front mounting screws on the front of the DARC Computer.
- 4.) Verify the DARC Computer is correct for the Site DAS Type. Refer to [Figure 2-2](#) for the color-coded sticker located on the front panel of the DARC Computer for DIP type information.
- 5.) Connect all rear DARC Computer cables. Refer to [Table 2-1](#) for cables to be connected. (Please see [Figure 2-2](#), for DARC Computer cable connection information.)

Reference Designator	DARC Connection	Connection
B or IG1	IG1 Ethernet port 1	IG 1 Node Ethernet port 0
A or IG2	IG2 Ethernet port 0	IG 2 Node Ethernet port 0
Network Sticker or IG3	IG3 Ethernet port 3	IG 3 Node Ethernet port 0
No sticker or Host	Host Ethernet port 2	Host Computer port C
DIP RX	DIP (verify not connected to TX)	Console Interface Bulkhead Assembly
DIP Scan Abort	J53 DIP Scan Abort	ICOM J53
RHARD	J54 RHARD	ICOM J54
Power	Power	Power Distribution Box

**Table 2-1 DARC Computer cable connection**



Figure 2-2 Jarrell DARC Computer

**NOTICE:** All cable connector screw-locks must be finger-tightened, only. Do NOT use any tool to tighten cable screw-locks.

- 6.) If present, verify the power switch located at the rear of the DARC Computer is in the ON position.

## Section 4.0 Install IG Computers


This section describes the procedure to install 3 IG computers in console.

- 1.) Remove the all the covers at front of console and retain the six (6) front mounting screws from the console chassis.
- 2.) Slide the first IG Computer into the Operator Console.
- 3.) Replace the two (2) front mounting screws on the front of the IG Computer.
- 4.) Verify the power switch (located at the rear of the IG Computer, if present) is in the ON position.
- 5.) Connect the power cable to the rear of IG. The Power cable is already in the chassis.
- 6.) Connect the Ethernet cable between DARC Ethernet (IG1 (eth1)) and IG Ethernet (eth0).
- 7.) Repeat the above procedures from Step 2 to step 6 to install IG 2 and IG 3 computer.  
For IG 2: Connect the Ethernet cable between DARC Ethernet [ IG 2 (eth0) ] and IG 2 Ethernet (eth0).  
For IG 3: Connect the Ethernet cable between DARC Ethernet [ IG 3 (eth3) ] and IG 3 Ethernet (eth0).

Refer to [Figure 2-3](#) and [Figure 2-4](#) for Ethernet cable connection.

**Figure 2-3 Rear of Jarrell IG Computer**

**Figure 2-4 Rear of Jarrell DARC Computer**


- 8.) After complete install 3 IG computers, please install the two metal plates to cover the blank slots at front of console.

## Section 5.0 Console Power on and verification

### 5.1 Power on Console

**NOTICE:** The Operating software and application software has been loaded on Host Computer and DARC Computer by manufacturer.

- 1.) Power ON the Operator Console at the front switch.
- 2.) At the Host Computer verify the front panel Power (PWR) LED is on. If not press the front panel power switch to apply power to the Host computer.
- 3.) At the DARC Computer verify the front panel Power (PWR) LED is on. If not press the front panel power switch (hold 3-10 seconds) to apply power to the DARC computer.
- 4.) At each IG Node, verify the front panel Power (PWR) LED is on. If not press the front panel power switch (hold 3-10 seconds) to apply power to each IG Node. Please refer to [Figure 2-5](#)

**Figure 2-5 Jarrell IG Computer**

- 5.) Stop Applications from auto-starting at the 5 second PINK Box.

### 5.2 DARC Computer verification

- 1.) At the Toolchest, with Application software down, open a Unix Shell.
- 2.) Ping the DARC Computer:  

```
{ctuser@hostname} ping darc
```

Verify ping is successful.  
Press <Ctrl+C> to stop ping process.  
Output similar to the following will appear:  

```
--- darc ping statistics ---  
4 packets transmitted, 4 received, 0% loss, time 3000ms  
rtt min/avg/max/mdev = 0.134/0.168/0.209/0.031 ms
```
- 3.) Remote Shell to the DARC Computer:  

```
{ctuser@hostname} rsh darc
```

Last login: Thu Oct 9 11:17:35 on ttyS1  

```
{ctuser@darc}
```

Verify `ctuser@darc` prompt is present.
- 4.) Verify the DARC Computer DAS Interface Processor Board (DIP) Diagnostics pass.

```
{ctuser@darc} DipDiag -x 1,2,3,4,5,6
```

**NOTICE:** Because tests 7 and 8 require the installation of a loopback cable, they are not being performed.

- 5.) Verify the scan data disks in DARC Computer:  
.  
{ctuser@darc} sudo gre-raid -z -c  
Verify the create partitions diagnostic passes and all output looks good.  
{ctuser@darc} sudo gre-raid -a  
Verify assemble diagnostic passes and all output looks good.  
{ctuser@darc} sudo gre-raid -q  
Verify query diagnostic passes and all output looks good.

### 5.3 IG Computers verification

- 1.) Open a Unix Shell  
{ctuser@hostname} rsh darc
- 2.) Ping IG Node #1:  
{ctuser@darc} ping ig1  
Verify ping is successful.  
Press <Ctrl+C> to stop ping process.
- 3.) Remote Shell to the IG Node #1:  
{ctuser@darc} rsh ig1  
Verify ctuser@ig1 prompt is present.
- 4.) Exit from the IG Node # 1 Remote Shell:  
{ctuser@ig1} exit
- 5.) Ping IG Node #2:  
{ctuser@darc} ping ig2  
Verify ping is successful.  
Press <Ctrl+C> to stop ping process.
- 6.) Remote Shell to the IG Node #2:  
{ctuser@darc} rsh ig2  
Verify ctuser@ig2 prompt is present.
- 7.) Exit from the IG Node # 2 Remote Shell:  
{ctuser@ig2} exit
- 8.) Ping IG Node # 3:  
{ctuser@darc} ping ig3  
Verify ping is successful.  
Press <Ctrl+C> to stop ping process.
- 9.) Remote Shell to the IG Node # 3:  
{ctuser@darc} rsh ig3  
Verify ctuser@ig3 prompt is present.
- 10.) Exit from the IG Node # 3 Remote Shell:  
{ctuser@ig3} exit
- 11.) Exit from the DARC Node Remote Shell:

```
{ctuser@darco} exit  
{ctuser@hostname}
```

12.) The IG Computer(s) RAC (Recon Accelerator Card) will be tested In the following test.

```
{ctuser@hostname} racdiags
```

**NOTICE:** The command `racdiags` will test up to three (3) IG Node RAC's.

13.) Verify rac diagnostic passes for each IG Node configured in the Operator Console.

**NOTICE:** The Unix Shell will be used again, so do not exit or close at this time.

## Section 6.0 Check the hardware setting in reconfig

```
ctuser@hostname: su -  
Password: #bigguy  
{root@hostname} reconfig
```

Click on the "Hardware", check the setting of console and IG.

DAS type: GDAS16

Console type: GOC3 or GOC4

IG # : 3

If no error found, click "accept" to exit reconfig GUI.


**NOTICE:** The console hardware configuration was setup in manufacture.

## Section 7.0 Restore system state

Restore the system state with the system state DVD-RAM shipped with system.

**NOTICE:** Do not restore the reconfig information.

- 1.) Insert the System State DVD-RAM into the SCSI Tower DVD RAM drive.
- 2.) Wait until the DVD drive is ready (i.e., until the front panel green LED is no longer lit).
- 3.) Select: Service icon to access the CSD (Common Service Desktop).
- 4.) Select: [Utilities].
- 5.) Select: [System State]

**Figure 2-6 System State GUI**

- 6.) Select all the tabs except "Reconfig info".
- 7.) Select [Restore]
- 8.) If applicable, select [Yes] when the following message appears: System State Media Status: Please insert a DVD into the drive and press Save again.

Figure 2-7 Save System State Ready

9.) When completed, select [CANCEL]

**NOTICE:** Do NOT select [DISMISS] more than once. Doing so will result in the GUI becoming hung, and the System State will be unusable. Should this occur, Re-Save the System State again.

10.) When completed, select [Dismiss]

11.) Close the *Service Desktop* window at the upper left corner of the screen.

## Section 8.0 Install Software Options

**NOTICE:** For Perfusion and Denta Scan:

If you loaded software in a language other than English, Please reconfig, changing the language to English, to load these two options.

Otherwise, the LFC procedure may failed.

**NOTICE:** Option must be loaded before Protocols.

1.) Insert the Options DVD in the SCSI Tower DVD-RAM drive.


**Note:** *(For GE Healthcare Employees Only)* If you do not have the Option DVD, then you can use the eLicence feature. Follow this link to the eLicence Page.

<http://elicense.gehealthcare.com>

2.) With the Application up, select the [SERVICE DESKTOP].

3.) Select [CONFIGURATION].

4.) Select [INSTALL OPTIONS], A blank software options screen appears:

- 5.) Select [INSTALL] . A Select Mechanism window opens (see below illustration). Select the mechanism through which you want to install Option Keys.


- 6.) Select [Permanent]. A Select Device window opens (see below illustration). Select the mechanism through which you want to install the Permanent Option Keys.

- 7.) Select the [MEDIA] button and insert the Options DVD or MOD.
- 8.) Select [OK] . Available options appear on the Software Options screen.

Example:  
Options screen  
with Permanent  
& FlexTrial  
options.

- 9.) Select either the DMPR (Direct Multi Planar Reformat) or Variviewer Option installations, if applicable, per Note:

**Note:** **IMPORTANT!**

**DMPR and Variviewer are mutually exclusive options with 05MW14.5 or later software. If you have DMPR (the full option), software will not allow you to install Variviewer, and vice-versa. For example, do not install Variviewer, if you plan to load DMPR.**

- 10.) Perform the following Option installations, if applicable:

**Note:** **The order of installation is critical.**

- A.) VolumeViewer
- B.) CardIQ2

- C.) CardEP
- D.) Conect pro (enable pps)
- E.) Exam Split (VES)
- 11.) Pick remaining options one at a time from the Available Options list and select [INSTALL] to update each selection and place it on the Installed Options List.
- 12.) When the process is complete, select [Quit] then [Quit] to close the window.

**Note:** After Select [Quit] , please wait for 5 minutes before perform any other operation.

## Section 9.0 Load Other Service Software

*(For GE Healthcare Employees or Licensed Customers Only)* If applicable, load Class C Service Software now.

*(For GE Healthcare Employees Only)* If applicable, load Class M Service Software now.

### 9.1 Load Class C

- 1.) Open a Unix Shell.
- 2.) Type: `su -`  
Password: `#bigguy<Enter>`
- 3.) Insert the Class C software CD-ROM into the HP Host CD drive. Wait for the disk read to STOP.
- 4.) Install the Class C software on the system:  
At the superuser prompt, type: `install.classC`  
When the installation process is complete, the Console reports the following:  
`Finished loadingSystem TypeClass C software`
- 5.) Verify that the Class C software is on the system:  
At the superuser prompt, type: `showprods | grep -i classc`  
Verify that the window displays a list of CLASS C file information and locations.

### 9.2 Load Class M

- 1.) Open a Unix Shell.
- 2.) Type: `su -`  
Password: `#bigguy<Enter>`
- 3.) Insert the Class M software CD-ROM into the HP Host CD drive. Wait for the disk read to STOP.
- 4.) Install the Class M software on the system:  
At the superuser prompt, type: `install.classM`  
When the installation process is complete, the Console reports the following:  
`Finished loadingSystem TypeClass M software`
- 5.) Verify that the Class M software is on the system:  
At the superuser prompt, type: `showprods | grep -i classm`  
Verify that the window displays a list of CLASS M file information and locations.

## Section 10.0 Retro Recon Test

The Reconstruction (not Scan) portion of the system can be tested at this point.

- 1.) Start up application  
{ctuser@hostname} st
- 2.) Verify the FSA Box is Idle (not Paused). If paused, click on Recon Management and Restart the Queue.
- 3.) Select [RETRO RECON] on the left Monitor.
- 4.) Select [RAT GOLD] series.
- 5.) Select [SELECT SERIES] .
- 6.) Select [CONFIRM].
- 7.) Verify the image(s) reconstruct and are displayed.
- 8.) Select [QUIT] .

## Section 11.0 Verify Recon Data Path Diagnostic

- 1.) Select the *Service* Icon.
- 2.) Select [DIAGNOSTICS]
- 3.) Select [Recon Data Path]
- 4.) Select [1] and [ALL]
- 5.) Select [RUN]
- 6.) Verify Recon Data Path tests pass.
- 7.) Select [Dismiss]

## Section 12.0 Re-save Reconfig Info File on System State

Save system state with a new DVD-RAM shipped with 16 FPS option. Label the DVD-RAM with "System State – 16 FPS configuration".

- 1.) Insert the System State DVD-RAM into the SCSI Tower DVD RAM drive.
- 2.) Wait until the DVD drive is ready (i.e., until the front panel green LED is no longer lit)
- 3.) Select: *Service icon* to access the CSD (Common Service Desktop)
- 4.) Select: [Utilities]

5.) Select: [System State]


Figure 2-8 System State GUI

- 6.) Select [All]
- 7.) Select [Save]
- 8.) If applicable, select [Yes] when the following message appears: System State Media Status: Please insert a DVD into the drive and press Save again.

Figure 2-9 Save System State Ready

- 9.) When completed, select [CANCEL]

**NOTICE:** Do NOT select [DISMISS] more than once. Doing so will result in the GUI becoming hung, and the System State will be unusable. Should this occur, Re-Save the System State again.

- 10.) When completed, select [Dismiss]
- 11.) Close the *Service Desktop* window at the upper left corner of the screen
- 12.) Install software options

## Section 13.0 Finalization

- 1.) Perform sanity scans (e.g., scout, axial and helical)

**NOTICE:** If the System will not scan or reconstruct verify the memory on the DARC2 Node is correct (Ihinv). Verify the DIP diagnostics pass Verify the RAC diagnostics pass.

- 2.) Verify the image(s) reconstruct and are displayed.
- 3.) Install the front and rear Operator Console covers.

# Chapter 3


## Final Check

### Section 1.0

#### Attaching The Rating Plate

Attach the provided rating plate to the following location (in the label area left to the rating plate of the OC).

OC Rear Side

[www.gehealthcare.com](http://www.gehealthcare.com)

**Imagination at work**