GE Healthcare

# Simple Clamp Replacement Manual

OPERATING DOCUMENTATION

5826675-1EN

Rev 1

© 2019 General Electric Company.

All rights reserved.

---

# **LEGAL NOTES**

---

## **TRADEMARKS**

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

## **COPYRIGHTS**

All Material Copyrighted (c) 2019 by the GE Company, All rights reserved.

This page intentionally left blank.

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## **DAMAGE IN TRANSPORTATION**

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 6.
- Fill out the GIQ workflow for any items missing, damaged, OBF/FOI for in process installs:  
[http://supportcentral.ge.com/ProcessMaps/form\\_new\\_request.asp?prod\\_id=268679&form\\_id=573167&node\\_id=1916016&map\\_id=&reference\\_id=&reference\\_type](http://supportcentral.ge.com/ProcessMaps/form_new_request.asp?prod_id=268679&form_id=573167&node_id=1916016&map_id=&reference_id=&reference_type)
- Contact your local service coordinator for more information on this process.

Rev. Nov. 10, 2017

## **CERTIFIED ELECTRICAL CONTRACTOR STATEMENT**

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

## **IMPORTANT...X-RAY PROTECTION**

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

## LITHIUM BATTERY CAUTIONARY STATEMENTS

### CAUTION

**Risk of Explosion.**

**Danger of explosion if battery is incorrectly replaced.**

**Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**

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## OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the Healthcare PQR Process to report all omissions, errors, and defects in this publication.

---

# *Revision History*

---

<b>Revision</b>	<b>Date</b>	<b>Reason for change</b>
1	July 11, 2019	Initial Release

This page intentionally left blank.

---

# List of Effected Pages

---

PAGES	REVISION	PAGES	REVISION	PAGES	REVISION
1 through 54	1				

This page intentionally left blank.

# Table of Contents

## Preface

**Publication Conventions ..... 19**

**Section 1.0**

**Safety & Hazard Information ..... 19**

1.1 Text and Character Representation..... 19  
1.2 Graphical Representation ..... 20

**Section 2.0**

**Publication Conventions ..... 21**

2.1 General Paragraph and Character Styles..... 21  
2.2 Page Layout..... 21  
2.3 Computer Screen Output/Input Text Character Styles ..... 22  
2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 22

## Chapter 1

**Introduction..... 23**

**Section 1.0 Purpose ..... 23**

**Section 2.0 Applicable Detector and corresponding FRU kit ..... 23**

**Section 3.0 Preliminary Requirements ..... 23**

3.1 Furnished Materials in FRU kit ..... 23  
3.2 Field Supplied Materials ..... 24  
3.3 Optional Tools..... 24

## Chapter 2

**Procedure..... 25**

**Section 1.0 Overview ..... 25**

**Section 2.0 Preparation for Replacement ..... 25**

**Section 3.0 Simple Clamp and Interposer Replacement ..... 27**

3.1 FRDM Removal ..... 30  
3.2 Interposer and Simple Clamp Removal ..... 31  
3.3 Cleaning..... 38  
3.4 Re-assembly Module ..... 42  
3.5 Interconnection check (plenum off, gantry covers off) ..... 48  
3.6 Light seal cover and Plenum Installation ..... 49  
3.7 Detector Testing ..... 49  
3.8 Gantry reassembly..... 49  
3.9 System Calibration..... 49

**Section 4.0 Finalization..... 49**

## Chapter A

**Torque Specification ..... 51**

---

# Preface

## Publication Conventions

---

Purpose: Standardize conventions for representing information in a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

## Section 1.0

### Safety & Hazard Information

#### 1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

#### EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

**DANGER  
SEVERE  
PERSONAL  
INJURY**

**DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

**WARNING  
SERIOUS  
PERSONAL  
INJURY**

**WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

 **CAUTION**  
**Minor Personal Injury**  
Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:

- Loss of critical patient data
- Crush or pinch points
- Sharp objects

 **NOTICE**  
**Property Damage Only**  
Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:

- Disk drive will be destroyed
- Internal mechanical damage, such as to the x-ray tube
- Coasting the rotor through resonance.

## 1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

<p><b>ELECTRICAL</b></p>  <p><b>LASER</b></p>  <p><b>LASER LIGHT</b></p>	<p><b>HAZARDS MECHANICAL</b></p>  <p><b>HEAT</b></p> 	<p><b>RADIATION</b></p>  <p><b>PINCH</b></p> 
<p><b>AVOID STATIC ELECTRICITY</b></p> 	<p><b>PROCEDURAL TAG AND LOCK OUT</b></p> 	<p><b>WEAR EYE PROTECTION</b></p>  <p><b>EYE PROTECTION</b></p>

# Section 2.0 Publication Conventions

## 2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

### EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


**Purpose:** Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

**Note:** Conveys information that should be considered important to the reader.

**Example:** Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

**Comment:** Represents "additional" information that may or may not be relevant to your situation.

## 2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document's part number and revision number appears in every header on every page. Odd

numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

## 2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators "greater than" and "less than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

Example:  
Fixed Output

This paragraph's fonts represents computer generated screen "fixed" output. It's output is fixed from the sense that it does not vary from application to application. It's the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

Example:  
Variable Output

*This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable\_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.*

Example:  
Fixed Input

**This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p**

Example:  
Variable Input

*This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable\_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as ypcat hosts | grep 3.45.120.3 without the greater-than and less-than operators.*

## 2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example:  
Hard Keys

A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example:  
Soft Keys

Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

**END OF SECTION**

# Chapter 1

## Introduction

### Section 1.0 Purpose

Purpose of this kit is to replace Simple Clamp on all detector modules in the detector listed below. Simple Clamp is one of parts which makes interconnection of analog signals between analog module and AD sandwich.

### Section 2.0 Applicable Detector and corresponding FRU kit

Det Part #	Detector Description	Applicable FRU kit part #
5126400	Sherlock Det	5826622SLHL
5196400	Halo Det	
5206400	VCT HD750 DET	
5202200	SAT Det	5826622HLST
5202200-2	SAT Det SDCB	
5202200-3	SAT Det SDCB RoHS	


### Section 3.0 Preliminary Requirements

#### 3.1 Furnished Materials in FRU kit

	GE parts #	5826622SLHL for Sherlock, Halo, HD	5826622HLST for Halo-Saturn
Interposer MP2	5108407	114 pcs	97 pcs
Top Clamp Assembly	5313650	57 pcs	40 pcs
Plastic Top Clamp	5202035	-	17 pcs
Bottom Clamp Assembly	5313653	57 pcs	57 pcs
Torx screw(T10)	5144267	180 pcs	180 pcs
Cloth wiper or cleaning paper Wiper	NA	1 box	1 box
Round Toothpick		1 box (50 pcs)	1 box (50 pcs)

### 3.2 Field Supplied Materials

Item	Description
Anhydrous alcohol to clean FPC/AD board surface	
ESD Mat and Wrist Strap	
Anti-static nitrile gloves	
Standard FE tool kit	
FE Torque Wrench Kit	46-268445G1 or 5135228 or equivalent
Torx Bit – T10  	
Torque driver capable of torquing from 3 to 22 lb-in(0.34Nm ~ 2.5Nm)	As an option, if you had four torque drivers, one set to 3 lb-in and another set to 7 lb-in the process to replace Simple Clamps and interposers would be more efficient. To install the module one torque drive set at 10lb-in(1.1Nm) and another set to 3.5 lb-in(0.4Nm) for the Digital Cable Screw.
HEX BIT 2.5MM with 0.25inch HEX Shank, 1.9inch length  	
HEX BIT 3.0MM with 0.25inch HEX Shank, 1.9inch length  	

### 3.3 Optional Tools

Item	Description
DETECTOR DIGITAL MODULE FIELD REPLACEMENT ACCESSORY COLLECTOR (5133782)	This contains the tools needed to remove HD and Halo detector modules, and is optional as many Field Engineers already have these tools.
DETECTOR DIGITAL MODULE FIELD REPLACE ACCESSORY COLLECTOR (5138154)	This contains the tools needed to remove Saturn detector modules, and is optional as many Filed Engineers already have these tools.
Module alignment gage (5172309) - order from Tool Depot when working with modules 27-31.	Not required reference only
Deionizing fan (2335059)	OPTIONAL

---

# Chapter 2

## Procedure

---

### Section 1.0 Overview

The procedure will include the following steps:

- Preparing for replacement. This step will direct you to prepare a clean, well-organized work area to replace the interposers. ESD mat, wrist strap, and Nitrile gloves are required.
- Opening up the detector to access the modules. You will be directed to the Service documentation for covers and plenum removal instructions.
- Pulling modules and replacing interposers and Simple clamp. If you need clean the Light seal comb in DAS card cage, several adjacent modules should be removed at same time to keep space for cleaning. Be CAREFUL not to damage collimator plate and alignment fixture. Be CAREFUL when pulling out and inserting the module to prevent damage to the digital cable if it rubs against the neighboring module. If three Service Engineers are available, one can pull modules and one remove Clamps and cleaning modules, and one re-assemble modules.
- It is CRITICAL to ensure the modules get returned to the same location.
- Retesting alignment, reference the Service documentation.
- Replacing covers, reference the Service documentation.
- Calibration and Finalization, reference the Service documentation.
- Ensure that you have the tools needed to remove modules. If you do not have the tools to remove the module, contact your regional support leader for the Detector Digital Module Field Replacement Accessory Collector.
- Estimated working hours are 12 hours with three persons.

### Section 2.0 Preparation for Replacement

The following section will prepare the system and a work space for clean interposer replacement.

**NOTICE**

**SAFETY GLASSES MUST BE WORN.**

As with all tasks that require tools, safety glasses must be worn at all times during the FMI.

**NOTICE**

Prevent permanent damage to the static-sensitive boards. Attach an antistatic wrist strap to your wrist and to a bare metal grounding point.

**NOTICE**

Wear Nitrile gloves to protect the detector module from contaminants.

- 1.) Prior to beginning the simple clamp and interposer replacement it is important to check current alignment to establish a baseline for comparison later. Follow the below steps to review the systems current alignment:

Note: If any modules are out of alignment, it is recommended to wait until simple clamp and interposers have been replaced to attempt alignment.

- a.) From the Common Service Desktop, Utilities Menu, Open a Shell and type: `zalignmentWizard` Enter.

Note: `zalignmentWizard` will NOT open from a regular unix shell, it must be opened from the Utilities Menu in CSD.

- b.) The alignment check will run 2 stationary scans and show plots as shown in `#concept_t4y_rln_13b/fig_ucc_ynn_13b`. Processing will take approximately 2 minutes after the scans are run.


Figure 2-1 Alignment check Display

- c.) Save as a copy of the output as a jpg in the `usr/g/scripts` directory.

Note: Refer to CT System Service Methods) FRDM Replacement Procedure, Section 3.7 step 3 for additional information regarding current alignment.

2.) Clean and Organized Workspace Preparation.

**Figure 2-2 ESD wrist strap and Nitrile gloves**

**NOTICE**

Clean workspace area required to reduce possible contaminants.

- a.) Locate a clean surface and place an ESD mat to use as a work area. The detector module is very sensitive to ESD and dust; it needs to be handled with care.
- b.) Clean all work surfaces, including the ESD mat, the work area must start off clean and organized. A comfortable chair and table should be used to ensure a quality interposer replacement.
- c.) Ensure gloves and ESD strap are on and properly connected.
- d.) Gather the following at the workspace.
  - \* T10 torx bit.
  - \* M3 button head screws from kit.
  - \* One Torque driver capable of 3 to 7 lb-in range (.34 to .8 N-m)
  - \* Alcohol Wipes from kit(s).
  - \* Interposers from kit.
  - \* Top Clamps and Bottoms Clamp from kit.

## Section 3.0 Simple Clamp and Interposer Replacement

**DANGER**

**POSSIBLE ELECTRIC SHOCK**

**UNEXPECTED GANTRY ROTATION CAN CAUSE SERIOUS INJURY.**

**NEVER SERVICE THE GANTRY WITH THE AXIAL DRIVE ENABLED. USE THE GANTRY ROTATION LOCK TO LOCK OUT GANTRY ROTATION. SEE SAFETY MENU OF SERVICE DOCUMENT GUI FOR APPROPRIATE PROCEDURES.**

**⚠ DANGER**

**CRUSH HAZARD**

**UNEXPECTED GANTRY ROTATION CAN CAUSE SERIOUS INJURY.**



**NEVER SERVICE THE GANTRY WITH THE AXIAL DRIVE ENABLED. USE THE GANTRY ROTATION LOCK TO LOCK OUT GANTRY ROTATION. SEE SAFETY MENU OF SERVICE DOCUMENT GUI FOR APPROPRIATE PROCEDURES.**

Please follow Equipment Service – Lockout-Tagout-PPE procedures (Direction 534353202EN Rev.11) as part of the FRDM Replacement Procedure.

In the following section, the FRDM replacement procedure in the Service Documents will be referenced to provide detailed instructions for removing and replacing FRDM's.

	<p>Diode and Flex Asm - with "A" SIDE marking,                  lower channel side                  (Halo example)</p>
	<p>A/D Sandwich Asm "A" Side, lower channel side                  No Digital Cable side.                  (Halo example).</p>
	<p>A/D Sandwich Asm "A" Side                  lower channel side                  No Digital Cable side                  Only for Saturn Outer</p>
	<p>Diode and Flex Asm – with "B" SIDE marking                  higher channel side                  (Halo example)</p>

**Table 2-1 Explanation of Components of FRDM**

	<p>A/D Sandwich Asm – “B” side Note Cable on top                  higher channel side                  Digital cable side                  (Halo example)</p>	
	<p>Bottom Clamp – goes on “B” Side</p>	
	<p>Top Clamp, goes on A Side.</p>	
	<p>Top Plastic Clamp, goes on A Side.                  Only for Saturn</p>	
	<p>Interposers</p>	
	<p>1</p>	<p>A Side. Lower channel side(no digital cable side)</p>
	<p>2</p>	<p>B Side. Higher channel side (digital cable side)</p>
	<p>3</p>	<p>Rear (Backplane) Side. (digital cable connection side)</p>
	<p>4</p>	<p>Front (Table) Side. (digital cable handle side)</p>
	<p>5</p>	<p>Rear, High corner</p>

**Table 2-1 Explanation of Components of FRDM**


## 3.1 FRDM Removal

Refer to Direction in 5338260-8EN (LightSpeed 7.X Service Methods) or 5343532-2EN (Revolution Discovery CT Service Methods) or 5409363-8EN (Optima CT660 Service Methods)

FDRM Replacement Procedure (Section 3.4 FRDM Removal) for the modules needing Interposer and or Simple Clamp replacement.

Note: ESD strap and Nitrile gloves are needed for FRDM removal.

- 1.) Follow Section "3.1 Special Tools Overview" in mentioned Service Methods above.
- 2.) Follow Section "3.2 Preparing the System" in mentioned Service Methods above.
- 3.) Follow Section "3.3 Air Plenum Removal" in mentioned Service Methods above.
- 4.) Mark the slot position number on each module, before removing FRDM. All modules must be returned to exactly the same slot.

**Figure 2-3 Numbering on Module**

- 5.) Follow Section "3.4 FRDM Removal" in mentioned Service Methods above.

### 3.2 Interposer and Simple Clamp Removal

- 1.) Move the FRDM module to the clean workspace prepared earlier. Place the FRDM module with the “A” side up. See Figure 4 on page 20, showing the orientation of a Module with the “A” side up.
- 2.) Using a T10 torx bit, loosen and remove three (3) M3 button head screws. Discard these screws. New screws are provided in the kit.

**Saturn outer module case:** Saturn module A side has Plastic Top Clamp. Remove Plastic Top Clamp and discard it.

①

②

1	Halo example
2	Saturn example, including Plastic Top clamp removal


**Figure 2-4 M3 Button Head Screws**

Note: Be careful not to open Flex too much. Handle Flex not to get out from the analog module, with supporting Flex root part.


Note: Be careful not to deform two alignment holes of Flex.

- 3.) Slowly lift up the top clamp and flex (they may stick together), remove the interposer and discard. Be careful not to bend the flex too far, or damage it during removal.

**Saturn outer module case:** Slowly pull down the bottom Clamp to release FPC alignment holes from alignment pins. Saturn module has no interposer at A side.

①

②


1	Halo example
2	Saturn example

**Figure 2-5 Top Clamp and Flex**


**Figure 2-6 Top Interposer Removal (except for Saturn Outer module)**

- 4.) While holding the diode and flex assembly, lift up on the A/D sandwich assembly straight up and remove.

①


②

1	Halo example
2	Saturn example, Top Plastic cover is already removed.

**Figure 2-7 A/D Sandwich Assembly Removal**

5.) Remove the interposer from the bottom clamp and discard, see illustration below.

①


②

1	Halo example
2	Saturn example

**Figure 2-8 Bottom Interposer**

6.) Separate the Flex A side from Top Clamp. If Flex is stickily attached on the clamp, insert tip of toothpick and rotate the toothpick. This operation makes separation easy. Discard the Top Clamp.

**Saturn outer module case:** Saturn module A side has Plastic Top Clamp without elastomer, instead of Top Clamp shown in [Figure 2-4](#). Also, Saturn module has no interposer at A side.

**Figure 2-9 Flex A side separated from Top Clamp**

Note: Be careful not to deform alignment holes of Flex B side. Deformation of the holes may cause open signals.


Figure 2-10 Alignment holes to be protected

- 7.) Flex remount for Area 1  
With using toothpick, separate the B side Flex half part which alignment pins do not exist. "1" of Illustration above.


Figure 2-11 Flex separated from Bottom Clamp, area without alignment pin

8.) Separate the near diode side. "2" of Figure 2-12.


Figure 2-12 Flex separated from Bottom Clamp, area without alignment pin


Figure 2-13 Keeping attach the area around two alignment pins. Area "3" of Figure 2-11

9.) Separate Flex alignment hole area from Bottom Clamp, so that these holes are not damaged.

With using fine stick like a tooth stick.


Figure 2-14 Flex long hole area separation


Figure 2-15 Flex separation with using tooth stick


Figure 2-16 Not to damage center alignment hole


### 3.3 Cleaning

- Note: Do not use disinfectant alcohol which contains water. It may cause leakage issue after reassembling modules. Must use Anhydrous alcohol.
- Note: While cleaning the electronics in the following steps, there may be debris from broken interposer contacts. Ensure to collect and remove any debris found. The gold particles are extremely conductive and could cause significant damage to the electronics.
- Note: Be careful NOT to touch your fingers directly to the contact pads.
- 1.) Carefully clean the contact pads of the A/D sandwich on both sides using the Anhydrous Alcohol and Waste Cloth. If melted elastomer material is attached, it should be cleaned.

Figure 2-17 A/D Sandwich contact pads cleaning (Halo example)

- 2.) If Light seal gasket is dirty with melted material, it should be also cleaned.

**Figure 2-18 Cleaning of Light Seal gasket (Halo example)**

- 3.) Carefully clean both outer sides of Flex using the Anhydrous Alcohol and Waste Cloth. If melted elastomer material is attached, it should be cleaned.

**Figure 2-19 Cleaning of outer side of Flex (Halo example)**


Note: Be careful not to bend (open) the flex too much, it can get damaged.

- 4.) Carefully clean both inner sides of contact pads on the flex using Anhydrous Alcohol and Waste Cloth. If melted elastomer material is attached, it should be cleaned.

**Figure 2-20 Flex pads cleaning**


- 5.) Wait one (1) minute for the alcohol to dry before proceeding.
- 6.) Clean the upper side of Light seal comb in DAS card cage, using the Anhydrous Alcohol and Waste Cloth. If melted elastomer material is attached, it should be cleaned.

**Figure 2-21 Cleaning of upper side of Light Seal comb**

**Note:** Be careful not to touch Collimator plates and Alignment fixture.

- 7.) Clean the lower side of Light seal comb in DAS card cage, using the Anhydrous Alcohol and Waste Cloth. If melted elastomer material is attached, it should be cleaned.


**Figure 2-22** Cleaning of lower side of Light Seal comb

## 3.4 Re-assembly Module

**Note:** It is important that the alcohol dries completely before proceeding. If the alcohol is not dry, it can cause leakage and noise failures.


- 1.) Carefully examine the flex contacts and A/D Sandwich contacts under a light source to ensure that no small particle or residue remain.
- 2.) Attach a new Bottom Clamp to B side of Diode & Flex assembly. Not to deform Flex alignment hole.

**Figure 2-23 Bottom Clamp at B side**

**Note:** Do Not touch the gold contacts on the interposer.


- 3.) Carefully remove one interposer from the package. Only handle the interposer by the ends and absolutely avoid touching the gold contacts on the interposer.

1	Handling area
2	Gold contacts - Do NOT touch.

**Note:** Be careful NOT to touch your fingers directly to the contact pads on the flex.

- 4.) Place the interposer on top of the bottom (B side) flex, the alignment pins on the bottom clamp will guide the interposer into the correct position.

**Figure 2-24 Bottom Interposer Install**

- 5.) Position the A/D sandwich with the ribbon cable (Digital cable) side down and place it on top of the bottom clamp and interposer.

Figure 2-25 AD Sandwich Aligned to Bottom Clamp

Figure 2-26 A/D Sandwich Installed

Note: Remind how to handle Interposer. Refer to 4-4.

- 6.) Carefully take another new interposer from the package and place it top of the A/D sandwich assembly.

**For Saturn:** Interposer at top side (A side) is not needed. Skip this step.

**Figure 2-27 Top Interposer Install**

- 7.) Place the A-side flex over the top Interposer. Not to deform Flex alignment holes.

**For Saturn:** Place the A-side flex over the AD sandwich. No Interposer on top side Saturn.

**Figure 2-28 Top Clamp Reinstalled (Halo example)**

- 8.) Place the Top Clamp over the top Interposer. Not to deform Flex alignment holes.  
**For Saturn:** Use Plastic Top Clamp instead of metal one.

**Figure 2-29 Top Clamp Reinstalled (Halo example)**

- 9.) Install three (3) new M3 button head screws provided in the kit onto the clamp, finger tighten only. Ensure the top and bottom clamps are evenly engaged.

**Figure 2-30 M3 Button Head Screws**

10.) Torque the three (3) M3 button head screws starting at the center (1), then top (2), and then bottom (3). See [Figure 2-31](#) and [Table 2-2](#).


**Figure 2-31 Torque Pattern**

Bottom (3)	Center (1)	Top (2)
3 lbin	7 lbin	3 lbin
0.34Nm	0.8 Nm	0.34Nm

**Table 2-2 Factory Torque Spec for Clamp Screws**

11.) Visually check correct assembling of Simple Clamp and Interposer.

1	Sharlock, Halo, HD
2	Saturn
3	Interposer

**Figure 2-32 Simple Clamp and Interposer layer structure (View from Table side)**

12.) Reinstall the detector module, refer to the FRDM Replacement Procedure (Section 3.5 FRDM Installation)”. Check the ribbon cable on the FRDM so that it is not sticking out too far prior to installing back into the detector.

13.) If there are other “at risk” modules, repeat Steps 1 to 12 above.

### 3.5 Interconnection check (plenum off, gantry covers off)

When all “at risk” modules are completed, check the connectivity of the modules. This can be done once, without having to reinstall the plenum.

Note: The following window will pop up to indicate X-Ray will be created. Ensure that the gantry is locked and Axial Switch is off, the click Confirm.

Note: Be sure to no leave the systems powered up for long without the plenum, this can cause the detector to overheat.

- 1.) Ensure that the rotational lock is engaged.
- 2.) Ensure that the front display is connected to satisfy e-stop.
- 3.) Apply the HVDC and the 120VAC power ONLY. Do NOT turn on the axial drive switch.
- 4.) Perform DasTools to ensure connectivity of analog signals.
- 5.) Power OFF the HVDC and the 120VAC power after tube being proper status for power off.

---

### 3.6 Light seal cover and Plenum Installation

- 1.) Install the light seal cover being careful around thermistor wires, using torx head screws. Torque screws to the value shown the Appendix.
- 2.) Install the Air plenum. Refer to refer to the FDRM Replacement Procedure (Section 3.6 Plenum Installation).

### 3.7 Detector Testing

Perform MiniCheck and others. Refer to the FDRM Replacement Procedure (Section 3.7 Detector Testing).

### 3.8 Gantry reassembly

- 1.) Reinstall the screw cover bands (both bands are identical) and torque the screws to the value shown in the Appendix.
- 2.) Install all gantry covers remembering to enable Axial Drive just before installing gantry right side cover.  
Refer to Replacement → Gantry → Enclosure → (Cover Removal Procedures).

Note: Gantry covers MUST be in place and detector up to proper temperature prior to running full calibration.

### 3.9 System Calibration

Complete the FRDM Replacement Procedure (Section 3.9 System Calibration).

## Section 4.0 Finalization

- 1.) Test scans done as part of FRDM Replacement Procedure.
- 2.) Perform a [Save State] to save the new calibration data.
- 3.) Remove all of the old parts from the site and dispose of it.

*This page intentionally left blank.*

# Appendix A

## Torque Specification

Factory Torque Spec for Detector module Clamp Screws (Flex Connection Screws).

Note: This is different than the specs in the Service Procedure.

Bottom (3)	Center (1)	Top (2)
3 lbin	7 lbin	3 lbin
<b>0.34Nm</b>	<b>0.8 Nm</b>	<b>0.34Nm</b>

1.1 Nm	0.8 lb-ft	10 lb-in	11.5 kg-cm
--------	-----------	----------	------------

**Table 2-3 Detector Module Screw**

0.4 Nm	0.3 lb-ft	3.5 lb-in	4 kg-cm
--------	-----------	-----------	---------

**Table 2-4 Digital Cable Screw**

2.5 Nm	1.8 lb-ft	22 lb-in	25.5 kg-cm
--------	-----------	----------	------------

**Table 2-5 Light Seal Screw**

0.8 Nm	0.6 lb-ft	7 lb-in	8.1 kg-cm
--------	-----------	---------	-----------

**Table 2-6 Screw Cover Band**

19.2 Nm	14.2 lb-ft	170 lb-in	196 kg-cm
---------	------------	-----------	-----------

**Table 2-7 Plenum top and side M8 screws**

8 Nm	6 lb-ft	71 lb-in	82 kg-cm
------	---------	----------	----------

**Table 2-8 Plenum bottom M6 screws**

8 Nm	6 lb-ft	71 lb-in	82 kg-cm
------	---------	----------	----------

**Table 2-9 Fan Plate Screws**

2.4 Nm	1.8 lb-ft	21.2 lb-in	24.5 kg-cm
--------	-----------	------------	------------

**Table 2-10 Cable Safety Covers over the CFC and DHC, for 2006 style Plenum**

*This Page Intentionally Blank.*

© 2019, General Electric Company. All Rights Reserved.

[www.gehealthcare.com](http://www.gehealthcare.com)

**Imagination at work**