# Optima CT540 Performix Tube Demotion Instruction

5779630-1EN

Rev 1

© 2017, GE Company,  
All Rights Reserved.

# **LEGAL NOTES**

---

## **TRADEMARKS**

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

## **COPYRIGHTS**

All Material, Copyright © 2017 by GE Company, Inc. All rights reserved.

*This page intentionally left blank.*

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### DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent write "Damage In Shipment" on ALL copies of the freight or express bill BEFORE delivery is accepted or "signed for" by a GE representative or hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

Fill out a report on <http://egems.med.ge.com/edq/home.jsp>

Contact the local service coordinator.

### CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

### IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

## LITHIUM BATTERY CAUTIONARY STATEMENTS

**CAUTION**  
Risk of  
Explosion

**Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**

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## OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

---

# *Revision History*

---

Revision	Date	Reason for change
1	07/11/2017	Initial release.

*This page intentionally left blank.*

# Table of Contents

---

<b>Chapter 1</b>	
<b>Upgrade Overview .....</b>	<b>17</b>
<b>Section 1.0</b>	
<b>Flowchart .....</b>	<b>17</b>
<b>Section 2.0</b>	
<b>Material List .....</b>	<b>18</b>
2.1    B-Cat Inventory.....	18
2.2    Demotion Kit List (B79582CB).....	18
<b>Chapter 2</b>	
<b>Upgrade Procedure .....</b>	<b>19</b>
<b>Section 1.0</b>	
<b>Record System INFO Information.....</b>	<b>19</b>
<b>Section 2.0</b>	
<b>Save System State .....</b>	<b>19</b>
<b>Section 3.0</b>	
<b>Console Upgrade.....</b>	<b>20</b>
3.1    NIO16 Console (Z820) Removal .....	20
3.2    OpenOC16 Console (Z840) Installation.....	20
3.2.1    Preparative .....	20
3.2.2    Unpack OpenOC16 Console .....	23
3.2.3    GPU Card Installation .....	24
3.2.4    OpenOC16 Console Cable Connections .....	25
<b>Section 4.0</b>	
<b>Tube Replacement .....</b>	<b>30</b>
4.1    Requirements .....	30
4.2    Safety.....	30
4.3    Overview.....	31
4.4    Procedure .....	32
4.4.1    Performix™ Plus Tube Removal.....	32
4.5    Performix™ Ultra Tube Installation .....	40
4.6    LFC Procedure .....	46
4.6.1    Hardware Configuration.....	46
4.6.2    Restore System State.....	46
4.7    Convert Protocol.....	47
4.8    Flash Download.....	47
4.9    Tube Replacement Finalization .....	48
<b>Section 5.0</b>	
<b>ASiR Option Installation .....</b>	<b>50</b>
5.1    Check GPU Card.....	50
5.2    Option Installation.....	50

---

<b>Section 6.0</b>	
<b>IQ Check .....</b>	<b>50</b>
<b>Section 7.0</b>	
<b>Checking and Attaching Labels and Emblems .....</b>	<b>51</b>
<b>Section 8.0</b>	
<b>Finalization .....</b>	<b>51</b>


# Chapter 1

## Upgrade Overview

**Performix Tube Demotion Kit** is only available for Optima CT540 with Performix Plus Tube (LB) and NIO16 Console Z820 Host Computer.

This upgrade instruction is to guide FE how to demote Performix Tube from Plus (LB) to Ultra (BB) and Operator Console from NIO16 with Z820 to OpenOC16 with Z840 on CT-only Optima CT540 system.

### Section 1.0 Flowchart

## Section 2.0 Material List

Verify that the contents of upgrade kit are correct according to the material list below.

### 2.1 B-Cat Inventory

Description	Number	
	ASiR	Non ASiR
IB Tube Demotion Kit	B79582CB	
CT540 TP/SW Collector	B79502CB	
OpenOC16 GEHW	B79532CB	
OpenOC16 ASiR Collector	B79842CB	N/A

**Table 1-1 B-Cat Inventory**

### 2.2 Demotion Kit List (B79582CB)

Mark	Number	Qty.	Description
1	D3885T	1.0	Performix Ultra Tube Assembly ROHS
2	U0203AC	20.0	CLAMP CV_200
3	U0204AC	10.0	CLAMP CV_300
4	5780028	1.0	Caution Label
5	5128609-7	1.0	HV Cable
6	5779630-1EN	1.0	Optima CT540 Performix Tube Demotion Instruction
7	2368768-2	1.0	JEDI TO POWER IF SMART ID CABLE ASSY

**Table 1-2 IB Tube Demotion Kit (5779244)**

# Chapter 2

## Upgrade Procedure

### Section 1.0

#### Record System INFO Information

Write down all of the system INFO information in Console Information Sheets, please refer to the Service Methods CD 5350500-2EN rev20 or later service manual.

Note: ORP type must be recorded in advance according to the following steps.

- 1.) Open a shell.
- 2.) Type: rlogin orp.
- 3.) Type: d 0x10000000
- 4.) The content of 0x10000000 to 0x1000000B is for ORP/SCORP/ORPSV board.

```
*****  
10000000: 3534 3337 3632 312d 3035 3032 ffff 0000 *5437621-0502....*  
10000010: 3534 3834 3639 322d 3003 ffff ffff ffff *5484692-0....*  
*****
```

- 5.) If the content is 5437621, the board is ORPSV, others are SCORP.

### Section 2.0

#### Save System State

Refer to the latest service methods (5350500-8EN) to understand the detailed Save State Procedure and other information saving procedure.

Note: If a Save State was just completed, it can be omitted here.

Note: When performing Save/Restore System State, make sure there is **ONLY** one USB storage device plugged in USB ports of both Console and Media Tower simultaneously, otherwise it may result in incorrect data storage.

- 1.) Insert the System State DVD-RAM media into the peripheral tower drive or USB media in any of the console's USB ports.
- 2.) Select [Service Desktop] -> [Utilities].
- 3.) Select [System State] for DVD-RAM media or [System State - USB] for USB media. The System State Save/Restore screen appears.
- 4.) Select [All] to select all calcs, characterization, etc.
- 5.) Select [Save].
- 6.) If DVD-RAM or USB is ready, select [Yes].
- 7.) When completed select [Dismiss].
- 8.) Remove System State DVD-RAM media from the peripheral tower or USB media from USB port.

## Section 3.0 Console Upgrade

### 3.1 NIO16 Console (Z820) Removal

**NOTICE** Make Sure to follow all Lockout/Tagout requirements while performing this procedure.

- 1.) Shut down the system and power off the Operator Console power switch.
- 2.) Switch OFF the breaker in the Power Distribution Box (hospital supply).

Note: If there is UPS, switch OFF the UPS output.

- 3.) Visually verify proper labeling of each cable and disconnect all console rear cables.
- 4.) Disconnect the cables between the old console and system (Gantry, PDU and local hospital network).
- 5.) Remove the old OC aside.


### 3.2 OpenOC16 Console (Z840) Installation

#### 3.2.1 Preparative

Table 2-1 System Interconnect Cables for Console

No.	Description	Long Cable	Short Cable
53	LVAC Power Cable - PDU to Operator Console	2343531	2343531-2
56	Ground, Raceway to Console	2371450-3	2371450-4
101	Signal Cable - Gantry to Console	5419981	5419981-2
102	Signal Cable (Ethernet) - Gantry to Console	5432020	2373436-3
103	Data Cable (Fiber Optic) - Gantry to Console	5432019	

Figure 2-1 System Interconnect Diagram (with Z840 Host Computer)

### 3.2.2 Unpack OpenOC16 Console

Please read details about OpenOC16 installation procedure in Service Manuals.

- 1.) Unpack console packaging.
  - a.) Prize up buckles by using screwdriver and remove the packaging covers.

**Figure 2-3 OpenOC16 Packaging Box**

- b.) Remove all packaging and discard per local requirement.


**Figure 2-4 Remove all packaging**

- 2.) Move and place the console to the operator room.

3.) Adjust console position, and then pedal four wheel brakes to prevent console from sliding.

**Figure 2-5 Wheel Brake**

### 3.2.3 GPU Card Installation

Note: If customer does not purchase ASiR option, please ignore this section.

Please refer to ASiR Upgrade Installation Manual (5771958-1EN) to install GPU card (shipped with B79842CB) in Z840 Host Computer.

### 3.2.4 OpenOC16 Console Cable Connections

Figure 2-6 Z840 Host Computer Connections



Figure 2-7 OpenOC16 Console Rear View

Number	Description
J1	Scan Monitor Power Connection
J2	Display Monitor Power Connection
J4	In-Room Monitor Connection
J5	Peripheral Media Tower Power Connection
J12	RPM Power Connection
J13	Injector Power Connection
J56	GSCB Power Connection


- 1.) Connect the Keyboard and Mouse (5763460, shipped with B79532CB) to Z840 Host Computer.  
 If the length of keyboard and mouse cables is not enough, add the following extension cables (shipped with B79532CB).

Description	PN#
Keyboard Ext Cable	5431909
Mouse Ext Cable	5458346

Table 2-2 keyboard and Mouse Extension Cables

- 2.) Connect GSCB cable connectors to OpenOC16 AC Box and Z840 Host Computer according to Figure 2-8 and Table 2-8.

**Figure 2-8 GSCB Cable**

#	ITEM	Description
1	DB9/F Black	Host PC DIP Serial Port
2	DB9/F Gray	Host PC RS232
3	Green Audio	Host PC Audio Out (Green)
4	Blue Audio	Host PC Audio In (Blue)
5	USB	Empty
6	DB9/M Black	AC Box J14
7	DB25/M Black	TGP Gantry Cable

**Table 2-3 GSCB Cable Connectors**

- 3.) Connect the media tower cables.
  - a.) Remove the old media power cable (5432953-2).
  - b.) Connect the new media power cable (5432953-6, shipped with B79532CB) to OpenOC16 AC Box J5.
  - c.) Connect the USB cables to Z840 Host Computer.

Description	PN#	Qty
USB Cable (PMT)	5315370	2
Power Cable, Tower to AC Box	5432953-6	1

**Table 2-4 Media Tower Cables**

- 4.) Connect the LCD Monitors for Z840 Host Computer.
  - a.) Remove the old video cables (5408703 and 5332107-2) from the LCD monitors.
  - b.) Remove the old power cables (5432953-4 and 5432953-3) from the LCD monitors.
  - c.) Connect the new video cables and power cables (shipped with B79532CB) to OpenOC16 AC Box and Z840 Host Computer.
    - Scan Monitor  
 Video Cable: from Console Host DP to Monitor DVI  
 Power Cable: from Console AC Box J1 to Monitor

**Figure 2-9 Video Cable and Power Cable for Scan Monitor**

- Image Monitor  
 Video Cable: from Console Host DVI-I to Monitor VGA  
 Power Cable: from Console AC Box J2 to Monitor

**Figure 2-10 Video Cable and Power Cable for Image Monitor**


Description	PN#
Scan Monitor Power Cable	5478299-5
Scan Monitor Video Cable	5408703-2
Image Monitor Power Cable	5478299-6
Image Monitor Video Cable	5366259

**Table 2-5 LCD Monitor Cables**


d.) Route LCD video cable and power cable.

**Figure 2-11 Cable Routing and Keeper**

- 5.) Plug cables into the Switch Hub on console if request.
- 6.) Connect other cables to Z840 Host Computer according to [Figure 2-6](#).

**Figure 2-12 Switch Hub Connections**

- 7.) Connect the OpenOC16 Console power cable and ground cable according to [Figure 2-7](#).

## Section 4.0 Tube Replacement

### 4.1 Requirements

#### Personnel

- 1 person

#### Tools and Test Equipment

- Tube hoist and boom
- HV bleeder kit

#### Consumables

- Transformer oil
- Dry Wipes or Paper Towel

#### Replacement Parts

- Performix Ultra Tube (PN D3885T/D3886T)  
D3885T is Performix Ultra Tube Unit - RoHS  
D3886T is Performix Ultra ECO Tube Unit - RoHS
- Anode HV Cable (5128609-7)
- Smart ID Cable (2368768-2)

### 4.2 Safety

 **WARNING** **POTENTIAL FOR ELECTRIC SHOCK.**  
**LETHAL VOLTAGES MAY BE PRESENT IN THE GANTRY.**  
**FOLLOW PROPER ELECTRICAL LOCK-OUT/TAG-OUT PROCEDURES.**

 **WARNING** **POTENTIAL FOR PERSONAL INJURY (PINCH, CRUSH, ENTANGLEMENT OR IMPACT INJURIES ARE POSSIBLE).**  
**GANTRY HAS ROTATING ASSEMBLY.**  
**USE THE GANTRY'S ROTATIONAL LOCK, TO PREVENT UNCOMMANDED MOTION OF THE GANTRY'S ROTATING ASSEMBLY.**

 **NOTICE** **Potential HV Cable Damage.**  
**Do not sharply bend HV cables, or damage the cables. The minimum allowed radius of bent HV cables is 6.35 cm (2.5 in.).**

 **NOTICE** **Follow ALL required safety and PPE procedures customary for your organization, when working on this product.**

### 4.3 Overview

Figure 2-13 Tube Removal Diagram



Figure 2-14 Performix Ultra Tube and Plus Tube Appearance Difference

Performix Plus Tube

Performix Ultra Tube

## 4.4 Procedure


### 4.4.1 Performix™ Plus Tube Removal

**⚠ DANGER ELECTROCUTION HAZARD  
HIGH VOLTAGES PRESENT.  
USE PROPER LOCKOUT/TAGOUT PROCEDURES BEFORE REPLACING  
COMPONENTS. SEE SAFETY MENU OF SERVICE DOCUMENT GUI FOR  
APPROPRIATE PROCEDURES.**

**⚠ DANGER CRUSH HAZARD.  
UNEXPECTED GANTRY ROTATION CAN CAUSE SERIOUS INJURY.  
USE THE GANTRY ROTATION LOCK TO LOCK OUT GANTRY ROTATION.  
SEE SAFETY MENU OF SERVICE DOCUMENT GUI FOR APPROPRIATE  
PROCEDURES.**


- 1.) Remove and set aside gantry side, top and front covers.
- 2.) Slide out the gantry rear cover, and rotate the top fan assembly outward.

**Figure 2-15 Top Fan Assembly**

- 3.) Remove the M12 screws holding support bracket for the right front gantry cover and set the assembly aside. It may be necessary to tilt the gantry back to remove the third bolt (not normally installed).

**Figure 2-16 Right Front Gantry Cover Support Bracket**

**NOTICE** If the gantry must be tilted, remember to tilt the gantry back to zero degrees.

- 4.) Ensure that Lockout/Tagout procedure has been applied, and that gantry power is removed. Ensure Axial Drive, HVDC, 120VAC are in OFF position.
- 5.) Unlock the Gantry rotation lock, and rotate the Gantry by hand until the tube reaches the 3 o'clock position.

Note: It may be easier to loosen the M12 tube mounting bolts with the tube at about the 2 o'clock position before locking the tube at the 3 o'clock position. Simply loosen the tube mounting bolts one half (1/2) turn. Do not remove the bolts yet.

**Figure 2-17 Tube Angle to Loosen and Torque M12 Screws**

**↑ WARNING**

**CRUSH POINT.**

**IF THE GANTRY IS NOT SECURELY LOCKED, UPON REMOVAL OF THE X-RAY TUBE, THE GANTRY WILL VIOLENTLY ROTATE DUE TO WEIGHT IMBALANCE. ENSURE THE GANTRY IS SECURELY LOCKED BY ATTEMPTING TO ROTATE THE GANTRY BY HAND.**

**↑ NOTICE**

**Make sure the tube is at 90 degrees so the tube hangs at the correct engagement angle for removal and installation.**

- 6.) Engage gantry rotational lock, and check that the gantry is securely locked by attempting to rotate the gantry by hand.

**↑** **WARNING**


**SEVERE INJURY POSSIBLE.**

**IF THE HOIST HOOK'S LOCK NUT IS LOOSE, THEN THE X-RAY TUBE MIGHT DROP FROM THE HOIST, WITH THE RESULTING IN SERIOUS INJURY TO SERVICE PERSONNEL.**

**BEFORE USING THE HOIST, CHECK TO MAKE CERTAIN THAT THE HOOK'S LOCK NUT IS NOT LOOSE. IF IT IS LOOSE, DO NOT USE THE HOIST. INSTEAD, ORDER A NEW HOIST.**

- 7.) Disconnect the 12 pin tube I.D. system cable, from the top of the tube unit J1.
- 8.) Disconnect the 4 pin mate-n-lock pump and fan power system cable, from the top of the tube unit J2.

**Figure 2-18 Tube J1 and J2 Location**


- 9.) Disconnect the ground strap from the top of the tube unit.
- 10.) Remove the anode and the cathode cables from the tube:
  - a.) Carefully cut tie-wraps securing HV cables. Note HV cable routing.
  - b.) Loosen each cable's locking ring with the spanner wrench.
  - c.) Pull each cable terminal out of its receptacle on the tube unit.
  - d.) Remove the anode cable and keep its terminals in JIGs.
  - e.) Ground the end of the cables to the Gantry frame.
  - f.) Wipe up any oil that drips from the cable terminal.
  - g.) Use paper towels to soak up any oil in the wells.

**CAUTION** **Crush Hazard.**  
**Remove the mounting bars in the following (lower/upper) order to lessen the risk of injury to your hand. Keep one hand under the bolt and pressure plate while unfastening it. This is to prevent them from falling into the fan that is attached to the tube.**

Note: It may be easier to tape the socket to the extension. This will prevent the socket from being dislodged on the tube radiator assembly.

11.) Disconnect smart ID cable (5339643-2) connector J1 from Tube.

**Figure 2-19 Smart ID Cable Connector**

12.) Insert the lifting post, boom and chain hoist.

13.) The XRT is attached to the Collimator with a Tube Mount Bracket Assembly. Remove the mounting plate & XRT from the Collimator by removing the four M12 cap screws and lock washers, and two load plates, with a hex socket driver. With your hand, reach behind the radiator to the screws from either side of the XRT center section while removing the bolts with two 12-inch extensions (24 inch length) on a ratchet. Throw these M12 bolts and washers away, as they should not be reused.


14.) Carefully swing the tube clear of the gantry.

Note: Be careful not to damage the fragile copper filter or lead shield in the mounting plate for the next step.

**WARNING** **SEVERE INJURY POSSIBLE.**  
**TUBE CAN SEPARATE FROM GANTRY.**  
**DO NOT REMOVE FACTORY INSTALLED MOUNTING PLATE.**

15.) If the replacement tube has a mounting plate attached (see [Figure 2-20](#)), DO NOT REMOVE IT. Skip step [16.](#)) and step [17.](#))

**Figure 2-20 Tube Shipped with Mounting Plate Attached**

**NOTICE Potential for IQ Artifacts.**

**When removing the mounting plate from the tube, be careful with the Copper Filter. It should be free of debris, scratches and dust. Particles create artifacts in the image by affecting the attenuation properties of the copper filter.**

- 16.) If the replacement tube does NOT have a mounting plate attached, remove the mounting plate by removing the four M10 (P/N 46-328416P20) hex head screws. Throw these bolts and washers away, as they should not be reused.
- 17.) If new tube has no mounting plate attached, inspect the copper filter. The copper filter should be clean, dent and scratch free.
  - If contamination is visible (see [Figure 2-21](#)), clean the collimator.
  - Discoloration is acceptable.

**Figure 2-21 Extremely Contaminated Copper Filter**

Note: Perform the following inspection before installing the new tube unit. Also look at the tube side of the copper filter when you are swapping the interposer plate.

The following tools are required for this inspection procedure:

- \* Phillips #2 Screwdriver
  - \* Flat blade Screwdriver
  - \* Bright Flashlight
- 18.) Inspect the Bowtie Filter. It is possible that the Bowtie Filter is contaminated and the Primary Copper Filter is not contaminated.
- a.) Remove the Collimator Control Board Sheet Metal Cover. See [Figure 2-22](#).

**Figure 2-22 Collimator Control Board Cover**


**⚠ WARNING POTENTIAL FOR SHOCK. VOLTAGE IS PRESENT. DO NOT ATTEMPT TO MOVE THE FILTER WITH POWER ON.**

**⚠ NOTICE Do not force the filter if it feels stuck. Damage to the limit switch can result.**

**Do not move the filter more than necessary for inspection. The filter can fall off the drive screw.**


- b.) With power removed from the gantry, use a flat blade screwdriver and position the bowtie filter assembly so it is visible through the input port or tube side of the collimator assembly.
  - CCW will move the filter into the beam. See [Figure 2-23](#).
  - Do not move the filter back to the home position

**Figure 2-23 Filter Position Adjuster**

- c.) Using a flashlight, inspect the bowtie filter for contamination. Reference [Figure 2-24](#). Look through the input port or tube side of the collimator. If contamination is visible, clean the collimator.

**Figure 2-24 Contaminated Bowtie Filter**

- 19.) Remove the plastic cover from the power I/F board.
- 20.) Disconnect the Smart ID cable (5339643-2) between the JEDI J7 and IF Board J2, and remove it.
- 21.) Get the Smart ID cable (2368768-2) from the upgrade kit, route and fix it, keep sufficient length for cable connection.
  - A1 to JEDI J7
  - A2 to Power IF Board J2

**Figure 2-25 Smart ID Cables Compare**



Figure 2-26 Smart ID Cable (2368768-2) Connection

## 4.5 Performix™ Ultra Tube Installation

**⚠ DANGER CRUSH HAZARD.**  
**UNEXPECTED GANTRY ROTATION CAN CAUSE SERIOUS INJURY.**  
**NEVER SERVICE THE GANTRY WITH THE AXIAL DRIVE ENABLED. USE THE GANTRY ROTATION LOCK TO LOCK OUT GANTRY ROTATION. SEE SAFETY MENU OF SERVICE DOCUMENT GUI FOR APPROPRIATE PROCEDURES.**

**⚠ DANGER ELECTROCUTION HAZARD.**  
**HIGH VOLTAGES PRESENT.**  
**USE PROPER LOCKOUT/TAGOUT PROCEDURES BEFORE REPLACING COMPONENTS. SEE SAFETY MENU OF SERVICE DOCUMENT GUI FOR APPROPRIATE PROCEDURES.**

**⚠ WARNING SEVERE INJURY TO PATIENT COULD OCCUR.**  
**USING THE WRONG BOLTS COULD PLACE STRESS ON BOLTS, CAUSING TUBE TO SEPARATE FROM GANTRY.**  
**USE PROPER BOLTS.**

- 1.) Allow the tube unit to warm to room temperature before you install it.

**↑** **WARNING**

**SEVERE INJURY POSSIBLE.  
TUBE CAN SEPARATE FROM GANTRY.  
DO NOT REMOVE FACTORY INSTALLED MOUNTING PLATE.**

- 2.) If the new tube has a factory installed mounting plate, DO NOT REMOVE IT. Skip to step 4.)
- 3.) Attach the mounting plate from the old tube using new M10 x 25mm bolts provided with the tube. Do NOT use Loctite.
  - Finger tighten all four (4) bolts
  - Torque all four (4) bolts to the following specification:


15 ft-lbs	20 N-m	180 in-lbs	210 kg-cm
-----------	--------	------------	-----------

This seats each bolt, enabling you to visually ensure that the mounting holes are not stripped while applying final torque.

Apply final torque on all four bolts:

30 ft-lbs	41 N-m	360 in-lbs	420 kg-cm
-----------	--------	------------	-----------

**Figure 2-27 Tube and Mounting Plate**

**↑** **NOTICE**

**Potential for IQ artifacts.**

**When attaching the mounting plate on the tube, be careful with the Copper Filter. It should be free of debris, scratches and dust. Particles create artifacts in the image. Use a lint free wipe to clean, if needed.**

- 4.) Re-check facility power and make sure it is off.
- 5.) Use the hoist to lift the new tube unit:
  - a.) Position the tube on the gantry: The “crosses” on the mounting plate and on the collimator should fit in perfectly when the tube is aligned properly.

- Note:
- To ease installation, fasten the top pressure plate to the rotating structure first. Then attach the bottom pressure plate.
  - Use new bolts and washers from tube crate. Make sure to select the proper bolts. There are instructions in the crate, and on the tube itself.
  - b.) Fasten the lower and upper and pressure plates to the rotating structure with the four M12


(50 mm) bolts, and set pre-load torque to:

25 ft-lbs	33 N-m	295 in-lbs	340 kg-cm
-----------	--------	------------	-----------

c.) Set final torque to:


49 ft-lbs	66 N-m	590 in-lbs	680 kg-cm
-----------	--------	------------	-----------

**Figure 2-28 Tube and Mounting Plate**


- 6.) Attach the tube I.D. cable to the 12 pin mate-N-Lock connector (J1) on top of the tube. (See [Figure 2-18](#))
  - 7.) Attach the tube pump and fan power cable to the 4 pin mate-N-Lock connector (J2). (See [Figure 2-18](#))
  - 8.) Fasten the grounding strap to the 1/4-20 ground stud on top of the tube unit.
  - 9.) Remove the plastic cap plug from both cable receptacles on the tube unit.
- Note:** **Take care not to lose the rubber quad rings for the High Voltage cables.**
- 10.) Lightly wet the new rubber quad ring with transformer oil (917).
  - 11.) Return the quad ring to its slot at the top of the receptacle retaining ring.
  - 12.) Pour transformer oil (917) into the receptacles to a depth of 10 mm (0.375 in).
  - 13.) Remove the anode HV cable (5429083-2) from HEMIT tank ST+ terminal.
    - a.) Loosen the anode cable's locking ring with the spanner wrench.

**Figure 2-29 Loosen the Anode Cable Locking Ring**

b.) Pull the anode cable terminal out of its receptacle and keep it in the protective CAP (see [Figure 2-30](#)). This CAP comes from the anode HV cable in demotion kit.

**Figure 2-30 Loosen the Anode Cable Locking Ring**

**Remove and Keep the Cable Terminal**


c.) Wipe up any oil that drips from the cable terminal.

d.) Replenish the well with transformer oil to a depth of 10mm (0.375in).

14.) Install/route the anode HV cable (P/N 5128609-7).

Note: The two anode HV cables have terminals with different orientations to prevent incorrect tube installation.

**Figure 2-31 Old and New Anode HV Cable Compare**

**NOTICE**  Incorrectly routed or secured HV cables will result in damage to the HV cables and/or other parts of the gantry.

15.) Be sure to route the HV cable as shown in [Figure 2-32](#).

**Figure 2-32 HV Cables Properly Routed and Secured**

16.) Align the cable terminal orienting key with the notch in the receptacle.

**NOTICE** Do not over tighten the locking ring. Over tightening can deform the cable plug sealing surfaces, break the oil seal between receptacle and housing, twist the receptacle, and disrupt internal wiring.

- 17.) Slowly insert the cathode and anode cables to engage the connector pins, and seat the cable in the well.
- Tighten the cable locking ring.

**Note:** Use service wrench to turn the threaded ring  $\frac{1}{4}$  to  $\frac{1}{2}$  of a full turn after the gap is closed and friction is felt by hand.

- Rotate the cable strain relief for a clean cable dress.

**NOTICE** IF YOU GET OIL ON YOUR HANDS, WASH THEM NOW.

- 18.) Carefully wipe up all excess oil.
- 19.) Secure the HV Cables using large tie-wraps as shown in [Figure 2-32](#).
- 20.) Disconnect hoist from tube and boom.
- 21.) Remove the post and boom from the gantry.
- 22.) Check for oil leaks:
- Wrap rags or paper towels around the cable horns, and tape them into place.
  - Manually rotate the tube to the 6 o'clock position.
  - Return the tube to the 3 o'clock position
  - Remove the toweling and wipe up all excess oil.
  - Wipe off the cable horns, locking rings and strain reliefs with an alcohol-dampened rag.
  - Repeat with a dry rag.
  - Wrap the cable strain reliefs and locking rings with a single layer of absorbent paper tissue. You can use two inch wide strips cut from a paper napkin.
  - Wrap the bottom edge of the paper around the top end of the cable horn, and tape it in place.
  - Extend the top edge of the paper over the top of the locking ring, and tape it to the plastic cable strain relief.
  - Remove paper after leak check.
  - If oils leaks are found, tighten the locking ring slowly until there is no leakage, paying attention to not over tighten.
- 23.) Bolts checklist.
- New mounting bolts and washers are used for mounting the tube; old mounting bolts and washers are discarded.
  - Proper torque specifications are followed for all fasteners, refer to Torque Wrench Information in 5350500-2EN rev20 or later service manual.
  - Tie-wraps and cables are in place.
- 24.) Restore system power at the main disconnect panel.
- 25.) Turn on gantry 120 VAC, HVDC ENABLE and AXIAL DRIVE ENABLE on the service switch panel. Wait at least 10 minutes to warm up the filament.

## 4.6 LFC Procedure

LFC procedure is needed, refer to LFC Procedure in 5350500-2EN rev20 or later Service Manual.

### 4.6.1 Hardware Configuration

Note: Select the appropriate **Hardware Configuration** for the system according to be recorded in [Section 1.0: Record System INFO Information](#).

Note: Select **[Tube Type]** to **Performix\_Ultra\_(BB)** in **Hardware Setting** window.

Note: Select the correct **ORP Type** according to be recorded in [Section 1.0: Record System INFO Information](#).

### 4.6.2 Restore System State

- 1.) Insert a previously saved system state DVD-RAM media into the peripheral tower drive or USB media in any of the console's USB ports.
- 2.) Wait until the media is ready.
- 3.) Select [Service] -> [Utilities].
- 4.) Select [System State] for DVD-RAM media or [System State - USB] for USB media. The System State Save/Restore screen appears.
- 5.) Select [All].
- 6.) Select [Restore]. The Restore System State box appears.
- 7.) Select [Yes]. A message at the end of the System State Log window should display: Restore System State Completed Successfully.
- 8.) When completed, select [Cancel].

Note: If Does Check option is installed, Dose Audit tool starts up when Reconfig Info or Protocol is selected. That shows configuration and Protocol summary. Press [Quit] to continue.

- 9.) Select [Yes] when the Scan Hardware Reset pop-up appears.
- 10.) When completed, select [Dismiss].
- 11.) Remove System State media and place in a safe location.
- 12.) Verify Restore System State.

## 4.7 Convert Protocol

It's necessary for protocol conversion because the protocol of CTT system is not available in SUSE system.

- 1.) Open a Unix Shell, type:  
{ctuser@hostname} **su -**  
Passwrod: **#bigguy**
- 2.) [root@hostname] # **/usr/g/scripts/conv\_proto** <Enter>
- 3.) Input **y** then press [Enter] to convert protocols.  
*Do you want to convert protocols? [Y/N]*
- 4.) Below information should show if the convert protocols successfully:  
*Success, process finished.*

Note: After finishing the user protocol convert, make sure every protocol can be opened by protocol management and save the current system state to another new media.

## 4.8 Flash Download

Note: The Flash Download takes anywhere from 5-30 minutes, depending on which subsystems require updates.


- 1.) Perform the Flash Download Utility found on the Common Service Desktop - Utilities Tab, select [Flash Download].

Figure 33 Common Service Desktop - Utilities Tab, Flash Download

- 2.) When the Flash Download Window opens, Select [Update].

Figure 34 Flash Download Window5

- 3.) Once the Gantry Hardware Flash Downloads successfully, select [Dismiss].
- 4.) Close the Common Service Desktop.

## 4.9 Tube Replacement Finalization

- 1.) Perform Tube Install Certification, refer to LFC Procedure in 5350500-2EN rev20 or later service manual.
- 2.) Perform tube alignment and calibration:

Note: For detailed procedures, refer to latest service manual for existing system with Performix™ Ultra and JEDI 60DC generator.

- If you have a system with the detector and collimator previously aligned before this tube change and there are no image quality issues, use the “Quick Tube Change” procedure below in the order described.
  - a.) Gantry Rotation Safety Check
  - b.) Reset TnT Data in JEDI Generator Tool
  - c.) Filament Calibration
  - d.) Tube Z-Align
  - e.) ISO Alignment Procedure
  - f.) Meter Verification
  - g.) HV Tank Feedback Resistor Verification (bleeder test)
  - h.) HHS Scans (Adv)  
Tube Usage Information is found on the Common Service Desktop following path:  
Error Logs>Tube Usage>Details.
  - i.) Gantry Balance Procedure, refer to Gantry Balance Procedure in 5350500-2EN rev20 or

- later Service Methods CD
- j.) Hot ISO Alignment
- k.) New Tube Configure
- l.) Z-Slope
- m.) DAS Gain Calibration
- n.) Collimator Calibration
- o.) Detailed Cals
- p.) N Number Check
- q.) System Scanning Test
- r.) Quality Assurance Test - (Need to do EA3 registration as DCadmin and AVexceeding users, detail information refer to 5350500-2EN rev20 or later Service Methods CD.)
- If you cannot follow the quick tube change procedure for any reason, perform the retest steps below in the order described.
  - a.) Install TubeInstallFlex option through eLicense, as needed.
  - b.) Gantry Rotation Safety Check
  - c.) Reset TnT Data in JEDI Generator Tool
  - d.) Filament Calibration
  - e.) Plane of Rotation (POR) Alignment
  - f.) Beam on Window (BOW) Alignment
  - g.) CBF/SAG Alignment Process
  - h.) ISO Alignment Procedure
  - i.) Meter Verification
  - j.) HV Tank Feedback Resistor Verification (bleeder test)
  - k.) HHS Scans (Adv)  
Tube Usage Information is found on the Common Service Desktop following path:  
Error Logs>Tube Usage>Details.
  - l.) Gantry Balance procedure, refer to Gantry Balance Procedure in 5350500-2EN rev20 or later service manual.
  - m.) Hot ISO Alignment
  - n.) Z-Slope
  - o.) DAS Gain Calibration
  - p.) Collimator Calibration
  - q.) Detailed Cals
  - r.) N Number Check
  - s.) System Scanning Test
  - t.) Quality Assurance Test - (Need to do EA3 registration as DCadmin and AVexceeding users, detail information refer to 5350500-2EN rev20 or later Service Methods CD.)

**NOTICE**

**Keep the Anode HV Cable (5429083-2) in safe place at customer side. This cable will be needed/used when this performix Ultra (BB) tube is removed and Performix Plus (LB) tube will be installed in the future.**

---

## Section 5.0

### ASiR Option Installation

Note: If customer does not purchase ASiR option, please ignore this section.

#### 5.1 Check GPU Card

Please refer to ASiR Upgrade Installation Manual (5344359-1EN) to check GPU card.

#### 5.2 Option Installation

- 1.) Un-install Neuro 3D Filter option if needed. (This option conflicts with ASiR Application.)
- 2.) Install the Recon Enhancement Option software (5344374-2, shipped with B79842CB) using Option Install menu.
- 3.) Install ASiR option software (5350779-3, shipped with B79842CB) using Option Install menu.
- 4.) Confirm that system started up normally.

## Section 6.0


### IQ Check

Image Quality Check may include the following:

- 1.) Alignment
  - a.) MTF with GE Performance phantom, both axial and helical images, comparing with spec.
  - b.) Aliasing test, both cold and hot.
  - c.) Z-Beam motion.
- 2.) Noise/Artifacts
  - a.) Noise measurement with 20 cm water phantom, comparing with spec.
  - b.) Cone beam artifact.
- 3.) Microphonics - Check image quality (lack of rings, streak and/or center artifacts).
- 4.) Clever DAS Gain - Confirm that the DAS gains used in the patient scans match those in the cals.
- 5.) CT Number Uniformity - Check both water and poly phantoms' CT number uniformity, comparing with spec.

## Section 7.0 Checking and Attaching Labels and Emblems

- 1.) Turn off Axial Drive and HVDC switches on gantry service panel.
- 2.) Check and attach the provided label.  
**Before attaching label, thoroughly clean the machine surface where to attach the label.**  
 Refer to the table below for proper location and attaching method.

Description	Part Number	Illustration	Location	Qty
CT540 Conversion Kit Locator Card	5780028		Find a proper place to attach the locator card on Gantry Base 	1

- 3.) Re-install Gantry covers.
  - Turn off 120VAC switch on Gantry Service Panel
  - Install Gantry front cover.
  - Install Gantry Left Side cover.
  - Install Gantry top covers.
  - Turn on HVDC, 120V
  - AC, and Axial Drive switches on the Gantry Service Panel. Press E-Reset.
  - Install Gantry Right Side cover.
- 4.) Update GIB with the Product Locator Card information. Go to the following address to access the site: [http://gib.gehealthcare.com/gib/gib\\_entry.jsp](http://gib.gehealthcare.com/gib/gib_entry.jsp)

## Section 8.0 Finalization

- 1.) Insite Check out, refer to 5350500-2EN rev20 or later Service Methods CD.
- 2.) Save System State
  - a.) Insert the Save System State DVD into the Peripheral Tower DVD-RW drive.
  - b.) Select [SERVICE DESKTOP].
  - c.) If reloading software, select [UTILITIES].
  - d.) Select [SYSTEM STATE].
  - e.) Select [ALL] to select all the cals, characterization, etc.
  - f.) Select [SAVE].
  - g.) If applicable, select [OK] when the following message appears: System State Media Status: Please insert a DVD into the drive and press Save again

- h.) When completed select [DISMISS].
  - i.) Close the Service Desktop window in the upper left corner of the screen.
- 3.) Install Gantry side and front covers and keep other components in the kit at customer safe place for further usage in site.

**Components to be kept in customer site:**

- Anode HV Cable
- 4.) Send the following information to [Performix.Ultra@ge.com](mailto:Performix.Ultra@ge.com):
- BB tube serial number
  - BB tube install date, BB tube kAs
  - BB tube exams
  - System ID

[www.gehealthcare.com](http://www.gehealthcare.com)