GE Healthcare

Console Long Cable Kit Installation Manual

OPERATING DOCUMENTATION

5456816-1EN

Rev 5

© 2013-2020 General Electric
Company. All rights reserved.

LEGAL NOTES

TRADEMARKS

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material Copyrighted (c) 2013-2020 by the GE Company, All rights reserved.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 6.
- Fill out the GIQ workflow for any items missing, damaged, OBF/FOI for in process installs:
http://supportcentral.ge.com/ProcessMaps/form_new_request.asp?prod_id=268679&form_id=573167&node_id=1916016&map_id=&reference_id=&reference_type
- Contact your local service coordinator for more information on this process.

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS

 **CAUTION**
Risk of Explosion **Danger of explosion if battery is incorrectly replaced. Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GEMS CQA Process to report all omissions, errors, and defects in this publication.

Revision History

Revision	Date	Reason for change
1	8/23/2013	Initial release
2	4/13/2016	Add Z820 information for NIO16/GOC6.6. Update mouse information for NIO16/GOC6.6 Z800.
3	6/26/2018	Add Z840 information for Open Console.
4	03/01/2019	Add Z8G4 information for Open Console
5	12/16/2020	According to SPRHCSDM00627199 and DOC2441968, update Overview to modify console environment to meet PIM.

Page Left Intentionally Blank

Table of Contents

Preface

Publication Conventions 17

Section 1.0

Safety & Hazard Information 17

1.1 Text and Character Representation..... 17

1.2 Graphical Representation 18

Section 2.0

Publication Conventions 18

2.1 General Paragraph and Character Styles..... 18

2.2 Page Layout..... 19

2.3 Computer Screen Output/Input Text Character Styles 19

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 20

Chapter 1

Introduction & Procedure 21

Section 1.0

Personnel Requirements 21

Section 2.0

Overview 21

Section 3.0

Preliminary Requirements 22

3.1 Applies To:..... 22

3.2 Material List 22

3.3 Safety..... 23

3.4 Required Conditions 23

Chapter 2

Installation Procedure 25

This Page Intentionally Blank.

Preface

Publication Conventions

Purpose: Standardize conventions for representing information in a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

Section 1.0 Safety & Hazard Information

1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

**DANGER
SEVERE
PERSONAL
INJURY**

DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

**WARNING
SERIOUS
PERSONAL
INJURY**

WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

**CAUTION
Minor
Personal
Injury**


Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:

- **Loss of critical patient data**
- **Crush or pinch points**
- **Sharp objects**

- NOTICE** Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:
- Property Damage Only**
- Disk drive will be destroyed
 - Internal mechanical damage, such as to the x-ray tube
 - Coasting the rotor through resonance.

1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

Section 2.0 Publication Conventions

2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


Purpose: Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

Note: Conveys information that should be considered important to the reader.

Example: Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

Comment: Represents “additional” information that may or may not be relevant to your situation.

2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document’s part number and revision number appears in every header on every page. Odd numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that’s either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators “greater than” and “less than” symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. italics) to indicate information specific to your system or site.

Example:
Fixed Output

This paragraph’s fonts represents computer generated screen “fixed” output. It’s output is fixed from the sense that it does not vary from application to application. It’s the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

Example: *This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.*

Example: **This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p**

Example: *This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as ypcat hosts | grep 3.45.120.3 without the greater-than and less-than operators.*

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example: A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example: Whereas the computer **MENU** button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

END OF SECTION

Chapter 1

Introduction & Procedure

Section 1.0 Personnel Requirements



Required Persons	Preliminary Reqs	Procedure	Finalization
1	Timing Not Applicable	30 mins	0 mins

Section 2.0 Overview

This option provides customer a capability to install the NIO / Open console and GOC6.6 console away from operator’s workspace up to approximately 3 meters.

NOTICE When install the console cabinet in the scan room, do not place it within the patient environment. Refer to [Figure 1-1](#) for the details.

Figure 1-1 Patient Environment

IEC 2513/2000

Only Scanning Gantry, Patient Table components and the following options can be placed in this area.

- Advantage 4D
- In room monitor
- SmartStep
- Extream Injector

Section 3.0 Preliminary Requirements

3.1 Applies To:

This option are applies to the following console (Please follow the proper system PCM configuration):

- GOC6.6 Console
- NIO Console
- Open Console

3.2 Material List

Mark	Part Number	Quantity	Description
1	5389625-2	1.0	Long Power Cable For PMT, NIO
2	5389625-4	1.0	Long Power Cable For Left Monitor, NIO
3	5389625-3	1.0	Long Power Cable For Right Monitor, NIO
4	5458246	1.0	Monitor Video Cable DVI-VGA
5	5431912 / 5431912-2	2.0	Monitor Video Cable From DP to DVI
6	5443978	1.0	USB Extension Cable For PMT
7	5431911	1.0	USB Extension Cable For Barcode Reader
8	5458346	1.0	3.5m USB Extension Cable For Mouse
9	5431909	1.0	USB Extension Cable For Keyboard
10	5431910	1.0	USB Extension Cable For Trackball
11	5449053	1.0	Extension Cable For GSCB
12	5443978-2	1.0	USB Extension Cable For EXT HDD
13	5404262	1.0	GSCB Cable
14	5456816-1EN	1.0	Console Long Cable Kit Installation Manual

Table 1-1 5456816 Material List

3.3 Safety

No general safety information.

3.4 Required Conditions

This kit must be installed by field engineer.

Page Left Intentionally Blank

Chapter 2

Installation Procedure

NOTICE Follow ALL Electro-Static Discharge procedures. Always use ESD Wrist Strap and grounding cords.

This section describes the installation procedure of the long cables to the NIO /Open console and GOC6.6 console.


1.) Shutdown the system and power OFF the Operator Console at the front panel switch.

Note: If GSCB with the new Cable, please skip the step 2.)

Figure 2-1 Old Cable and New Cable for GSCB

Old GSCB Cable Without Extension
(5404262 Rev.3)

New GSCB Cable With Extension
(5404262 Rev.4)

- 2.) Replace the GSCB Cable
 - i.) Remove GSCB metal plate by losing the four screws.

Figure 2-2 Metal Plate Removal

- ii.) Remove the 8 screws on GSCB rear cover.

Figure 2-3 Remove Screws from GSCB Rear Cover

- iii.) Carefully open GSCB rear cover, disconnect E-stop cable on PWA J31 and remove the rear cover.


Figure 2-4 Disconnect E-Stop Cable Connector

- iv.) Remove the cable clamp by unscrewing its two screws. (refer to [Figure 2-7](#))


- v.) Carefully use flat-head screw driver to unlock two connectors of GSCB cable on J2 and J38, and unplug them, the recommended method shown as below.

Figure 2-5 Disconnect Cable Connectors J2 and J38


- vi.) Connect new GSCB cable 5404262 (included in the Long Cable Kit 5456816) Conn 1 and Conn 2 (shown as below) to PWA J2 and J38.

Figure 2-6 GSCB Cable 5404262

vii.) Adjust the GSCB cable, let GSCB front cover edge is exactly in the slot on cable boot.

Figure 2-7 Fix the GSCB Cable

viii.) Use cable clamp and the two screws to fix GSCB cable.

NOTICE: Cable clamp should press on cable foil.

ix.) Check the setting of switch S1 and S2 shown as below **System Settings**. See [Figure 2-8](#) and [Table 2-1](#)

Figure 2-8 Switches Setting


Table 2-1 Switches Setting

Switch	Pin	Pin Name	Description	Switch Setting		Note
				Default	Requested on system	
S1	P1	Not Used	N/A			
	P2	BEEPER_ON / OFF	BEEPER_ON: enable beeper function	BEEPER_OFF	BEEPER_OFF	X- ray ON sound can be turned off / on by using S1-2 switch
			BEEPER_OFF: disable beeper function			
	P3	232 / USB	232: Host communicate with GSCB through RS232, for "Gantry Reset" command.	232	232	Wrong setting will cause system reset failure.
USB: Host communicate with GSCB through USB, for "Gantry Reset" command.						
P4	Not Used	N/A				
S2		E-RESET / TILT2	E-RESET: enable GSCB E-RESET function	E-RESET	For system using TGPG E-RESET.	
			TILT2: disable GSCB E-RESET function		For system with other TGPs (TGPHD/ TGPU, etc.): TILT2.	

- x.) Connect GSCB rear cover E-stop connector on PWA J31.

Figure 2-9 Reconnect E-stop Connector

- xi.) Install rear cover by the 8 screws.

Figure 2-10 Rear Cover Installation

- xii.) Install the metal plate by screwing four screws and fix the ground cable.

Figure 2-11 Metal Plate Installing and Ground Cable Connecting

- 3.) Connect the cables for **NIO** console according to the [Table 2-2](#) and [Figure 2-12/Figure 2-13](#).

Connect the cables for **Open** console according to the [Table 2-3](#) and [Figure 2-14](#).
Connect the cables for **GOC6.6** console according to [Table 2-5](#) and [Figure 2-16](#).

Figure 2-12 Cable Connections (NIO16 with Z800)


Figure 2-13 Cable Connections (NIO16 with Z820)


Table 2-2 Cable Connections List for NIO Console

MARK	FROM	LONG CABLE KIT	TO	MARKUP
1	Image Monitor	power cable (5389625-3)	AC Box J9	Replace the old cable with the new one
2	Image Monitor	Video cable (5458246) for Z800	Host Computer DVI	Replace the old cable with the new one
		Video cable (5431912) for Z820	Host Computer DP0	Replace the old cable with the new one
3	Scan Monitor	power cable (5389625-4)	AC Box J10	Replace the old cable with the new one
4	Scan Monitor	Video cable (5431912) for Z800	Host Computer DP	Replace the old cable with the new one
		Video cable (5431912) for Z820	Host Computer DP1	Replace the old cable with the new one
5	PMT Tower	USB ext cable for EXT HDD (5443978-2)	Host Computer USB	Replace the old cable with the new one
6	PMT Tower	USB ext cable for DVD-RW (5443978)	Host Computer USB	Replace the old cable with the new one
7	PMT Tower	PMT power cable (5389625-2)	AC Box J11	Replace the old cable with the new one
8	Keyboard	3.5m USB ext cable for keyboard (5431909)	Host Computer USB	Add the long new cable

MARK	FROM	LONG CABLE KIT	TO	MARKUP
9	Mouse	3.5m PS/2 ext cable (5450275) for Z800	Host Computer PS2	Not in NIO/GOC6.6 console long cable kit (5456816)
		3.5m USB ext cable (5458346) for Z820	Host Computer USB	Add the long new cable

NOTE for Z800:

- 1, PS2 extension cable (5450275) and PS2 mouse (5132092) are shipped with NIO console assembly.
- 2, Customer can order PS2 extension cable (5450275) and PS2 mouse (5132092) to replace USB extension cable (5458346) and USB mouse (5143798-6).
- 3, If Customer has installed PS2 mouse and extension cable, the USB extension cable (5458346) shipped with NIO/GOC6.6 Console Long Cable Kit (5456816) is out of use.

10	Trackball	USB ext cable for trackball (5431910)	Host Computer USB	Add the long new cable
11	Barcode Reader	USB ext cable for barcode reader (5431911)	Host Computer USB	Add the long new cable

NOTE:

For Barcode Reader installation, two connection methods can be used:

- 1.) If single bay tower PMT-20/21/22 (5270510-20/21/22) is used, directly connect the Barcode Reader to the tower USB port, the Barcode Reader extension cable is unused.
- 2.) For other towers, connect the Barcode Reader extension cable between the host computer USB port and the Barcode Reader.

12	GSCB	GSCB ext cable (5449053)	Host Computer, TGP, AC Box	Add the long new cable between Conn 3 and Conn 4 of GSCB Cable 5404262 (Refer to Figure 2-17)
----	------	--------------------------	----------------------------	--

Figure 2-14 Cable Connections (NIO64 with Z820 or NIO/OpenOC with Z840)


Table 2-3 Cable Connections List for NIO64 with Z820 or NIO/OpenOC with Z840

MARK	FROM	LONG CABLE KIT	TO	MARKUP
1	Image Monitor	Power cable (5389625-3)	AC Box J2	Replace the old cable with the new one
2	Image Monitor	Video cable (5458246)	Host Computer DVI	Replace the old cable with the new one
3	Scan Monitor	Power cable (5389625-4)	AC Box J1	Replace the old cable with the new one
4	Scan Monitor	Video cable (5431912-2)	Host Computer DP	Replace the old cable with the new one
5	PMT Tower	USB ext cable for EXT HDD (5443978-2)	Host Computer USB	Replace the old cable with the new one
6	PMT Tower	USB ext cable for DVD-RW (5443978)	Host Computer USB	Replace the old cable with the new one
7	PMT Tower	PMT power cable (5389625-2)	AC Box J5	Replace the old cable with the new one
8	Keyboard	3.5m USB ext cable for keyboard (5431909)	Host Computer USB	Add the long new cable
9	Mouse	3.5m USB ext cable (5458346)	Host Computer USB	Add the long new cable
10	Trackball	USB ext cable for trackball (5431910)	Host Computer USB	Add the long new cable

MARK	FROM	LONG CABLE KIT	TO	MARKUP
11	Barcode Reader	USB ext cable for barcode reader (5431911)	Host Computer USB	Add the long new cable
12	GSCB	GSCB ext cable (5449053)	Host Computer, TGP, AC Box	Add the long new cable between Conn 3 and Conn 4 of GSCB Cable 5404262 (Refer to Figure 2-17)

NOTE:

For USB port detailed connection, pls refer to Port Label on the host computer rear cover.

Figure 2-15 Cable Connections (OpenOCG2 with Z8G4)


Table 2-4 Cable Connections List for OpenOCG2 with Z8G4

MARK	FROM	LONG CABLE KIT	TO	MARKUP
1	Image Monitor	Power cable (5389625-3)	AC Box J2	Replace the old cable with the new one
2	Image Monitor	Video cable (5431912-2)	Host Computer DP	Replace the old cable with the new one
3	Scan Monitor	Power cable (5389625-4)	AC Box J1	Replace the old cable with the new one
4	Scan Monitor	Video cable (5431912-2)	Host Computer DP	Replace the old cable with the new one
5	PMT Tower	USB ext cable for EXT HDD (5443978-2)	Host Computer USB	Replace the old cable with the new one

MARK	FROM	LONG CABLE KIT	TO	MARKUP
6	PMT Tower	USB ext cable for DVD-RW (5443978)	Host Computer USB	Replace the old cable with the new one
7	PMT Tower	PMT power cable (5389625-2)	AC Box J5	Replace the old cable with the new one
8	Keyboard	3.5m USB ext cable for keyboard (5431909)	Host Computer USB	Add the long new cable
9	Mouse	3.5m USB ext cable (5458346)	Host Computer USB	Add the long new cable
10	Trackball	USB ext cable for trackball (5431910)	Host Computer USB	Add the long new cable
11	Barcode Reader	USB ext cable for barcode reader (5431911)	Host Computer USB	Add the long new cable
12	GSCB	GSCB ext cable (5449053)	Host Computer, TGP, AC Box	Add the long new cable between Conn 3 and Conn 4 of GSCB Cable 5404262 (Refer to Figure 2-17)

NOTE:

For USB port detailed connection, pls refer to Port Label on the host computer rear cover.

Figure 2-16 Cable Connections (GOC6.6 Console)


Table 2-5 Cable Connections List for GOC6.6

MARK	FROM	LONG CABLE KIT	TO	MARKUP
1	Image Monitor	power cable (5389625-3)	Rear Bulkhead J2	Replace the old cable with the new one
2	Image Monitor	Video cable (5458246)	Host Computer DVI	Replace the old cable with the new one
3	Scan Monitor	power cable (5389625-4)	Rear Bulkhead J3	Replace the old cable with the new one
4	Scan Monitor	Video cable (5431912)	Host Computer DP	Replace the old cable with the new one
5	PMT Tower	USB ext cable for EXT HDD (5443978-2)	Host Computer USB	Replace the old cable with the new one
6	PMT Tower	USB ext cable for DVD-RW (5443978)	Host Computer USB	Replace the old cable with the new one
7	PMT Tower	PMT power cable (5389625-2)	Rear Bulkhead J5	Replace the old cable with the new one
8	Keyboard	3.5m USB ext cable for keyboard (5431909)	Host Computer USB	Add the long new cable
9	Mouse	3.5m PS/2 ext cable (5450275) for Z800	Host Computer PS2	Not in NIO/GOC6.6 console long cable kit (5456816)
		3.5m USB ext cable (5458346) for Z820	Host Computer USB	Add the long new cable

MARK	FROM	LONG CABLE KIT	TO	MARKUP
------	------	----------------	----	--------

NOTE for Z800:

- 1, PS2 extension cable (5450275) and PS2 mouse (5132092) are shipped with GOC6.6 console peripheral collector.
- 2, Customer can order PS2 extension cable (5450275) and PS2 mouse (5132092) to replace USB extension cable (5458346) and USB mouse (5143798-6).
- 3, If Customer has installed PS2 mouse and extension cable, the USB extension cable (5458346) shipped with NIO/ GOC6.6 Console Long Cable Kit (5456816) is out of use.

10	Trackball	USB ext cable for trackball (5431910)	Host Computer USB	Add the long new cable
11	Barcode Reader	USB ext cable for barcode reader (5431911)	Host Computer USB	Add the long new cable

NOTE:

For Barcode Reader installation, two connection methods can be used:

- 1.) If single bay tower PMT-20/21/22 (5270510-20/21/22) is used, directly connect the Barcode Reader to the tower USB port, the Barcode Reader extension cable is unused.
- 2.) For other towers, connect the Barcode Reader extension cable between the host computer USB port and the Barcode Reader.

12	GSCB	GSCB ext cable (5449053)	Host Computer, TGP, AC Box	Add the long new cable between Conn 3 and Conn 4 of GSCB Cable 5404262 (Refer to Figure 2-17)
----	------	--------------------------	----------------------------	--

Figure 2-17 GSCB Extension Cable Connection

- 4.) Power ON the Operator Console at the console front panel switch.
- 5.) Confirm the functionality GSCB, Monitor and Peripheral Media Tower.

Page Left Intentionally Blank

www.gehealthcare.com

Imagination at work