GE Healthcare

# Xtream Injector Option Installation Manual

OPERATING DOCUMENTATION

5317408-1EN

Rev 9

© 2019 General Electric Company.

All rights reserved.

---

# **LEGAL NOTES**

---

## **TRADEMARKS**

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

## **COPYRIGHTS**

All Material Copyrighted (c) 2019 by the GE Company, All rights reserved.

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## **DAMAGE IN TRANSPORTATION**

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://us44hdd21/sctq/InstallFulfill/InstalFulfillment.htm>
- Contact your local service coordinator for more information on this process.

Rev. Jan. 5, 2005

## **CERTIFIED ELECTRICAL CONTRACTOR STATEMENT**

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

## **IMPORTANT...X-RAY PROTECTION**

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

## LITHIUM BATTERY CAUTIONARY STATEMENTS

### CAUTION

**Risk of Explosion.**

**Danger of explosion if battery is incorrectly replaced.**

**Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**

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## OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the Healthcare PQR Process to report all omissions, errors, and defects in this publication.

# Revision History

Revision	Date	Reason for change
9	April 16, 2019	Added information Medtron Injector
8	November 30, 2015	Added information of E-Box Option.
7	September 17, 2014	Added urlich injector information. Updated Medrad interconnect information.
6	May 17, 2012	Added Brivo CT information.
5	October 24, 2011	Updated Table 1-1 (CiA425 Compliant Injector certified)
4	May 16, 2011	RoHS Flag conversion
3	05/18/10	Added Medrad information
2	03/24/09	Updated Section 1.0, 2.0, 3.0 and 6.0
1	08/28/08	Initial Release

---

# List of Effected Pages

---

PAGES	REVISION	PAGES	REVISION	PAGES	REVISION
1 through 44	9				

# Table of Contents

## Preface

**Publication Conventions ..... 19**

### Section 1.0

**Safety & Hazard Information ..... 19**

1.1 Text and Character Representation..... 19  
1.2 Graphical Representation ..... 20

### Section 2.0

**Publication Conventions ..... 21**

2.1 General Paragraph and Character Styles..... 21  
2.2 Page Layout..... 21  
2.3 Computer Screen Output/Input Text Character Styles ..... 22  
2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 22

## Chapter 1

**Xtream Injector Option Installation ..... 23**

### Section 1.0

**Introduction ..... 23**

1.1 Overview..... 23  
1.2 Prerequisites..... 23

### Section 2.0

**Material List ..... 24**

### Section 3.0

**Interconnect Diagram ..... 25**

3.1 Nemoto Injector ..... 25  
3.2 MEDRAD Injector ..... 26  
3.3 Ulrich Injector..... 28  
3.4 Medtron Injector..... 29

### Section 4.0

**Installation Procedure..... 30**

4.1 Install Injector..... 30  
4.2 Install cables between CT and Injector ..... 30  
    4.2.1 Preparation ..... 30  
        4.2.1.1 Cover Removal (except for Brivo CT) ..... 30  
        4.2.1.2 Cover Removal for Brivo CT ..... 31  
    4.2.2 Install Cables ..... 31  
4.3 Restore Power to the Gantry ..... 32  
4.4 Install Option Software Key ..... 32

### Section 5.0

**Functional Test..... 33**

5.1 Verify Connection between CT and Injector ..... 33  
5.2 Functional Test for Xtream Injector Option ..... 33  
5.3 Functional Test for Enhanced Xtream Injector Option..... 36

---

<b>Section 6.0</b>	
<b>System Level Safety Test.....</b>	<b>39</b>
<b>Section 7.0</b>	
<b>Finish Installation .....</b>	<b>39</b>
7.1 Process Product Locator Card .....	39
<b>Chapter A</b>	
<b>Injector Status Buttons .....</b>	<b>41</b>
1.1 Xtream Injector Status Icons .....	41
1.2 Enhanced Xtream Injector Status Icons .....	42

# Preface

## Publication Conventions

---

Purpose: Standardize conventions for representing information in a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

## Section 1.0

### Safety & Hazard Information

#### 1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

#### EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

**DANGER  
SEVERE  
PERSONAL  
INJURY**

**DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

**WARNING  
SERIOUS  
PERSONAL  
INJURY**

**WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

 **CAUTION**  
**Minor Personal Injury**  
Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:

- Loss of critical patient data
- Crush or pinch points
- Sharp objects

 **NOTICE**  
**Property Damage Only**  
Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:

- Disk drive will be destroyed
- Internal mechanical damage, such as to the x-ray tube
- Coasting the rotor through resonance.

## 1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

<p><b>ELECTRICAL</b></p>  <p><b>LASER</b></p>  <p><b>LASER LIGHT</b></p>	<p><b>HAZARDS MECHANICAL</b></p>  <p><b>HEAT</b></p> 	<p><b>RADIATION</b></p>  <p><b>PINCH</b></p> 
<p><b>AVOID STATIC ELECTRICITY</b></p> 	<p><b>PROCEDURAL TAG AND LOCK OUT</b></p> 	<p><b>WEAR EYE PROTECTION</b></p>  <p><b>EYE PROTECTION</b></p>

# Section 2.0 Publication Conventions

## 2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

### EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


**Purpose:** Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

**Note:** Conveys information that should be considered important to the reader.

**Example:** Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

**Comment:** Represents "additional" information that may or may not be relevant to your situation.

## 2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document's part number and revision number appears in every header on every page. Odd

numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

## 2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators "greater than" and "less than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

Example:  
Fixed Output

This paragraph's fonts represents computer generated screen "fixed" output. It's output is fixed from the sense that it does not vary from application to application. It's the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

Example:  
Variable Output

*This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable\_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.*

Example:  
Fixed Input

**This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p**

Example:  
Variable Input

***This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable\_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as ypcat hosts | grep 3.45.120.3 without the greater-than and less-than operators.***

## 2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example:  
Hard Keys

A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example:  
Soft Keys

Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

**END OF SECTION**

# Chapter 1

## Xtream Injector Option Installation

### Section 1.0 Introduction

#### 1.1 Overview

The Xtream Injector is a start synchronization of CT scanner and Injector. Pressing the Start Scan button makes the scan and injector start simultaneously. This feature is only available with the injector, which supports CANopen Injector Interface Standard (CiA425) and is verified by GEHC.

There are 2 classes in Xtream Injector.

- Xtream Injector, which is the same as Class1 in CiA425, allows only ON/OFF.
- Enhanced Xtream Injector, which is the same as Class4 in CiA425, allows ON/OFF and injector parameters setting from the CT scanner.

CT scanner and injector are operated independently after start of scanning.

#### 1.2 Prerequisites

The following hardware must be installed prior to this option installation.

- Gantry I/F Panels with IPC (Integrated Peripheral Control) board
- CiA425 Compliant Injector certified by GEHC (see [Table 1-1](#))

Type (Continued)	Manufacturer/Model (Continued)
Patient contrast injector: For Xtream Injector option	Nemoto Dual Shot Alpha (CiA425 Class I) / GE 5328194 Nemoto Dual Shot Alpha (CiA425 Class IV) / GE 5328195 Nemoto Dual Shot Alpha 7 (CiA425 Class IV) Nemoto Dual Shot GX (CiA425 Class IV) Nemoto Dual Shot GX 7(CiA425 Class IV) Nemoto A-800 (CiA425 Class IV) / GE 5421942 Nemoto W1000 (CiA425 Class IV) Medrad ISI900 (for Stellant D) (CiA 425 Class IV)/GE 5335919 ulrich ohio M (XD2004) (CiA425 Class I) ulrich CT motion (XD8000) (CiA425 Class I)
Patient contrast injector: For Enhanced Xtream Injector option	Nemoto Dual Shot Alpha (CiA425 Class IV) / GE 5328195 Nemoto Dual Shot Alpha 7 (CiA425 Class IV) Nemoto Dual Shot GX (CiA425 Class IV) Nemoto Dual Shot GX 7(CiA425 Class IV) Nemoto A-800 (CiA425 Class IV) / GE 5421942 Nemoto W1000 (CiA425 Class IV) Medrad ISI900 (for Stellant D) (CiA 425 Class IV)/GE 5335919 ulrich CT motion (XD8000) (CiA425 Class IV) Medtron Accutron CT-D (CiA425 Class IV)

Table 1-1 CiA425 Compliant Injector certified by GEHC (as of April, 2019)

- Note:
- For compatibility information with CT system, refer to the Technical Reference Manual of your product.
  - MEDRAD provides an option to upgrade from previously installed injector to integrated version. Installation of additional hardware and firmware upgrade of injector must be completed by MEDRAD before installation of this option.
  - It is also recommended that both GE service and Injector vendor to be present for system checkout to confirm the equipment work properly.

## Section 2.0 Material List

Part Number	Quantity	Description
5317258	1.0	Injector Power Cable for Nemoto and Ulrich
5506484	1.0	Injector Power Cable for Medrad and Nemoto /w E-Box
	1.0	Option Key DVD (Xtream Injector or Enhanced Xtream Injector)
5317408-1EN	1.0	Instruction
5169456	1.0	Injector Control Cable

**Note:** For VCT and HD, Injector control cable is included in system cable collector.

**Note:** 5506484 will be included in the option kit after 2H of 2014.

# Section 3.0 Interconnect Diagram

## 3.1 Nemoto Injector



Figure 3-1 Nemoto Interconnect Diagram

(\*) Not use one of the three AC output. Band extra outlet somewhere safely.

Figure 3-2 Nemoto Interconnect Diagram (with E-Box option)

### 3.2 MEDRAD Injector


Figure 3-3 Interconnect Diagram (Medrad Stellant D)




Figure 3-4 Interconnect Diagram (Medrad Stellant D with GE Supplied Power Cable)

**Figure 3-5 Interconnect Diagram (Medrad Stellant D (Certegra WS) with GESupplied Power Cable)**

### 3.3 Ulrich Injector

**Note:** ohio M does not have equipotential stud.


**Figure 3-6 Interconnect Diagram with Ulrich supplied Cable**

**CAUTION**

**Model “CT motion” can be applied to patient with and without battery charge.  
Model “ohio M” can only be applied to patient without battery charge.  
Equipotential wire connection must be required when Injector is used with wall power connection for battery charge.**

### 3.4 Medtron Injector

## Section 4.0 Installation Procedure

### 4.1 Install Injector

Installation of the injector is done by Injector vendor or local distributors. If GE service is responsible for the installation of the injector in your country, follow the installation manual shipped with injector.

### 4.2 Install cables between CT and Injector


**DANGER RISK OF ELECTRIC SHOCK. USE PROPER LOCKOUT/TAGOUT PROCEDURES.**

#### 4.2.1 Preparation

##### 4.2.1.1 Cover Removal (except for Brivo CT)



- 1.) Use an 8mm Hex wrench to unlatch the side cover from the front cover.
- 2.) Remove the right side cover by lifting it upward to release the two (2) latches, located on the top edge of the cover. Once removed, the Service Switch Panel should be exposed.
- 3.) Turn OFF all 3 switches (Axial Drive, HVDC, 120 VAC) on the Service Switch Panel.
- 4.) If installed, remove both rear base covers, to gain access to Sub D15 connector on IPC Bracket, as shown in [Figure 4-6](#).

**Figure 4-6 Sub D15 Connector (example of BrightSpeed)**

### 4.2.1.2 Cover Removal for Brivo CT

- 1.) Remove the gantry left side cover.
- 2.) Turn OFF all three switches (HVDC, Axial Drive and 120VAC) on the service switch panel.
- 3.) Remove the right, top and rear covers from the gantry.
- 4.) Gain access to Sub D15 connector on IPC bracket, as shown in [Figure 4-7](#).

Sub D15 Connector

For Brivo CT385, remove the cover on the IPC panel, then connect the Injector DB15 cable here.

Figure 4-7 Sub D15 Connector (for Brivo CT)

### 4.2.2 Install Cables

- 1.) Install Injector power cable and Sub-D Cable 5169456 according to “Interconnect Diagram” ([Figure 3-1](#) through [Figure 3-6](#)).

**WARNING**

- To insure operator and patient safety, Medrad requires the use of earth equipotential cables as shown in the diagram ([Figure 3-3](#)). (This is not necessary if GE supplied power cable (5506484) is used. See [Figure 3-4](#))

**WARNING**

- Connect power cables directly to the console (Nemoto) or the wall outlets (Medrad) as shown. Use of multi-tap powerstrips is prohibited.

### 4.3 Restore Power to the Gantry

- 1.) At the Service Switch Panel. Ensure the 3 switches, (Axial Drive, HVDC, 120VAC) are set to ON.
- 2.) Install all covers.

### 4.4 Install Option Software Key

Install option software key (Xtream Injector or Enhanced Xtream Injector) from Install Options GUI in Common Service Desktop.

For Enhanced Xtream Injector option, the Injector Setup screen appears during the installation. Choose the same pressure limit unit with Injector.

**Figure 4-8 Xtream Injector Set Up Screen**

Note: Xtream Injector and Enhanced Xtream Injector option key cannot be installed together.

## Section 5.0

# Functional Test


Note: The following tests need to be done with injector service personnel.

### 5.1 Verify Connection between CT and Injector

- 1.) Turn on the injector.
- 2.) Confirm the following message is logged in gesyslog,  
`IPC : Information : Validated Injector is connected.`
  - If connection to injector is OK but non-GE certified injector is connected, the following message will be logged.  
`IPC : Information : Connected injector is nothing in validated list.`
  - If none of the above message appears, it means non-compatible injector is connected or there may be hardware problems on Injector or interface (IPC or cable).

### 5.2 Functional Test for Xstream Injector Option

- 1.) Turn on the Injector and select an injection protocol on the Injector Monitor. Keep injector "NOT arming state".
- 2.) On CT Console, select a scan protocol with Xstream Injector OFF in New Patient. Set scan parameters as Axial or Helical on view edit screen.
- 3.) After [Confirm] button becomes selectable, click [Xstream Injector (blue)] button.  
[Xstream Injector (blue)] button for injector should turn to [Xstream Injector (red)].  
Then [Confirm] button is disabled.

**Figure 5-1 View Edit Screen with Xstream Injector Button**

- 4.) Press arming button\* on Injector. [Xstream Injector (red)] button on CTUI should turn to [Xstream Injector (white)]. Then [Confirm] button becomes selectable.

\*On Nemoto DualShot injector, arming button means “Air Check” button.


Figure 5-2 View Edit Screen with Xstream Injector Button

- 5.) Modify injector parameters on Injector Monitor.  
 [Xstream Injector (white)] button on CTUI should turn to [Xstream Injector (red)].  
 Then [Confirm] button is disabled.



Figure 5-3 View Edit Screen with Xstream Injector Button

- 6.) Click [Xtream Injector (red)] button on CTUI.  
 [Xtream Injector (red)] button on CTUI should turn to [Xtream Injector (blue)] and [Confirm] button should become selectable.

**Figure 5-4 View Edit Screen with Xtream Injector Button**

- 7.) Click [Xtream Injector (blue)].  
 [Xtream Injector (blue)] turns to [Xtream Injector (red)].
- 8.) Press arming button on Injector.  
 [Xtream Injector (red)] turns to [Xtream Injector (white)] and [Confirm] button becomes active.
- 9.) Click [Confirm] button on view edit screen.  
 [Scan Start] on SCIM starts to blink.
- 10.) Press [Start] button on Injector.  
 Injection is not started.
- 11.) Press [Scan Start] button on SCIM. (If Scan Start gets Time Out, then click [Resume] button.)  
 Scan and Injection are started simultaneously. (Wait until scan and injection complete.)
- 12.) This completes the Xtream Injector function test.

### 5.3 Functional Test for Enhanced Xtream Injector Option

- 1.) Turn on the Injector.
- 2.) On CT Console, select a scan protocol with Xtream Injector OFF in New Patient.
- 3.) Set scan parameters as Axial or Helical on view edit screen and click [Xtream injector (blue)].


Figure 5-5 View Edit Screen with Xstream Injector Button

4.) The Enhanced Xstream Injector screen appears. Click [Off] button to turn to [On].




Figure 5-6 Enhanced Xstream Injector Screen

5.) Default injector parameters are shown. Modify injector parameters on the screen.

**Figure 5-7 Enhanced Xstream Injector Screen**

- 6.) Click [Accept] button. This will download the modified parameters to Injector. [Xstream injector] button becomes [Xstream Injector (Not Ready)](red).
- 7.) Modify injector parameters on Injector Monitor. Then click [Xstream Injector (red)] button on CT Console.
- 8.) Attention pop-up appears.

**Figure 5-8 Attention Pop-up**

- 9.) Click [OK] on Attention pop-up and then click [Get Current] button on the Enhanced Xstream Injector screen. This will upload the modified parameters on injector to CT.

- 10.) Modify injector parameters on the Enhanced Xtream Injector screen. Click [Accept] button.  
Parameters on CT are downloaded to injector again.  
[Xtream injector] button becomes [Xtream Injector (Not Ready)](red).
- 11.) Press arming button on Injector, and then click [Confirm] and press [Scan Start] on SCIM.
- 12.) Scan and Injection are started simultaneously. (Wait until scan and injection complete.)  
After injection completed, actual data ([Injection Actual Data]) is logged into :  
/usr/g/service/log/gesys\_<system name>.log
- 13.) This completes Enhanced Xtream Injector function test.

## Section 6.0

### System Level Safety Test

Complete the enclosure leakage (patient touch) test according to the Enclosure Leakage (Patient Touch) and System Chassis Leakage Tests procedure in system service methods CD-ROM.

## Section 7.0

### Finish Installation

#### 7.1 Process Product Locator Card

- 1.) Process the product locator card according to the instruction in your country.

# Appendix A

## Injector Status Buttons

There are several status of Xstream Injector.

Note: Definition of “Xstream Injector (red)” is different between Xstream Injector option and Enhanced Xstream Injector option.

### 1.1 Xstream Injector Status Icons





Button Status	Condition
	Xstream Injector is turned off
	Xstream Injector is turned on and injector is ready to start injection
	Xstream Injector is turned on and injector is not connected. Xstream Injector is turned on and injector is not controllable from scanner.
	Xstream Injector is turned on and injector is controllable from scanner, but not ready for injection.

Table A-2 Xstream Injector Status Icons

## 1.2 Enhanced Xtream Injector Status Icons

Button Status	Condition
	Xtream Injector is turned off
	Xtream Injector is turned on and injector is ready to start injection in New Patient
	Xtream Injector is turned on and injector is not connected. Xtream Injector is turned on and injector is not controllable from scanner.
	Xtream Injector is turned on and the parameter injector has and scanner has are different. (Both Ready / Not ready to inject) This condition also happens when setting parameter is in progress.
	Xtream Injector is turned on, parameter injector has and scanner has match, Not ready to inject.

**Table A-3 Enhanced Xtream Injector Status Icons**

© 2019, General Electric Company. All Rights Reserved.

[www.gehealthcare.com](http://www.gehealthcare.com)

Imagination at work