



System ID:  
SR#:  
Date:

### PM Report Schedule W\_HHS

Customer Name: \_\_\_\_\_  
 Service Record #: \_\_\_\_\_ System ID or Serial #: \_\_\_\_\_  
 Customer Equipment ID: \_\_\_\_\_ Service Manual Doc #: \_\_\_\_\_  
 Software Version: \_\_\_\_\_ PM Frequency: \_\_\_\_\_  
 Form Status: Not Complete Start Date: \_\_\_\_\_  
 System Type: \_\_\_\_\_ Country: \_\_\_\_\_

#### Preliminary Tasks - Last PM Inspection:

Last Inspection Date :  Last Schedule Completed:   
 List any Schedule W\_HHS tasks that were completed earlier   
 PM A     PM B     PM C     PM W     Optional HHS Scan Data

### PM Schedule Sys W

#### Initial Procedures

#### E-Stops and Patient Tilt Sensor Inspection

Task	Results	Comments
Test E-stops		
Emergency Off Switch (on wall)		
Test Patient Tilt Sensors		

#### General Console Cleaning & Inspection

Task	Results	Comments
Clean Component Filters, Fans & Grills		
Seismic Anchor Check (if necessary)		

#### Slip Ring & Brush Block Inspection and Maintenance

Task	Results	Comments
Remove Slip Ring Brush Debris		
Slip Ring Tracks Inspected		



System ID:  
SR#:  
Date:

Task	Results	Comments
CleanBrush Blocks		
Brush Tip Inspection		

**Gantry Filter Cleaning**

Task	Results	Comments
Clean Gantry Heater Filter		
Inspect Rotating Path Components & Cables		

**Tube Heat Exchanger & Jedi Fan Inspection**

Task	Results	Comments
Inspect/Clean Tube Heat Exchanger		
Inspect/Clean JEDI Inverter Fan		

**General Table Cleaning & Inspection (for GT1700V & VT2000)**

Task	Results	Comments
Inspect/Clean Table Pan		
Inspect Table Covers		
Inspect Cradle Accessory Attachment		
Check Tape Switches		

**NGPDU**

Task	Results	Comments
Inspect Seismic Anchor Bolts		

**System Options Inspection**

Task	Results	Comments
Nemoto Injector Inspection		

**Gantry Verifications (Only as required by Region or Customer, Use Optional HHS Scan Schedule)**



System ID:  
SR#:  
Date:

Task	Results	Comments
HV Tank Resistor Verification		
Meter Verification		
HHS Scans		

**Quality Assurance Test**

Task	Results	Comments
Scan Window Inspection		
QA Phantom Inspection		

Task	GE Phantom Serial #	Results	Comments
Quality Assurance Test			

Task	Results	Comments
X-ray On Indicators Inspection		

**Finishing Up**

Task	Results	Comments
Save State Completed		
Complete PM Paper Work		

**Appendix: HHS Data**

**mA Shunt Register Accuracy**

**Shunt A - Shunt B = 5Ω +/- 2%**

Task	Readings	Results	Comments
Shunt register value			



System ID:  
SR#:  
Date:

The difference between mA measured (= [measured voltage (mV)] / [shunt resistor value] ) and the displayed mA on console should be within  $\pm 2\%$  of the requested mA.

Requested mA	Measured	Displayed	Delta	Pass/Fail	Comments
50 mA					
200 mA					

**mA Test Points Accuracy**

Calculate difference between TP mA Measured and TP mA Displayed. The Delta should be within  $\pm 4\%$  of the requested mA.

Requested mA	Measured	Displayed	Delta	Pass/Fail	Comments
50 mA					
200 mA					

**kV Test Points Accuracy**

TP Measured kV must be within  $\pm 3\%$  of requested kV.

TP Measured kV and Displayed KV must be within  $\pm 2\%$  of requested kV.

Requested kV	Measured	Displayed	Delta	Pass/Fail	Comments
80 kV					
100 kV					
120 kV					
140 kV					

**HV Tank Feedback Resistor Verification**

To pass the kV check, all four of the following conditions must be met:

- 1) The Measured kV must be within +/- 3% of the Requested kV.
- 2) The kV mA Tool Screen reading must be within +/- 3% of the Requested kV.
- 3) Both measurements (scope measured & screen) must be within +/- 2% of the Requested kV.
- 4) The internal scan timer measurement should be within +/-4%.

Requested kV	On Screen kV	Measured kV	Pass/Fail	Comments
80 kV				
100 kV				
120 kV				
140 kV				
Scan timer value	N/A			

**HHS Scan Data**

Small Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		
80							
80							
80							
80							
80							
80							
80							
80							
80							

**HHS Scan Data**

Small Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		



System ID:  
SR#:  
Date:

100							
100							
100							
100							
100							
100							
100							
100							

**HHS Scan Data**

Small Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		
120							
120							
120							
120							
120							
120							
120							
120							

**HHS Scan Data**

Small Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		



System ID:  
SR#:  
Date:

140							
140							
140							
140							
140							
140							
140							
140							

**HHS Scan Data**

Large Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		
80							
80							
80							
80							
80							
80							
80							
80							

**HHS Scan Data**

Large Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		



System ID:  
SR#:  
Date:

100							
100							
100							
100							
100							
100							
100							
100							

**HHS Scan Data**

Large Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		
120							
120							
120							
120							
120							
120							
120							
120							

**HHS Scan Data**

Large Focal Spot		On Screen Data				Results	Comments
kV	mA	kV	mA	5ms mA	Time		



System ID:  
SR#:  
Date:

140							
140							
140							
140							
140							
140							
140							
140							

Is follow-up work required?

Service Scheduled Date:

This inspection passed and meets all GE Healthcare PM specifications

**Comments:**

Characters Remaining:

**Tools Used:**

Description	Serial Number	Bar Code/ Asset Tag	Cal Due Date	
				<input type="button" value="Add Tool"/> <input type="button" value="Remove Tool"/>



System ID:  
SR#:  
Date:

**GE Representative:**

Number of Signers

Name:

SSO #:

SR #:

Signature

**Customer Signature:**  Yes  Not Required

Date Complete

Form Complete

Save