MAR Option Installation Instruction

Service Documentation

OPERATING DOCUMENTATION

5486358-1EN
Rev 3
2015 General Electric Company. All
rights reserved.

LEGAL NOTES

TRADEMARKS

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material Copyrighted (c) 2015 by the GE Company, All rights reserved.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a GE representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://us44hdd21/sctq/InstallFulfill/InstalFulfillment.htm>
- Contact your local service coordinator for more information on this process.

Rev. Jan. 5, 2005

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The GE Company will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the GE Company, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment. Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS

CAUTION

Risk of Explosion.

Danger of explosion if battery is incorrectly replaced.

Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the PQR Process to report all omissions, errors, and defects in this publication.

Revision History

Revision	Date	Reason for change
3	12/01/15	Chapter 1. Section 3.0: New Memory 6450000-108 instead of old memory 6450000-101
2	09/13/13	Added Preliminary Requirements on Chapter 1
1	06/28/13	Initial Release

Table of Contents

Preface

Publication Conventions 17

Section 1.0

Safety & Hazard Information 17

1.1 Text and Character Representation..... 17

1.2 Graphical Representation 18

Section 2.0

Publication Conventions 19

2.1 General Paragraph and Character Styles..... 19

2.2 Page Layout..... 19

2.3 Computer Screen Output/Input Text Character Styles 20

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 20

Chapter 1

Introduction..... 21

Section 1.0

Overview 21

Section 2.0

Preliminary Requirements..... 21

Section 3.0

Material List 22

Chapter 2

MAR Option Installation Procedure 23

Section 1.0

Memory Replacement Procedure 23

1.1 Preparation 23

1.2 Power-Off (Shutdown) 23

1.3 Memory Replacement..... 24

1.4 Verification 28

Section 2.0

MAR Option Installation..... 29

Section 3.0

Finalization..... 30

Preface

Publication Conventions

Purpose: Standardize conventions for representing information is a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

Section 1.0

Safety & Hazard Information

1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

 **DANGER
SEVERE
PERSONAL
INJURY**

DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

 **WARNING
SERIOUS
PERSONAL
INJURY**

WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

 **CAUTION**
Minor Personal Injury
Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:

- Loss of critical patient data
- Crush or pinch points
- Sharp objects

 **NOTICE**
Property Damage Only
Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:

- Disk drive will be destroyed
- Internal mechanical damage, such as to the x-ray tube
- Coasting the rotor through resonance.

1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

<p>ELECTRICAL</p>  <p>LASER</p>  <p>LASER LIGHT</p>	<p>HAZARDS MECHANICAL</p>  <p>HEAT</p> 	<p>RADIATION</p>  <p>PINCH</p> 
---	---	---

<p>AVOID STATIC ELECTRICITY</p> 	<p>PROCEDURAL TAG AND LOCK OUT</p> 	<p>WEAR EYE PROTECTION</p>  <p>EYE PROTECTION</p>
--	---	---

Section 2.0 Publication Conventions

2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


Purpose: Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

Note: Conveys information that should be considered important to the reader.

Example: Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

Comment: Represents "additional" information that may or may not be relevant to your situation.

2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document's part number and revision number appears in every header on every page. Odd

numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators "greater than" and "less than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

Example:
Fixed Output

This paragraph's fonts represents computer generated screen "fixed" output. It's output is fixed from the sense that it does not vary from application to application. It's the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

Example:
Variable Output

This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.

Example:
Fixed Input

This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p

Example:
Variable Input

This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as ypcat hosts | grep 3.45.120.3 without the greater-than and less-than operators.

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example:
Hard Keys

A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example:
Soft Keys

Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

END OF SECTION

Chapter 1

Introduction

Section 1.0

Overview

This procedure describes the MAR Option Installation procedure, includes:

- Memory Replacement
- MAR Option Installation

Section 2.0

Preliminary Requirements

MAR Option Installation only available on SUSE Operating System.

Please confirm the system information before install the MAR Option:

- If the OS is **GEHC-SLES-11-SP1**, please follow this instruction for MAR option installation.
- If the OS is not **GEHC-SLES-11-SP1**, Please follow up IB Upgrade Instruction (5490735-1EN) to upgrade CTT system to SUSE system first, then follow this instruction for MAR option installation.

Section 3.0

Material List

Refer to Table 1-1 for listing of the MAR Option Kit for GEHW.

MARK	PART NUMBER	QUANTITY	DESCRIPTION
1	5479514	1.0	MAR option, GEHW
2	5486307	1.0	LABEL FOR MAR OPTION
3	6450000-108	8.0	4GB DDR3 ECC unbuffered DIMM 1600MHZ or higher frequency
4	5486358-1EN	1.0	MAR Option Installation Instruction

Table 1-1 MAR Option Kit for GEHW (5486358)

Refer to Table 1-2 for listing of the MAR Option Kit for WSO.

MARK	PART NUMBER	QUANTITY	DESCRIPTION
1	5479514-2	1.0	MAR option, WSO
2	5486307-2	1.0	LABEL FOR MAR OPTION
3	6450000-108	8.0	4GB DDR3 ECC unbuffered DIMM 1600MHZ or higher frequency
4	5486358-1EN	1.0	MAR Option Installation Instruction

Table 1-2 MAR Option Kit for WSO (5486358-2)

Chapter 2

MAR Option Installation Procedure

Section 1.0

Memory Replacement Procedure

1.1 Preparation

Confirm that System installation has been completed.

Tools: Field ESD Kit is required

1.2 Power-Off (Shutdown)

- 1.) Select one of the following methods to Power OFF the Operator Console if the Power is ON.
 - A.) If Applications are up, click on the [Shut Down] button and select [Shutdown].
 - B.) If Applications are down, open a Unix Shell from pop-up menu. Enter (type) the following command, at the indicated prompts.

```
ctuser@hostname: su -
```

```
Password: #bigguy
```

```
root@hostname: halt
```

- 2.) Power OFF the console at the front panel switch when 'System Halted' message is displayed.

1.3 Memory Replacement

NOTICE Follow ESD Management and Device Handling procedure.

Figure 2-1 Side Cover

- 1.) Remove the console left side cover by unscrewing its four screws.
 - a.) Remove the side access panel by unscrewing its three screws shown as [Figure 2-2](#).

Figure 2-2 Side Access Panel

b.) Pull out the latch to release the left side cover shown as [Figure 2-3](#) .

Figure 2-3 Left Side Cover

2.) Remove the Airflow Guide and Expansion Card Support shown as below.:

a.) Pull the two handle of right sides shown as [Figure 2-4](#).

Airflow Guide

Expansion Card Support

Figure 2-4 Remove the Protect Covers

- 3.) Remove the Memory Fan as shown [Figure 2-5](#).

Press green clamps together by hand and pull out the fan from the memory socket


Figure 2-5 Memory Fan removal

- 4.) Release two clamps both side of memory, and carefully pull the memory out of the slot (total 2Gx6), as shown [Figure 2-6](#)


Figure 2-6 Memory removal

- 5.) Install the memory with new 4Gx8. Align the memory keyway with the slot key, and firmly seat the memory in the slot as shown [Figure 2-7](#).


Figure 2-7 New Memory Installation

The memory sequence on mainboard is refer to [Figure 2-8](#) drawing on left cover.




Figure 2-8 Memory sequence

1.4 Verification

- 1.) Startup OC, press F10 entering BIOS.
- 2.) Select "System Information" below "File" menu.

- 3.) Confirm that memory size (total) is 32768 MB

- 4.) Select "Save Changes and Exit" to exit the BIOS.
- 5.) Reinstall all removed components, and reconnect any cables that have been disconnected.
- 6.) Reinstall the Airflow Guide and Expansion Card Support as shown [Figure 2-4](#).
- 7.) Reinstall the computer left side cover and the side access panel. Reinstall the console left side cover by screwing its four screws.

Section 2.0 MAR Option Installation

Confirm that System started up normally.

Then install the option using Configuration->Install Option menu, and check the option is installed from Configuration->Verify Install Option menu

Section 3.0 Finalization

- 1.) Execute sanity scan and confirm the scan completed and image is generated without error.
- 2.) Take the label (5486307 or 5486307-2) from the kit and paste the label near the rating plates on the rear cover of NIO console.


Figure 2-9 Paste Label

www.gehealthcare.com