IB Upgrade from CTT to SUSE Instruction

Service Documentation

OPERATING DOCUMENTATION

5490735-1EN
Rev 6
© 2013-2016 GE Company. All
rights reserved.

LEGAL NOTES

TRADEMARKS

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

COPYRIGHTS

All Material Copyrighted (c) 2016 by the GE Company, All rights reserved.

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DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a GE representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://us44hdd21/sctq/InstallFulfill/InstalFulfillment.htm>
- Contact your local service coordinator for more information on this process.

Rev. Jan. 5, 2005

CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The GE Company will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the GE Company, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

LITHIUM BATTERY CAUTIONARY STATEMENTS

CAUTION

Risk of Explosion.

Danger of explosion if battery is incorrectly replaced.

Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.

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OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GE PQR Process to report all omissions, errors, and defects in this publication.

Revision History

Revision	Date	Reason for change
6	04/05/16	Added Service Pack information on section 4.4.3 on Chapter 2
5	01/30/14	Updated Section 4.3 Configure System Information on SUSE Added Section 4.4.3 for SP1.5 on Chapter 2
4	11/27/13	Added additional procedure on Chapter 2 section 4.4.1 and 4.4.2 due to convert tool can't convert all user's protocol of CTT system
3	10/17/13	Updated Convert Protocol procedure per SP1.3 Added Finalization after IB Upgrade
2	09/10/13	Add the convert protocol procedure on chapter 2.0, section 4.4
1	08/29/13	Initial Release

Table of Contents

Preface

Publication Conventions 17

Section 1.0

Safety & Hazard Information 17

1.1 Text and Character Representation..... 17
 1.2 Graphical Representation 18

Section 2.0

Publication Conventions 19

2.1 General Paragraph and Character Styles..... 19
 2.2 Page Layout..... 19
 2.3 Computer Screen Output/Input Text Character Styles 20
 2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 20

Chapter 1

Introduction..... 21

Section 1.0

System Requirement..... 21

Section 2.0

Overview and Flowchart..... 21

Chapter 2

Upgrade Procedure 23

Section 1.0

Preparative..... 23

Section 2.0

Save System State 23

Section 3.0

Record System Information 28

3.1 Record PPS Server Information 28
 3.2 Record All Remote Host Network Information 29
 3.3 Record Product Network Filters (PNF) Information 31
 3.4 Record AW IP for Direct Information 33
 3.5 Record Option Information Licenses 34

Section 4.0

Load From Cold Procedure..... 35

4.1 LFC Procedure 35
 4.2 Restore System State..... 35
 4.3 Configure System Information on SUSE 36
 4.4 Convert Protocol..... 39
 4.4.1 Convert Protocol on Service Pack 1.2 39
 4.4.2 Convert Protocol on Service Pack 1.3 or 1.4..... 40

4.4.3 Convert Protocol on Service Pack 1.5 or later	41
Section 5.0	
Finalization	41

Preface

Publication Conventions

Purpose: Standardize conventions for representing information is a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

Section 1.0

Safety & Hazard Information

1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

**DANGER
SEVERE
PERSONAL
INJURY**

DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

**WARNING
SERIOUS
PERSONAL
INJURY**

WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

 **CAUTION**
Minor Personal Injury
Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:

- Loss of critical patient data
- Crush or pinch points
- Sharp objects

 **NOTICE**
Property Damage Only
Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:

- Disk drive will be destroyed
- Internal mechanical damage, such as to the x-ray tube
- Coasting the rotor through resonance.

1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

<p>ELECTRICAL</p>  <p>LASER</p>  <p>LASER LIGHT</p>	<p>HAZARDS MECHANICAL</p>  <p>HEAT</p> 	<p>RADIATION</p>  <p>PINCH</p> 
<p>AVOID STATIC ELECTRICITY</p> 	<p>PROCEDURAL TAG AND LOCK OUT</p> 	<p>WEAR EYE PROTECTION</p>  <p>EYE PROTECTION</p>

Section 2.0 Publication Conventions

2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


Purpose: Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

Note: Conveys information that should be considered important to the reader.

Example: Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

Comment: Represents "additional" information that may or may not be relevant to your situation.

2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document's part number and revision number appears in every header on every page. Odd

numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators "greater than" and "less than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

Example:
Fixed Output

This paragraph's fonts represents computer generated screen "fixed" output. It's output is fixed from the sense that it does not vary from application to application. It's the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

Example:
Variable Output

This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.

Example:
Fixed Input

This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p

Example:
Variable Input

This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as ypcat hosts | grep 3.45.120.3 without the greater-than and less-than operators.

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example:
Hard Keys

A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example:
Soft Keys

Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

END OF SECTION

Chapter 1

Introduction

Section 1.0

System Requirement

IB upgrade from CTT to SUSE kits are available for the following system configuration:

- Discovery CT590 RT with NIO Console
- Optima CT580 with NIO Console


Section 2.0

Overview and Flowchart

This upgrade instruction describes the procedure of CTT upgrade to SUSE. Manually record the following information due to this information is not carried forward by restoring system state.

- PPS Server Info
- All remote Host Network Information
- PNF Information
- AW IP for direct connect
- Option Information Licenses

The illustration below shows the main flow of the upgrade procedures.

Chapter 2

Upgrade Procedure

Section 1.0

Preparative

Note: Refer to the latest service methods (5366639-8EN) about LFC procedure to understand detailed save state procedure and other information saving procedure.

Section 2.0

Save System State

Before proceeding with Save System State:

- Make sure to have formatted DVD-RAM disk. Make sure the DVD-RAM disk is not write-protected.

Or

- Have a FAT32 Formatted USB media (either Memory Stick or USB External Drive, minimum of 2 GB).


Note: When performing Save/Restore System State, check to ensure there is no two or more USB storage devices (SSA Key or other mobile storage devices) plugged in Console/Tower at the same time which may result in incorrect data storage.

Note: System State is a collection of system files saved in a directory structure. The system uses a default directory name (*service_mod_data*) on System State backup media. Whenever System State is saved to DVD-RAM or USB media, the existing default directory (*service_mod_data*) is erased when Save *All* is selected. To avoid loss of a saved System State, on USB media for example, rename the default directory with a unique name. (ex: System ID). When a specific System State backup is being accessed for restoring System State, the directory must be renamed back to the default directory name of *service_mod_data*.

Due to the default directory name and the directory size, it will not be practical to store multiple System States on a single USB Flash Memory Stick. For emergency backup, it's recommended that System State backups be transferred to a laptop or external USB Hard Drive. Use USB Flash Memory Stick only for the immediate saving and restoring of System State during Load from Cold.

- 1.) Insert the System State DVD-RAM Media into the Peripheral Tower DVD-RAM Drive or USB Media in the console's USB ports.
- 2.) Wait until the DVD-RAM Drive or USB Media is ready (i.e., front panel DVD-RAM Drive LED or USB Device Access LED is no longer flashing).
- 3.) Select: [Service] icon to access the Common Service Desktop (CSD).
- 4.) Click [OK] as required.
- 5.) Select: [Utilities] tab.

Figure 2-1 Common Service Desktop (CSD), Utilities Tab


- 6.) Select: [System State] for DVD-RAM Media or [System State - USB] for USB Media. The System State Save/Restore screen appears.

Figure 2-2 Save System State Window (DVD-RAM Media)

Figure 2-3 Save System State Window (USB Media)

- 7.) Select [All].
- 8.) Select [Save]. The System State Media Ready pop-up window appears.
- 9.) Select [Yes] in the System State Media Ready pop-up window.

Figure 2-4 System State DVD Media Ready Pop-up Window




Figure 2-5 System State USB Media Ready Pop-up Window

- 10.) Verify that the "Save" of System State was successful. A message at the end of the System State Message Log Window should state: Save/Restore System State: Completed Successfully.

Figure 2-6 System State Message Log Window

- 11.) When completed, select [Cancel], then [Dismiss].
- 12.) Close the Common Service Desktop window at the upper left corner of the screen.
- 13.) Remove System State DVD-RAM disk from the DVD Peripheral Tower or USB Media from USB port. Label System State Backup Media appropriately and place in a safe location.

Section 3.0


Record System Information

3.1 Record PPS Server Information


- 1.) Open a Unix Shell:
`{ctuser@hostname} installhisris`
- 2.) Press <ENTER>, the following illustration appeared, record the data then press [Accept].

- 3.) Press [Yes].


4.) Record PPS data then press [Accept].


3.2 Record All Remote Host Network Information


1.) Select [ImageWorks] from the desktop.

2.) Select [Network]

- 3.) Click [Select remote host] on the Network list menu
- 4.) Select Remote host from the Remote host list, then click [update].

5.) Record the Network information as below.

The screenshot shows a dialog box titled "Remote host parameters". It contains several input fields and options:

- Host name: bj05
- Network address: 3.36.231.115
- Network protocol: DICOM
- Port number: 4006
- AE Title: bj05
- Comment: No Comment


Below these fields are three sections of options:

- Archive Node: Auto Yes No
- Access to local host:
 - send images? Yes No
 - query/retrieve images? Yes No
- Custom search? On Off


At the bottom of the dialog are two buttons: "Save" and "Cancel".

3.3 Record Product Network Filters (PNF) Information

- 1.) Select the [SERVICE DESKTOP].
- 2.) Select in sequence: Configuration → Configure PNF


3.) Record PNF information.


3.4 Record AW IP for Direct Information

- 1.) Select the [SERVICE DESKTOP].
- 2.) Select in sequence: Configuration → Configure AW DirectConnect and record AW IP address


- 3.) Select in sequence: Configuration → Configure AW Server Client and record AW server INFO

3.5 Record Option Information Licenses

- 1.) Select the [SERVICE DESKTOP].
- 2.) Select in sequence: Configuration → Install Options

- 3.) View Option List and check whether include below 5 options, these five options need apply the new e-licence when CTT upgrade to SUSE.
 - a.) Advanced Vessel Analysis Xpress on OC
 - b.) AutoBone Xpress on OC
 - c.) Perf 4 - Multi-Organ on the OC
 - d.) Perf 4 - Neuro on Xstream OC
 - e.) Advantage CTC Pro for OC

Option Name on CTT system	Option Name on SUSE system
Advanced Vessel Analysis Xpress on OC	AVA_Xpress
AutoBone Xpress on OC	AutoBone_Xpress
Perf 4 - Multi-Organ on the OC	CT_Perfusion_4D_MultiOrgan
Perf 4 - Neuro on Xstream OC	CT_Perfusion_4D_Neuro
Advantage CTC Pro for OC	CT_Colono_Pro3D_EC

Table 2-1 Compare Option Name between CTT and SUSE system

Section 4.0

Load From Cold Procedure

4.1 LFC Procedure

Detail LFC procedure refer to service methods (5366639-8EN)->Software->LFC procedure (13BW16.x).

Restore System State, please refer to section [4.2](#)

4.2 Restore System State

- 1.) Insert a previously saved System State DVD-RAM Media into the Peripheral Tower DVD-RAM Drive or USB Media in the USB ports
- 2.) Wait until the DVD-RAM Drive or USB Media is ready (i.e., front panel DVD-RAM Drive LED or USB Device Access LED is no longer flashing).
- 3.) Select the [Service] icon to access the CSD (Common Service Desktop).
- 4.) Select: [Utilities].
- 5.) Select: [System State] for DVD-RAM Media or [System State - USB] for USB Media. The System State Save/Restore screen appears.
- 6.) Select [All].
- 7.) Select [Restore]. The Restore System State box appears.
- 8.) Select [Yes].

Note: If the site installed above five options, refer to [section 3.5 on page 34](#). Error Message will be pop-up during this time, refer to below illustration, please click [ok] to continue. Option Name different between CTT and SUSE system, refer to [Table 2-1](#).

- 9.) Verify that the "Restore" of System State was successful. If not, restore the System State again. A message at the end of the System State Log window should display: Restore System State Completed Successfully.
- 10.) When completed, select [Cancel].

Note: Manually install the e-licence for five options.

Note: If Dose check option is installed, Dose audit tool starts up when Reconfig Info or Protocol is selected. That shows configuration and Protocol summary. Press Quit to continue.

- 11.) Select [Yes] when the Scan Hardware Reset pop-up appears.
- 12.) When completed, select [Dismiss].
- 13.) Remove System State DVD-RAM disk from the DVD Peripheral Tower or USB Media from USB port and place in a safe location.


- 14.) Reboot the Operators Console by selecting [Shutdown] Desktop, then select [Restart] in the Attention Window.


- 15.) Verify Restore System State
 - a.) Visually verify that the system has the correct configuration and preference settings by opening a Terminal Window.
 - b.) Execute the system configuration script "reconfig".

4.3 Configure System Information on SUSE


- 1.) Configure PPS server information. Please refer to [section 3.1 on page 28](#)
- 2.) Configure all remote network information
 - a.) Select [ImageWorks] from the desktop.


- b.) Select in sequence: Tools-> Network Configuration


c.) Select [Add]

d.) Input the Host Name then press enter.

- 3.) Configure PNF information
 - a.) Select the Service Desktop
 - b.) Select in sequence: Configuration -> Configure PNF

- 4.) Configure AW IP for direct information, Please refer to [section 3.4 on page 33](#)

4.4 Convert Protocol


Upgrade system from CTT to SUSE is done, it's necessary for FE to do convert protocol because SUSE system can't use the protocol from CTT.

Check the software version string by following methods:

- A.) Open a Unix Shell and type the following command:

```
{ctuser@hostname} swhwinfo
```
- B.) Look at the Common Service Desktop Home Page.

Figure 2-7 Common Service Desktop

- If the system has installed the SP1.2, please follow the step 4.4.1.
- If the system has installed SP1.3 or SP1.4, please follow the step 4.4.2.
- If the system has installed the SP1.5 or later, please follow the step 4.4.3

4.4.1 Convert Protocol on Service Pack 1.2

Insert convert tool CD into Host PC driver, and follow up below steps:

- 1.) In the Unix Shell, become root:

```
{ctuser@hostname} su -
```

 Password: #bigguy
- 2.)

```
[root@hostname]# mount /dev/sr0 /mnt/cdrom <Enter>
```
- 3.)

```
[root@hostname]# cd /mnt/cdrom <Enter>
```
- 4.)

```
[root@hostname cdrom]# ./conv_proto <Enter>
```

- 5.) Input **y** then press **Enter** to convert protocols
Do you want to convert protocols ?[Y/N] (default: N)_Y
- 6.) Below information should show if the convert protocols successfully:
Success, process finished.
- 7.) Below command is for ejecting CD from Host PC driver:
[root@hostname cdrom]# cd ~ <Enter>
[root@hostname]#umount /mnt/cdrom <Enter>
[root@hostname]#eject
- 8.) Open a NEW shell, then execute following commands:
 - a.) {ctuser@hostname} cd /usr/g/protocols/scanRx <Enter>
 - b.) {ctuser@hostname}sed -i '/frameOfReferenceUID/d' *.* <Enter> (no blank space between * and .)
 - c.) {ctuser@hostname}cd Service-4 <Enter>
 - d.) {ctuser@hostname}sed -i '/frameOfReferenceUID/d' *.* <Enter> (no blank space between * and .)

After finishing the user protocol convert, make sure every protocol can be opened by protocol management and save current system state to another new disk. The convert tool or command above will never be used again after upgrading from CTT to SUSE.

Note: If protocol convert is failed, restore previous protocol of system state and perform the steps from [4.4.1](#).

4.4.2 Convert Protocol on Service Pack 1.3 or 1.4

- 1.) In the Unix Shell, become root:
{ctuser@hostname} su -
Password: #bigguy
- 2.) [root@hostname]# /usr/g/scripts/conv_proto <Enter>
- 3.) Input **y** then press **Enter** to convert protocols
Do you want to convert protocols ?[Y/N] (default: N)_Y
- 4.) Below information should show if the convert protocols successfully:
Success, process finished.
- 5.) Open a NEW shell, then execute following commands:
 - a.) {ctuser@hostname} cd /usr/g/protocols/scanRx <Enter>
 - b.) {ctuser@hostname}sed -i '/frameOfReferenceUID/d' *.* <Enter> (no blank space between * and .)
 - c.) {ctuser@hostname}cd Service-4 <Enter>
 - d.) {ctuser@hostname}sed -i '/frameOfReferenceUID/d' *.* <Enter> (no blank space between * and .)

After finishing the user protocol convert, make sure every protocol can be opened by protocol management and save current system state to another new disk. The convert tool or command above will never be used again after upgrading from CTT to SUSE.

Note: If protocol convert is failed, restore previous protocol of system state and perform the steps from [4.4.2](#).

4.4.3 Convert Protocol on Service Pack 1.5 or later

- 1.) In the Unix Shell, become root:

```
{ctuser@hostname} su -  
Password: #bigguy
```
- 2.)

```
[root@hostname]# /usr/g/scripts/conv_proto <Enter>
```
- 3.) Input **y** then press **Enter** to convert protocols

```
Do you want to convert protocols ?[Y/N] (default: N)_Y
```
- 4.) Below information should show if the convert protocols successfully:

```
Success, process finished.
```

After finishing the user protocol convert, make sure every protocol can be opened by protocol management and save current system state to another new disk. The convert tool or command above will never be used again after upgrading from CTT to SUSE.

Note: If protocol convert is failed, restore previous protocol of system state and perform the steps from [4.4.3](#).

Section 5.0 Finalization

Perform *Quality Assurance Test* of Functional Checks on Service Methods, make sure all of them pass.

www.gehealthcare.com