# OC Upgrade from TIO to NIO for Discovery CT590RT and Optima CT580

## Service Documentation

OPERATING DOCUMENTATION

5782217-1EN  
Rev 2  
© 2017 GE Company.  
All rights reserved.

---

# **LEGAL NOTES**

---

## **TRADEMARKS**

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

## **COPYRIGHTS**

All Material Copyrighted (c) 2017 by the GE Company, All rights reserved.

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## **DAMAGE IN TRANSPORTATION**

All packages should be closely examined at time of delivery. If damage is apparent, have notation "damage in shipment" written on all copies of the freight or express bill before delivery is accepted or "signed for" by a GE representative or a hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 6.
- Fill out the GIQ workflow for any items missing, damaged, OBF/FOI for in process installs: [http://supportcentral.ge.com/ProcessMaps/form\\_new\\_request.asp?prod\\_id=268679&form\\_id=573167&node\\_id=1916016&map\\_id=&reference\\_id=&reference\\_type](http://supportcentral.ge.com/ProcessMaps/form_new_request.asp?prod_id=268679&form_id=573167&node_id=1916016&map_id=&reference_id=&reference_type)
- Contact your local service coordinator for more information on this process.

## **CERTIFIED ELECTRICAL CONTRACTOR STATEMENT**

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

## **IMPORTANT...X-RAY PROTECTION**

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The GE Company will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the GE Company, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

## LITHIUM BATTERY CAUTIONARY STATEMENTS

### CAUTION

**Risk of Explosion.**

**Danger of explosion if battery is incorrectly replaced.**

**Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**

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## OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the GE PQR Process to report all omissions, errors, and defects in this publication.

*This Page Intentionally Blank.*

---

# Revision History

---

Revision	Date	Reason for change
1	08/23/2017	Initial release.
2	10/10/2017	Update Section 3.6 Record Option Information Licenses. Update Patch SW version from SP1.2 to SP1.3. Update Section 7.6 Old Console Disposal.

*This Page Intentionally Blank.*

# Table of Contents

## Preface

**Publication Conventions ..... 17**

**Section 1.0**

**Safety & Hazard Information ..... 17**

1.1 Text and Character Representation..... 17

1.2 Graphical Representation ..... 18

**Section 2.0**

**Publication Conventions ..... 19**

2.1 General Paragraph and Character Styles..... 19

2.2 Page Layout..... 19

2.3 Computer Screen Output/Input Text Character Styles ..... 20

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 20

## Chapter 1

**Introduction..... 21**

**Section 1.0**

**System Requirement..... 21**

**Section 2.0**

**Overview ..... 21**

**Section 3.0**

**Flowchart ..... 22**

**Section 4.0:**

**Material List ..... 23**

4.1 B-Cat Inventory..... 23

## Chapter 2

**Upgrade Procedure ..... 25**

**Section 1.0:**

**Preparative..... 25**

**Section 2.0:**

**Save System State ..... 25**

**Section 3.0:**

**Record System Information ..... 30**

3.1 Check TGP Board Type..... 30

3.2 Record PPS Server Information ..... 30

3.3 Record All Remote Host Network Information ..... 32

3.4 Record Product Network Filters (PNF) Information ..... 34

3.5 Record AW IP for Direct Information ..... 36

3.6 Record Option Information Licenses ..... 37

<b>Section 4.0:</b>		
<b>TIO Console Removal and NIO Console Installation.....</b>		<b>38</b>
4.1	Old TIO Console Removal .....	38
4.2	New NIO Console (Z840) Installation.....	38
4.2.1	Unpack NIO16 Console .....	38
4.2.2	GPU Card Installation .....	40
4.2.3	Memory Installation .....	40
4.2.4	Console Cable Connections.....	40
4.2.5	Console Cover Installation .....	48
<b>Section 5.0</b>		
<b>TGP Board Replacement.....</b>		<b>49</b>
<b>Section 6.0:</b>		
<b>Load From Cold Procedure .....</b>		<b>49</b>
6.1	LFC Procedure.....	49
6.1.1	LFC .....	49
6.1.2	Patch Installation.....	51
6.2	Restore System State .....	52
6.3	Convert Protocol .....	53
6.4	Configure System Information on SUSE.....	54
<b>Section 7.0:</b>		
<b>Finalization .....</b>		<b>56</b>
7.1	Flash Download .....	56
7.2	GPU Card Checking.....	58
7.3	System Checking .....	58
7.4	Save System State.....	58
7.5	Label Attaching .....	59
7.6	Old Console Disposal.....	59

# Preface

## Publication Conventions

Purpose: Standardize conventions for representing information in a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

## Section 1.0

### Safety & Hazard Information

#### 1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

#### EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

**DANGER  
SEVERE  
PERSONAL  
INJURY**

**DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

**WARNING  
SERIOUS  
PERSONAL  
INJURY**

**WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

 **CAUTION**  
**Minor Personal Injury**  
**Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:**

- Loss of critical patient data
- Crush or pinch points
- Sharp objects

 **NOTICE**  
**Property Damage Only**  
**Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:**

- Disk drive will be destroyed
- Internal mechanical damage, such as to the x-ray tube
- Coasting the rotor through resonance.

## 1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

<p><b>ELECTRICAL</b></p>  <p><b>LASER</b></p>  <p><b>LASER LIGHT</b></p>	<p><b>HAZARDS MECHANICAL</b></p>  <p><b>HEAT</b></p> 	<p><b>RADIATION</b></p>  <p><b>PINCH</b></p> 
<p><b>AVOID STATIC ELECTRICITY</b></p> 	<p><b>PROCEDURAL TAG AND LOCK OUT</b></p> 	<p><b>WEAR EYE PROTECTION</b></p>  <p><b>EYE PROTECTION</b></p>

# Section 2.0 Publication Conventions

## 2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

### EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


**Purpose:** Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

**Note:** Conveys information that should be considered important to the reader.

**Example:** Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

**Comment:** Represents "additional" information that may or may not be relevant to your situation.

## 2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document's part number and revision number appears in every header on every page. Odd

numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

## 2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators "greater than" and "less than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

Example:  
Fixed Output

This paragraph's fonts represents computer generated screen "fixed" output. It's output is fixed from the sense that it does not vary from application to application. It's the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

Example:  
Variable Output

*This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable\_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.*

Example:  
Fixed Input

**This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: cd /usr/3p**

Example:  
Variable Input

*This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable\_input>. In these cases, the (<>) operators would be dropped prior to input. For example: ypcat hosts | grep <3.45.120.3> would be typed into the computer as ypcat hosts | grep 3.45.120.3 without the greater-than and less-than operators.*

## 2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

Example:  
Hard Keys

A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

Example:  
Soft Keys

Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

**END OF SECTION**

---

# Chapter 1

## Introduction

---

### Section 1.0

#### System Requirement

**The Upgrade Kit** is available for Discovery CT590 RT and Optima CT580 with TIO Console. SmartView option and SmartView 2D/3D option (replacement of SmartView option after the upgrade) is also NOT available to APP 16BW32.x with SP1.2. If the SmartView is already installed before the upgrade, please DON'T purchase and install this upgrade kit.

### Section 2.0

#### Overview

The upgrade instruction is to guide FE how to upgrade Console from TIO to NIO.


If the site need to upgrade from CTT to SUSE, please manually record the following information due to this information is not carried forward by restoring system state.

- PPS Server Info
- All remote Host Network Information
- PNF Information
- AW IP for direct connect
- Option Information Licenses

# Section 3.0 Flowchart

The illustration below shows the main flow of the upgrade procedures.

**Figure 1-1 Flowchart**

# Section 4.0: Material List

Verify that the contents of upgrade kit are correct according to the material list below.

## 4.1 B-Cat Inventory

Number	Description
B79422CB	NIO16 GEHW with Z840
B79442CB	GPU Fluoro for Z840
B70112DA	Z840 Upgrade Collector
B75362CB	Media Tower
B79432CB	TP_SW_Option Collector Z840
B77052CA	GSCB

**Table 1-1 B-Cat Inventory (B70102DA)**

Number	Description
5782217-1EN	IB Upgrade Instruction for RT Freedom
5410711	Signal Cable, Fiber Optical
5370000	TGPHD Board Assy
2265406-2	FIBER-OPTIC COUPLING RECEPTACLE KIT, 2.5MM BAYONET, AMP OPTIMATE--RoHS
5445679	ARMADILLO LOOPBACK MULTIMODE 62.5 125um LCPC at 850NM VCSEL LASER
5575700-299	RT Innovation II_User Interface Quick Guide_Multi-language

**Table 1-2 Z840 Upgrade Collector (B70112DA)**

*This Page Intentionally Blank.*

---

# Chapter 2

## Upgrade Procedure

---

### Section 1.0: Preparative

Refer to the latest service methods (5366639-8EN) to understand detailed hardware and software installation procedure.

### Section 2.0: Save System State

Before proceeding with Save System State:


- Make sure to have formatted DVD-RAM disk. Make sure the DVD-RAM disk is not write-protected.

Or

- Have a FAT32 Formatted USB media (either Memory Stick or USB External Drive, minimum of 2 GB).
- 1.) Insert the System State DVD-RAM Media into the Peripheral Tower DVD-RAM Drive or USB Media in the console's USB ports.
  - 2.) Wait until the DVD-RAM Drive or USB Media is ready (i.e., front panel DVD-RAM Drive LED or USB Device Access LED is no longer flashing).
  - 3.) Select: [Service] icon to access the Common Service Desktop (CSD).
  - 4.) Click [OK] as required.


- 5.) Select: [Utilities] tab.

**Figure 2-1 Common Service Desktop (CSD), Utilities Tab**


- 6.) Select: [System State] for DVD-RAM Media or [System State - USB] for USB Media. The System State Save/Restore screen appears.

**Figure 2-2 Save System State Window (DVD-RAM Media)**


**Figure 2-3 Save System State Window (USB Media)**


- 7.) Select [All].
- 8.) Select [Save]. The System State Media Ready pop-up window appears.
- 9.) Select [Yes] in the System State Media Ready pop-up window.

**Figure 2-4 System State DVD Media Ready Pop-up Window**

**Figure 2-5 System State USB Media Ready Pop-up Window**

- 10.) Verify that the "Save" of System State was successful. A message at the end of the System State Message Log Window should state: Save/Restore System State: Completed Successfully.

Figure 2-6 System State Message Log Window

- 11.) When completed, select [Cancel], then [Dismiss].
- 12.) Close the Common Service Desktop window at the upper left corner of the screen.
- 13.) Remove System State DVD-RAM disk from the DVD Peripheral Tower or USB Media from USB port. Label System State Backup Media appropriately and place in a safe location.

## Section 3.0: Record System Information


### 3.1 Check TGP Board Type

- 1.) Open a Unix Shell.
- 2.) Type: `rlogin tgp`, [Enter].
- 3.) Type: `boardShow`, [Enter].
- 4.) TGP board type will be displayed, **TGPU** or **TGPHD** (example *TGPHD* as below shown).


```
{ctuser@ct29}rlogin tgp
-> boardShow
Board: TGPHD
value = 6 = 0x6
-> █
```

### 3.2 Record PPS Server Information


- 1.) Open a Unix Shell:  
`{ctuser@hostname} installhisris`
- 2.) Press <ENTER>, the following illustration appeared, record the data then press [Accept].

3.) Press [Yes].

- 4.) Record PPS data then press [Accept].


### 3.3 Record All Remote Host Network Information


- 1.) Select [ImageWorks] from the desktop.

2.) Select [Network]

- 3.) Click [Select remote host] on the Network list menu
- 4.) Select Remote host from the Remote host list, then click [update].

5.) Record the Network information as below.

Remote host parameters

Host name	bj05
Network address	3.36.231.115
Network protocol	DICOM
Port number	4006
AE Title	bj05
Comment	No Comment

Archive Node:  Auto  Yes  No

Access to local host:

send images?  Yes  No

query/retrieve images?  Yes  No

Custom search?  On  Off


Save Cancel

### 3.4 Record Product Network Filters (PNF) Information

- 1.) Select the [SERVICE DESKTOP].
- 2.) Select in sequence: Configuration → Configure PNF


3.) Record PNF information.


## 3.5 Record AW IP for Direct Information

- 1.) Select the [SERVICE DESKTOP].
- 2.) Select in sequence: Configuration → Configure AW DirectConnect and record AW IP address


- 3.) Select in sequence: Configuration → Configure AW Server Client and record AW server INFO

### 3.6 Record Option Information Licenses

- 1.) Select the [SERVICE DESKTOP].
- 2.) Select in sequence: Configuration → Install Options.

- 3.) View Option List and check whether include below 5 options, these five options need apply the new e-licence when system platform upgrade from CTT to SUSE.
  - AutoBone Xpress on OC
  - Advanced Vessel Analysis Xpress on OC
  - Perf 4 Neuro on Xstream OC
  - Perf 4 Multi Organ on Xstream OC
  - Advantage CTC Pro for OC
  - SmartView Multi-Slice

B-CAT	Option Name (Installed in Site)	Flex String (CTT System)	Flex String (SUSE System)
B7870JD	AutoBone Xpress on OC	AutoBone	AutoBone_Xpress
B7870JC	Advanced Vessel Analysis Xpress on OC	AdvVesselAnalysis	AVA_Xpress
B7870HN	Perf 4 Neuro on Xstream OC	CTPerfusion4Neuro	CT_Perfusion_4D_Neuro
B7870HM	Perf 4 Multi Organ on Xstream OC	CTPerfusion4	CT_Perfusion_4D_MultiOrgan
B7864KL	Advantage CTC Pro for OC	CTColonoPro	CT_Colono_Pro3D_EC
B7868FM (Smartview 2D)	SmartView Multi-Slice	RealTimeCTFluoro	Enhanced_Real_Time_CT_Fluoro

**Table 2-1 Compare Flex String Between CTT and SUSE System**

## Section 4.0: TIO Console Removal and NIO Console Installation

### 4.1 Old TIO Console Removal

- 1.) Shut down the system and power off the Operator Console power switch.
- 2.) Switch OFF the breaker in the Power Distribution Box (hospital supply).

Note: If there is UPS, switch OFF the UPS output.

- 3.) Visually verify proper labeling of each cable and disconnect all console rear cables.
- 4.) Remove the old Console aside.

### 4.2 New NIO Console (Z840) Installation

#### 4.2.1 Unpack NIO16 Console

Please read the details about NIO16 installation procedure in Service Manuals.

- 1.) Unpack console packaging.
  - a.) Remove all items from the console.
  - b.) Remove all packaging and discard per local requirement.
  - c.) Place the step-board under the front edge of the skid and step on it to raise the front edge of the skid as in [Figure 2-7](#).

**Figure 2-7 Step-Board Used to Raise Front Edge of Skid**

d.) Remove the two front cushions from the bottom of the skid.

**Figure 2-8 Cushion on Botton of Skid**


e.) Lift up on the strap on the front of the step-board (Figure 2-7) to lower the skid. Remove the step-board.

f.) Ensure the console stabilizers are in line with thee notched portion at the front of the skid. This will allow enough clearance to smoothly roll the console down the ramps.

2.) Move and place the console to the operator room.


3.) Adjust two rubber footer to fasten the console for preventing slide.

**Figure 2-9 Cushion on Botton of Skid**

- 4.) Remove the Host PC support bracket from left side of Console.

**Figure 2-10 Host PC Support Bracket**

#### 4.2.2 GPU Card Installation

Please refer to *GPU Fluoro Prerequisite Installation Manual (5344368-1EN\_r5 or later)* to install GPU Card (8770000-126, shipped with B79442CB) in Z840 Host Computer.

#### 4.2.3 Memory Installation

Note: If the customer has purchased MAR Option and installed in old system, please order MAR Option Kit (B79602CB) for new system.

Please refer to MAR Option Installation Manual (5767680-1EN) to install Memory (8770000-100, shipped with B79602CB) in Z840 Host Computer.

#### 4.2.4 Console Cable Connections

NO.	Description	Long Cable	Short Cable
53	LVAC Power Cable - PDU to Operator Console	2343531	2343531-2
56	Ground, Raceway to Operator Console	2371450-3	2371450-4
101	Signal Cable - Gantry to Operator Console	2333150	2333150-2
102	Signal Cable (Ethernet) - Gantry to Operator Console	2352714-2	2352714-3
103	Data Cable (Fiber Optic) - Gantry to Console	2117848-2	2117848-7

**Table 2-2 System Cables for OC Connections**

Figure 2-11 System Interconnect Diagram (with Z840 Host Computer)


Figure 2-12 Console Interconnect (with Z840 Host Computer)




Figure 2-13 Z840 Port Label

**Figure 2-14 NIO16 Console AC Box Assignments**


Number	Description
J9	Display Monitor Power Connection
J10	Scan Monitor Power Connection
J11	Peripheral Media Tower Power Connection
J12	In-Room Monitor Connection
J13	Injector Power Connection
J14	RPM Power Connection
J56	GSCB Power Connection

1.) Connect the Keyboard and Mouse to NIO16 Console Z840 Host Computer.

DESCRIPTION	PART NUMBER	QTY
USB Ext Cable (Mouse)	5458346	1
USB Ext Cable (Keyboard)	5431909	1

- 2.) Remove the old SCIM from the site, and take out the GSCB from the upgrade kit. Connect GSCB cable to NIO16 Console AC Box and Z840 Host Computer according to [Figure 2-15](#) and [Table 2-15](#).

**Figure 2-15 GSCB Cable**

#	Item	Description
1	GSCB - DB9/F Black	Host Computer DIP Serial Port
2	GSCB - DB9/F Gray	Host Computer RS232
3	GSCB - Green Audio	Host Computer Audio Out (Green)
4	GSCB - Blue Audio	Host Computer Audio In (Blue)
5	GSCB - USB	Not Used
6	GSCB - DB9/M Black	AC Box J56
7	GSCB - DB25/M Black	TGP Gantry Cable

**Table 2-3 GSCB Cable Connection**

- 3.) Remove the old Media Tower from the site and install the new Media Tower (5270510-22) shipped with Upgrade Kit. Connect Media Tower cables to NIO16 Console AC Box and Z840 Host Computer.

**Figure 2-16 Media Tower Connections**

DESCRIPTION	PART NUMBER	QTY
USB Cable (PMT)	5315370	2
PMT Power Cable	5432953-2	1

**Table 2-4 Media Tower Cables**

- 4.) Connect the LCD Monitors for Z840 Host Computer, all signal and power cables are shipped with upgrade kit.
- Scan Monitor  
 Video Cable: from Console Host DP to Monitor DVI  
 Power Cable: from Console AC Box J10 to Monitor

**Figure 2-17 Video Cable and Power Cable for LCD Monitor**


- Image Monitor  
Video Cable: from Console Host DVI to Monitor D-SUB (VGA)  
Power Cable: from Console AC Box J9 to Monitor

**Figure 2-18 Video Cable and Power Cable for LCD Monitor**

g.) Route LCDs video cable and power cable.


**Figure 2-19 Cable Routing and Keeper**

5.) Plug cables into the Switch Hub on Console if request.


- 6.) Connect other cables to Z840 Host Computer according to [Figure 2-13](#).

**Figure 2-20 Switch Hub Connection**

- 7.) Connect the Fiber Optical Cable (5410711) to the Z840 Host Computer, then connect it with the System Fiber Optical Cable.

**Figure 2-21 Fiber Optical Cable Connection**

- 8.) Connect the Power Cable and Ground Cable to the NIO16 Console (See [Figure 2-14](#)).

#### 4.2.5 Console Cover Installation

Please refer to **Service Methods (5366638-2EN\_r19 or later) -> Replacement -> NIO16 Console -> Console Cover Removal and Installation** to install console covers.

---

## Section 5.0

### TGP Board Replacement

- If TGP type is **TGPU** checked in [3.1 Check TGP Board Type](#), follow **Service Methods (5366638-2EN\_r19 or later) -> Replacement -> Gantry -> ORP and TGPU/TGPG** to replace TGP board from TGPU to TGPHD. TGPHD board is included in upgrade kit (**B70112DA**), do not order replacement board.
- If TGP type is **TGPHD** checked in [3.1 Check TGP Board Type](#), skip to [Section 6.0 : Load From Cold Procedure](#).

## Section 6.0:

### Load From Cold Procedure

#### 6.1 LFC Procedure


##### 6.1.1 LFC

Detail LFC procedure refer to **Service Methods (5366638-2EN\_r19 or later) -> Software -> LFC Procedure (13BW16.x/15BW19.x/16BW32.x)**.

Note: Please check and select the correct **Gantry Sub Type** in CT Software Configuration -> Hardware after loading system INFO.


- For Discovery CT590 RT  
**MUST** toggle Product Type - Discovery RT [OFF] and select “LightSpeed\_Discovery\_CT590RT” as Gantry Sub Type.

**Figure 2-22 Hardware Configuration Window (Discovery CT590 RT)**

- For Optima CT580  
**MUST** toggle Product Type - Discovery RT [OFF] and select "LightSpeed\_Optima\_CT580" as Gantry Sub Type.

**Figure 2-23 Hardware Configuration Window (Optima CT580)**

### 6.1.2 Patch Installation

Please confirm that Service Pack/Patch version 1.3 or later is shipped with the upgrade kit. SP1.3 is a mandatory version for all upgrade sites.

## 6.2 Restore System State

**NOTICE** Options need apply new e-licence through e-licence website (<http://elicense.gehealthcare.com/elicense/>) except 10 standard options when OC upgrade to Z840 computer. Please refer to option list recorded in [3.6 Record Option Information Licenses](#).

**Apply new option e-licence in e-licence website, please select “Discovery 590\_Optima 580” for Model selection.**

- 1.) Insert a previously saved System State DVD-RAM Media into the Peripheral Tower DVD-RAM Drive or USB Media in the USB ports
- 2.) Wait until the DVD-RAM Drive or USB Media is ready (i.e., front panel DVD-RAM Drive LED or USB Device Access LED is no longer flashing).
- 3.) Select the [Service] icon to access the CSD (Common Service Desktop).
- 4.) Select: [Utilities].
- 5.) Select: [System State] for DVD-RAM Media or [System State - USB] for USB Media. The System State Save/Restore screen appears.
- 6.) Select [All].
- 7.) Select [Restore]. The Restore System State box appears.
- 8.) Select [Yes].
- 9.) Verify that the "Restore" of System State was successful. If not, restore the System State again. A message at the end of the System State Log window should display: Restore System State Completed Successfully.


**Note:** Please ignore any SwOption Error information during restoring the System State.

- 10.) When completed, select [Cancel].

**Note:** If Dose check option is installed, Dose audit tool starts up when Reconfig Info or Protocol is selected. That shows configuration and Protocol summary. Press Quit to continue.

- 11.) Select [Yes] when the Scan Hardware Reset pop-up appears.
- 12.) When completed, select [Dismiss].
- 13.) Remove System State DVD-RAM disk from the DVD Peripheral Tower or USB Media from USB port and place in a safe location.
- 14.) Reboot the Operators Console by selecting [Shutdown] Desktop, then select [Restart] in the Attention Window.

15.) Verify Restore System State.

- a.) Visually verify that the system has the correct configuration and preference settings by opening a Terminal Window.
- b.) Execute script "reconfig" and manual reset system info.


## 6.3 Convert Protocol

- 1.) Open a Unix Shell, type:  
`{ctuser@hostname} su -`  
Password: **#bigguy**
- 2.) `[root@hostname] /usr/g/scripts/conv_proto <Enter>`
- 3.) Input **y** then press [Enter] to convert protocols.  
*Do you want to convert protocols? [Y/N]*
- 4.) Below information should show if the convert protocol successfully:  
*Success, process finished.*


Note: After finishing the user protocol convert, make sure every protocol can be opened by protocol management and save the current system state to another new media.

## 6.4 Configure System Information on SUSE


- 1.) Configure PPS server information. Please refer to [section 3.2 on page 30](#).
- 2.) Configure all remote network information
  - a.) Select [ImageWorks] from the desktop.
  - b.) Select in sequence: Tools-> Network Configuration

- c.) Select [Add]


d.) Input the Host Name then press enter.

3.) Configure PNF information.

a.) Select the Service Desktop.

b.) Select in sequence: Configuration -> Configure PNF.

4.) Configure AW IP for direct information, Please refer to [section 3.5 on page 36](#).

## Section 7.0: Finalization

### 7.1 Flash Download

The Flash Download takes anywhere from 5-30 minutes, depending on which subsystems require updates.


- 1.) Perform the Flash Download Utility found on the Common Service Desktop - Utilities Tab, select [Flash Download].

**Figure 24 Common Service Desktop - Utilities Tab, Flash Download**

2.) When the Flash Download Window opens, Select [Update].

**Figure 25 Flash Download Window5**


3.) Once the Gantry Hardware Flash Downloads successfully, select [Dismiss].

4.) Close the Common Service Desktop.

## 7.2 GPU Card Checking

Open a shell, then type: `ls /proc/driver/nvidia/gpus | wc -l`

A screenshot of a terminal window titled 'ctuser on bjj02: /usr/g/ctuser'. The terminal has a menu bar with 'File', 'Edit', 'View', 'Terminal', and 'Help'. The command `ls /proc/driver/nvidia/gpus | wc -l` is entered and executed, resulting in the output `2`. A red rectangular box highlights the command and its output. The prompt `[ctuser@bjj02]` is visible before and after the command execution.

- If “2” displays, the Recon GPU card is installed.
- If “1” displays, the Recon GPU card is un-installed.

## 7.3 System Checking

Perform *Quality Assurance Test* of Functional Checks on Service Methods, make sure all of them pass.

## 7.4 Save System State

Refer to [Section 2.0 : Save System State](#) to save the new system state after system function is normal.

## 7.5 Label Attaching

Attach the GPU Fluoro label beside the OC Rating Plate on the rear cover of console.

## 7.6 Old Console Disposal

The old console has components that are needed to cover FRU shortages, so please **RETURN** the de-installed console to GPRS Harvest team. Follow the normal local return process to return the console.

- 1.) Package the old console with the shipping material.
- 2.) Mark type of the old console on the outside of the box.
- 3.) Transporter ships the de-installed console to:
  - **For AM**  
GE Repair Operations Center Attn:  
Asset Recovery Leader, Harvest 120 West Opus Drive Oak Creek, WI 53154
  - **For Europe**  
Stick the harvest return sticker (see below) on the OC box  
Geodis Logistics – Harvest Bat EVL2 Quai 48 – ZI la pièce de la remise Route de Corbeil –  
CD26 – 91090 Lisses – France

<h1>HARVEST Upgrade</h1>	
Upgrade Name	<input type="text"/>
Country	<input type="text"/>
System ID	<input type="text"/>
Order/SO	<input type="text"/>
FSE SSO	<input type="text"/>
<b>Ship to Address</b>	
Geodis Logistics/Harvest	
Bat EVL2 - Porte A1 Quai 9 à 13 - ZI de la remise	
Route de Corbeil CD26 - 91090 Lisses – FRANCE	
<b>Contact : JP Lechat +33 6 72 18 60 32</b>	

- **For ASIA**  
GE Healthcare Japan Corporation Attn:  
Harvest, Service 4-7-127, Asahigaoka, Hino-shi-Tokyo-191-8503 Japan
- **For LATAM**
  - Brazil: Av. Dr. Marcos Pentead de Ulhoa Rodrigues, 690-Galpão 8A Pq.Jubran Cep-06460-040 Barueri-SP Phone: +55 (11) 3206-733
  - Argentina: OCASA Soluciones Logísticas Echeverría 1262/64 C1428DQN – C. A. de Buenos Aires (54 11) 4785-8585

[www.gehealthcare.com](http://www.gehealthcare.com)