# GPU Fluoro Prerequisite Installation Manual

OPERATING DOCUMENTATION

5771958-1EN

Rev 2

@ 2020 General Electric Company.

All rights reserved.

# **LEGAL NOTES**

---

## **TRADEMARKS**

Adobe, the Adobe logo, Acrobat, the Acrobat logo, Exchange, and PostScript are trademarks of Adobe Systems Incorporated or its subsidiaries and may be registered in certain jurisdictions.

Microsoft is a registered trademark and Windows is a trademark of Microsoft Corporation.

All other products and their name brands are trademarks of their respective holders.

## **COPYRIGHTS**

All Material Copyrighted (c) 2020 by the GE Company, All rights reserved.

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## DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent write "Damage In Shipment" on ALL copies of the freight or express bill BEFORE delivery is accepted or "signed for" by a GE representative or hospital receiving agent. Whether noted or concealed, damage MUST be reported to the carrier immediately upon discovery, or in any event, within 14 days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this 14 day period.

To file a report:

- Call 1-800-548-3366 and use option 6.
- Fill out the GIQ workflow for any items missing, damaged, OBF/FOI for in process installs:  
[http://supportcentral.ge.com/ProcessMaps/form\\_new\\_request.asp?prod\\_id=268679&form\\_id=573167&node\\_id=1916016&map\\_id=&reference\\_id=&reference\\_type](http://supportcentral.ge.com/ProcessMaps/form_new_request.asp?prod_id=268679&form_id=573167&node_id=1916016&map_id=&reference_id=&reference_type)
- Contact your local service coordinator for more information on this process.

## CERTIFIED ELECTRICAL CONTRACTOR STATEMENT

All electrical Installations that are preliminary to positioning of the equipment at the site prepared for the equipment shall be performed by licensed electrical contractors. In addition, electrical feeds into the Power Distribution Unit shall be performed by licensed electrical contractors. Other connections between pieces of electrical equipment, calibrations and testing shall be performed by qualified GE Medical personnel. The products involved (and the accompanying electrical installations) are highly sophisticated, and special engineering competence is required. In performing all electrical work on these products, GE will use its own specially trained field engineers. All of GE's electrical work on these products will comply with the requirements of the applicable electrical codes.

The purchaser of GE equipment shall only utilize qualified personnel (i.e., GE's field engineers, personnel of third-party service companies with equivalent training, or licensed electricians) to perform electrical servicing on the equipment.

## IMPORTANT...X-RAY PROTECTION

X-ray equipment if not properly used may cause injury. Accordingly, the instructions herein contained should be thoroughly read and understood by everyone who will use the equipment before you attempt to place this equipment in operation. The General Electric Company, Medical Systems Group, will be glad to assist and cooperate in placing this equipment in use.

Although this apparatus incorporates a high degree of protection against x-radiation other than the useful beam, no practical design of equipment can provide complete protection. Nor can any practical design compel the operator to take adequate precautions to prevent the possibility of any persons carelessly exposing themselves or others to radiation.

It is important that anyone having anything to do with x-radiation be properly trained and fully acquainted with the recommendations of the National Council on Radiation Protection and Measurements as published in NCRP Reports available from NCRP Publications, 7910 Woodmont Avenue, Room 1016, Bethesda, Maryland 20814, and of the International Commission on Radiation Protection, and take adequate steps to protect against injury.

The equipment is sold with the understanding that the General Electric Company, Medical Systems Group, its agents, and representatives have no responsibility for injury or damage which may result from improper use of the equipment.

Various protective materials and devices are available. It is urged that such materials or devices be used.

## LITHIUM BATTERY CAUTIONARY STATEMENTS

### CAUTION

**Risk of Explosion.**

**Danger of explosion if battery is incorrectly replaced.**

**Replace only with the same or equivalent type recommended by the manufacturer. Discard used batteries according to the manufacturer's instructions.**

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## OMISSIONS & ERRORS

Customers, please contact your GE Sales or Service representatives.

GE personnel, please use the Healthcare PQR Process to report all omissions, errors, and defects in this publication.

---

# Revision History

---

Revision	Date	Reason for change
2	05/13/2020	Chapter 2, Section 2.2.2: Add GPU Card Installation for Z8G4 Host Computer Chapter 2, Section 3.0: Add check GPU ECC state
1	05/09/2017	Initial release

*Page Left Intentionally Blank*

# Table of Contents

---

## Preface

**Publication Conventions ..... 17**

**Section 1.0**

**Safety & Hazard Information ..... 17**

1.1 Text and Character Representation..... 17

1.2 Graphical Representation ..... 18

**Section 2.0**

**Publication Conventions ..... 19**

2.1 General Paragraph and Character Styles..... 19

2.2 Page Layout..... 19

2.3 Computer Screen Output/Input Text Character Styles ..... 20

2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)..... 20

## Chapter 1

**Introduction..... 21**

**Section 1.0**

**Flowchart ..... 21**

**Section 2.0**

**Material List ..... 22**

## Chapter 2

**Installation Procedure ..... 23**

**Section 1.0**

**Overview ..... 23**

**Section 2.0**

**Procedure..... 23**

2.1 Preparation ..... 23

2.2 GPU Card Installation ..... 24

2.2.1 GPU Card Installation for OpenOC with Z840 Host Computer ..... 24

2.2.2 GPU Card Installation for OpenOC with Z8G4 Host Computer ..... 31

2.3 Check GPU Card ..... 37

2.3.1 GPU Card Inspection ..... 37

2.4 Power On..... 39

2.5 Install option (if required) ..... 39

**Section 3.0**

**Finalization..... 39**

---

# Preface

## Publication Conventions

---

Purpose: Standardize conventions for representing information is a uniform way of communicating information to a reader in a consistent manner. Conventions are used so that the reader can easily recognize the actions or decisions that must be made. There are a number of character and paragraph styles used in this publication to accomplish this task. Please become familiar with them before proceeding forward.

## Section 1.0

### Safety & Hazard Information

#### 1.1 Text and Character Representation

Within this publication, different paragraph and character styles have been used to indicate potential hazards. Paragraph prefixes, such as hazard, caution, danger and warning, are used to identify important safety information. Text (Hazard) styles are applied to the paragraph contents that is applicable to each specific safety statement. Words describe the type of potential hazard that may be encountered and are placed immediately before the paragraph it modifies. Safety information will normally include:

- Type of potential Hazard
- Nature of potential injury
- Causative condition
- How to avoid or correct the causative condition

#### EXAMPLES OF HAZARD STATEMENTS USED

It's important that you read and understand hazard statements, and not just ignore them.

**DANGER  
SEVERE  
PERSONAL  
INJURY**

**DANGER IS USED WHEN A HAZARD EXISTS WHICH WILL CAUSE SEVERE PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **ELECTROCUTION**
- **CRUSHING**
- **RADIATION**

**WARNING  
SERIOUS  
PERSONAL  
INJURY**

**WARNING IS USED WHEN A HAZARD EXISTS WHICH COULD OR CAN CAUSE SERIOUS PERSONAL INJURY OR DEATH IF INSTRUCTIONS ARE IGNORED. THEY CAN INCLUDE:**

- **POTENTIAL FOR SHOCK**
- **EXPOSED WIRES**
- **FAILURE TO TAG AND LOCKOUT SYSTEM POWER COULD ALLOW FOR UN-COMMANDED MOTION.**

 **CAUTION**  
Minor Personal Injury  
Caution is used when a hazard exists which can or could cause minor injury to self or others if instructions are ignored. They include for example:

- Loss of critical patient data
- Crush or pinch points
- Sharp objects

 **NOTICE**  
Property Damage Only  
Notice is used when a hazard is present that can cause property damage but has absolutely no personal injury risk. They can include:

- Disk drive will be destroyed
- Internal mechanical damage, such as to the x-ray tube
- Coasting the rotor through resonance.

## 1.2 Graphical Representation

Important information will always be preceded by the exclamation point  contained within a triangle, as seen throughout the Preface chapter. In addition to text, several different graphical icons (symbols) may be used to make you aware of specific types of hazards that could possibly cause harm.

<p><b>ELECTRICAL</b></p>  <p><b>LASER</b></p>  <p><b>LASER LIGHT</b></p>	<p><b>HAZARDS MECHANICAL</b></p>  <p><b>HEAT</b></p> 	<p><b>RADIATION</b></p>  <p><b>PINCH</b></p> 
<p><b>AVOID STATIC ELECTRICITY</b></p> 	<p><b>PROCEDURAL TAG AND LOCK OUT</b></p> 	<p><b>WEAR EYE PROTECTION</b></p>  <p><b>EYE PROTECTION</b></p>

# Section 2.0 Publication Conventions

## 2.1 General Paragraph and Character Styles

Prefixes are used to highlight important non-safety related information. Paragraph prefixes (such as purpose, example, comment or note) are used to identify important but non-safety related information. Text styles are also applied to text within each paragraph modified by the specific prefix.

### EXAMPLES OF PREFIXES USED FOR GENERAL INFORMATION


**Purpose:** Introduces and provides meaning as to the information contained within the chapter, section or subsection, Such as used at the beginning this chapter for example.

**Note:** Conveys information that should be considered important to the reader.

**Example:** Used to make the reader aware that the paragraph(s) that follow are examples of information possibly stated previously.

**Comment:** Represents "additional" information that may or may not be relevant to your situation.

## 2.2 Page Layout

The current section and its title are always shown in the footer of the left (even) page.

An exclamation point in a triangle is used to indicate important information to the user.

Paragraphs preceded by **Alphanumeric** characters (e.g. numbers) contain information that must be followed in a **specific order**.

The current chapter and its title are always shown in the footer of the right (odd) page.

Paragraphs preceded by a **symbol** (e.g. bullets) contain information that has **no specific order**.

Headers and footers in this publication are designed to allow you to quickly identify your location. The document's part number and revision number appears in every header on every page. Odd

numbered page footers indicate the current chapter, its title and current page number. Even page footers show the current section and its title, as well current page number.

## 2.3 Computer Screen Output/Input Text Character Styles

Within this publication, mono-spaced character styles (fonts) are used to indicate computer text that's either screen input and output. Mono spaced fonts, such as courier, are used to indicate text direction. When you type at your keyboard, you are generating computer input. Occasionally you will see the math operators "greater than" and "less than" symbols used to indicate the start and finish of variable output. When reading text generated by the computer, you are reading it as computer generated output. In addition to direction, characters are italicized (e.g. *italics*) to indicate information specific to your system or site.

**Example:**  
Fixed Output      This paragraph's fonts represents computer generated screen "fixed" output. It's output is fixed from the sense that it does not vary from application to application. It's the most commonly used style used to indicate filenames, paths and text that does not change from system to system. The character style used is a fixed width such as courier.

**Example:**  
Variable Output      *This paragraph's font represents computer screen output that is "variable". Its used to represent output that varies from application to application or system to system. Variable output is sometimes found placed between greater than and lesser than operators for clarification. For example: <variable\_ouput> or <3.45.120.3>. In both cases, the < and > operators are not part of the actual input.*

**Example:**  
Fixed Input      This paragraph's font represents fixed input. It's computer input that is typed-in via the keyboard. Typed input that does not vary from application to application or system to system. Fixed text the user is required to supply as input. For example: `cd /usr/3p`

**Example:**  
Variable Input      *This paragraph's font represents computer input that can vary from application to application or system to system. With variable text, the user is required to supply system dependent input or information. Variable input sometimes is placed between greater-than and less-than operators. For example: <variable\_input>. In these cases, the (<>) operators would be dropped prior to input. For example: `yocat hosts | grep <3.45.120.3>` would be typed into the computer as `yocat hosts | grep 3.45.120.3` without the greater-than and less-than operators.*

## 2.4 Buttons, Switches and Keyboard Inputs (Hard & Soft Keys)

Different character styles are used to indicate actions requiring the reader to press either a hard or soft button, switch or key. Physical hardware, such as buttons and switches, are called hard keys because they are hard wired or mechanical in nature. A keyboard or on/off switch would be a hard key. Software or computer generated buttons are called soft keys because they are software generated. Software driven menu buttons are an example of such keys. Soft and hard keys are represented differently in this publication.

**Example:**  
Hard Keys      A power switch **ON/OFF** or a keyboard key like **ENTER** is indicated by applying a character style that uses both over and under-lined bold text that is bold. This is a hard key.

**Example:**  
Soft Keys      Whereas the computer MENU button that you would click with your mouse or touch with your hand uses over and under-lined regular text. This is a soft key.

**END OF SECTION**

---


# Chapter 1

## Introduction

---

### Section 1.0

### Flowchart

## Section 2.0 Material List

MARK	NUMBER	QTY	DESCRIPTION/NAME
1	5767682	1	Label for Option on OC Collector
2	5771958-1EN	1	GPU Fluoro Prerequisite Installation Manual
3	8770000-126	1	GPU PCIe Nvidia M5000

**Table 1-1 GPU FLUORO for Z840 requisite (5772853)**

MARK	NUMBER	QTY	DESCRIPTION/NAME
1	5848913	1	Label for Option on OC Collector
2	5771958-1EN	1	GPU Fluoro Prerequisite Installation Manual
3	8780000-113	1	GPU PCIe Nvidia RTX5000

**Table 1-2 GPU FLUORO for Z8G4 requisite (5845300)**

---

# Chapter 2

## Installation Procedure

---

### Section 1.0 Overview

This procedure describes the GPU Fluoro Prerequisite Installation. This hardware works as a prerequisite with the following options.

- Enhanced Recon Option
- Smartview Fluoro Option

### Section 2.0 Procedure

#### 2.1 Preparation

- 1.) Select one of the following methods to Power OFF the Operator Console if the Power is ON.
  - a.) If Applications are up, click on the [Shut Down] button and select [Shutdown].
  - b.) If Applications are down, open a Unix Shell from pop-up menu. Enter (type) the following command, at the indicated prompts.

```
ctuser@hostname: su -  
Password: #bigguy  
root@hostname: halt
```

- 2.) Power OFF the console at the front panel switch when 'System Halted' message is displayed.

## 2.2 GPU Card Installation

### 2.2.1 GPU Card Installation for OpenOC with Z840 Host Computer

Follow ESD Management and Device Handling procedure.

**Figure 3-1 Host Computer Side Cover**

- 1.) Remove the side access panel by unscrewing its three screws shown as [Figure 3-2](#).

**Figure 3-2 Side Access Panel**

- 2.) Pull out the latch to release the Host PC left side cover..
- 3.) Remove the Expansion Card Support and Airflow Guide as shown below:
  - i.) Pull the handle of right side of the Expansion Card Support.
  - ii.) Hold two handles of Airflow Guide, Pull out vertically.

**Expansion Card Support**

**Airflow Guide**


**Figure 3-3 Remove the Protect Covers**



**Figure 3-4 Z840 Host Computer**

- 4.) Open PCI card retention clamp as follows:
  - a.) Press lever of the retention clamp to unhook.
  - b.) Rotate the retention clamp, and open it shown as [Figure 3-5](#).


**Figure 3-5 Open PCI Card Retention clamp**

c.) Remove two PCI slot cover from back side shown as [Figure 3-6](#)


**Figure 3-6 PCI Slot Cover Removal**

5.) Install the NVIDIA Quadro M5000 GPU Card to the slot #6 as follows:


a.) Remove the power connector from the bracket to free the second slot as follows:

**Figure 3-7 Remove Power Connector**

b.) Carefully push the GPU card in the slot#6 and avoid cables interference.

The Second Slot

Figure 3-8 Insert GPU Card in the Slot


Figure 3-9 Card Slot Release Lever

- c.) Connect the power cable disconnected in Step a.) to the GPU card.

Power Cable Connector

**Figure 3-10 Connect Power Cable to GPU Card**

- 6.) Close the PCI card retention clamp.

Note: If the PCI card retention clamp levers do not close, make sure that all cards are properly seated, and then try again.

- 7.) Reinstall the Expansion Card Support and Airflow Guide shown as [Figure 3-3](#).  
8.) Reinstall the computer left side covers and the side access panel.

## 2.2.2 GPU Card Installation for OpenOC with Z8G4 Host Computer

Follow ESD Management and Device Handling procedure.

Figure 4-11 Host Computer Side Cover

1.) Remove the side access panel by unscrewing its three screws shown as [Figure 4-12](#).


Figure 4-12 Side Access Panel

- 2.) Pull out the latch to release the host computer left side cover. shown as [Figure 4-13](#).



**Figure 4-13 Side Access Panel Removal**

- 3.) Open and remove the module small side cover from left side access panel.

**Figure 4-14 Remove module small side cover**

- 4.) Remove the fourth and fifth PCI slot cover from back side shown as [Figure 4-15](#)

**Figure 4-15 PCI Slot Cover Removal**

5.) Install the NVIDIA Quadro RTX5000 GPU Card to the slot #4 as follows:

a.) Remove two power connectors from the bracket, they will be used to power the GPU card.


Figure 4-16 Remove Power Connector

b.) Carefully push the GPU card in the slot#4 and avoid cables interference.


Figure 4-17 Insert GPU Card in the Slot



Figure 4-18 Card Slot Release Lever

c.) Connect the power cable disconnected in Step a.) to the GPU card.


**Figure 4-19 Connect Power Cable to GPU Card**

- 6.) Reinstall all removed components and reconnect any cables that have been disconnected.
- 7.) Reinstall the module small side cover.
- 8.) Reinstall the host computer left side cover.

## 2.3 Check GPU Card

### 2.3.1 GPU Card Inspection

- 1.) Confirm Host Computer powers up when console power is turned on.
- 2.) Check the GPU card installed. Open a shell, then type:  
{ctuser@hostname} `ls /proc/driver/nvidia/gpus | wc -l`. Check the result as below:

```
ctuser on bj02: /usr/g/ctuser
File Edit View Terminal Help
{ctuser@bj02}ls /proc/driver/nvidia/gpus | wc -l
2
{ctuser@bj02}
```

**Figure 4-20 GPU Check Result**

- If “2” displays, the Recon GPU card is installed.
  - If “1” displays, the Recon GPU card is un-installed.
- 3.) Ensure GPU ECC state is ON:
    - a.) Open a Unix shell and log on as root:
    - b.) Type: `su -` [ENTER]
    - c.) Type the root password [ENTER]
    - d.) Type: `nvidia-smi` [ENTER]

e.) Check GPU ECC status as below:

- \* If the GPU ECC is **ON**, below is what the output would look like (boxed in green):  
 (The screenshot below is just an example)

```

+-----+
| NVIDIA-SMI 361.28      Driver Version: 361.28      |
+-----+-----+
| GPU  Name          Persistence-M| Bus-Id        Disp.A | Volatile Uncorr. ECC |
| Fan  Temp  Perf    Pwr:Usage/Cap|      Memory-Usage | GPU-Util  Compute M. |
+-----+-----+
|   0   Quadro M5000        Off   | 0000:03:00.0  Off   |          0%         0 |
| 38%   32C   P8     12W / 150W | 109MiB / 7679MiB |          0%         Default |
+-----+-----+
|   1   Quadro K620        Off   | 0000:04:00.0  On    |          0%         N/A |
| 34%   35C   P8     1W / 30W  | 142MiB / 2047MiB |          0%         Default |
+-----+-----+

+-----+
| Processes:                      GPU Memory |
| GPU      PID  Type  Process name      Usage |
+-----+-----+
|   0      6177  G    /usr/bin/Xorg      9MiB |
|   0     13469  C    iungo_ig11        78MiB |
|   1      6177  G    /usr/bin/Xorg     133MiB |
+-----+
    
```

- \* If the GPU ECC is **OFF**, below is what the output would look like (boxed in red):  
 (The screenshot below is just an example)

```

+-----+
| NVIDIA-SMI 361.28      Driver Version: 361.28      |
+-----+-----+
| GPU  Name          Persistence-M| Bus-Id        Disp.A | Volatile Uncorr. ECC |
| Fan  Temp  Perf    Pwr:Usage/Cap|      Memory-Usage | GPU-Util  Compute M. |
+-----+-----+
|   0   Quadro M5000        Off   | 0000:03:00.0  Off   |          0%         Off |
| 38%   32C   P8     12W / 150W | 109MiB / 8191MiB |          0%         Default |
+-----+-----+
|   1   Quadro K620        Off   | 0000:04:00.0  On    |          0%         N/A |
| 34%   35C   P8     1W / 30W  | 142MiB / 2047MiB |          0%         Default |
+-----+-----+

+-----+
| Processes:                      GPU Memory |
| GPU      PID  Type  Process name      Usage |
+-----+-----+
|   0      6132  G    /usr/bin/Xorg      9MiB |
|   0     12267  C    iungo_ig11        78MiB |
|   1      6132  G    /usr/bin/Xorg     133MiB |
+-----+
    
```

f.) How to turn ECC back on

- \* Type: `nvidia-smi -g 0 --ecc-config=1` [ENTER]
- \* A message will show that ECC is enabled and a reboot is required:

```

Enabled ECC support for GPU 0000:03:00.0.
All done.
Reboot required.
    
```

- \* After reboot, check that the ECC is **ON** according to previous steps

## 2.4 Power On

Power On the system and confirm both of LCD display show normal screens, and start up without error.

## 2.5 Install option (if required)

Install the option which requires GPU Fluoro Option.

## Section 3.0 Finalization

- 1.) Execute sanity scan and confirm the scan completed and image is generated without error.
- 2.) Execute GPU Diagnostics in Console Diag. and confirm it's completed successfully.
- 3.) Attach the name plate near the rating plates on the right side cover of console. See [Figure 4-21](#)


Figure 4-21 Attach GPU Rating Plate

[www.gehealthcare.com](http://www.gehealthcare.com)