



GE Medical Systems

Technical Publications

Direction 2304882

Revision 2

Signa[®] TwinSpeed Block Diagrams & Supplemental Schematics

**Copyright© 2001, 2002, 2003 by General Electric Company, Inc.
All Rights Reserved**

Operating Documentation

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

Call Traffic and Transportation, Milwaukee, WI (414) 785-5052/8*323-5052 **immediately** after damage is found. At this time be ready to supply name of carrier, delivery date, consignee name, freight or express bill number, item damaged and extent of damage.

Complete instructions regarding claim procedure are found in Section "S" of the Policy & Procedure Bulletins.

3/12/92



GE Medical Systems

*GE Medical Systems: Telex 3797371
P.O. Box 414, Milwaukee, Wisconsin 53201 U.S.A.
(Asia, Pacific, Latin America, North America)*

*GE Medical Systems — Europe: Telex 261794
Shortlands, Hammersmith, London W6 8BX U.K.*

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

AVERTISSEMENT

- CE MANUEL DE MAINTENANCE N'EST DISPONIBLE QU'EN ANGLAIS.
- SI LE TECHNICIEN DU CLIENT A BESOIN DE CE MANUEL DANS UNE AUTRE LANGUE QUE L'ANGLAIS, C'EST AU CLIENT QU'IL INCOMBE DE LE FAIRE TRADUIRE.
- NE PAS TENTER D'INTERVENTION SUR LES ÉQUIPEMENTS TANT QUE LE MANUEL SERVICE N'A PAS ÉTÉ CONSULTÉ ET COMPRIS.
- LE NON-RESPECT DE CET AVERTISSEMENT PEUT ENTRAÎNER CHEZ LE TECHNICIEN, L'OPÉRATEUR OU LE PATIENT DES BLESSURES DUES À DES DANGERS ÉLECTRIQUES, MÉCANIQUES OU AUTRES.

WARNUNG

- DIESES KUNDENDIENST-HANDBUCH EXISTIERT NUR IN ENGLISCHER SPRACHE.
- FALLS EIN FREMDER KUNDENDIENST EINE ANDERE SPRACHE BENÖTIGT, IST ES AUFGABE DES KUNDEN FÜR EINE ENTSPRECHENDE ÜBERSETZUNG ZU SORGEN.
- VERSUCHEN SIE NICHT, DAS GERÄT ZU REPARIEREN, BEVOR DIESES KUNDENDIENST-HANDBUCH NICHT ZU RATE GEZOGEN UND VERSTANDEN WURDE.
- WIRD DIESE WARNUNG NICHT BEACHTET, SO KANN ES ZU VERLETZUNGEN DES KUNDENDIENSTTECHNIKERS, DES BEDIENERS ODER DES PATIENTEN DURCH ELEKTRISCHE SCHLÄGE, MECHANISCHE ODER SONSTIGE GEFAHREN KOMMEN.

AVISO

- ESTE MANUAL DE SERVICIO SÓLO EXISTE EN INGLÉS.
- SI ALGÚN PROVEEDOR DE SERVICIOS AJENO A GEMS SOLICITA UN IDIOMA QUE NO SEA EL INGLÉS, ES RESPONSABILIDAD DEL CLIENTE OFRECER UN SERVICIO DE TRADUCCIÓN.
- NO SE DEBERÁ DAR SERVICIO TÉCNICO AL EQUIPO, SIN HABER CONSULTADO Y COMPRENDIDO ESTE MANUAL DE SERVICIO.
- LA NO OBSERVANCIA DEL PRESENTE AVISO PUEDE DAR LUGAR A QUE EL PROVEEDOR DE SERVICIOS, EL OPERADOR O EL PACIENTE SUFRAN LESIONES PROVOCADAS POR CAUSAS ELÉCTRICAS, MECÁNICAS O DE OTRA NATURALEZA.

ATENÇÃO

- ESTE MANUAL DE ASSISTÊNCIA TÉCNICA SÓ SE ENCONTRA DISPONÍVEL EM INGLÊS.
- SE QUALQUER OUTRO SERVIÇO DE ASSISTÊNCIA TÉCNICA, QUE NÃO A GEMS, SOLICITAR ESTES MANUAIS NOUTRO IDIOMA, É DA RESPONSABILIDADE DO CLIENTE FORNECER OS SERVIÇOS DE TRADUÇÃO.
- NÃO TENDE REPARAR O EQUIPAMENTO SEM TER CONSULTADO E COMPREENDIDO ESTE MANUAL DE ASSISTÊNCIA TÉCNICA.
- O NÃO CUMPRIMENTO DESTE AVISO PODE POR EM PERIGO A SEGURANÇA DO TÉCNICO, OPERADOR OU PACIENTE DEVIDO A' CHOQUES ELÉTRICOS, MECÂNICOS OU OUTROS.

AVVERTENZA

- IL PRESENTE MANUALE DI MANUTENZIONE È DISPONIBILE SOLTANTO IN INGLESE.
- SE UN ADDETTO ALLA MANUTENZIONE ESTERNO ALLA GEMS RICHIEDE IL MANUALE IN UNA LINGUA DIVERSA, IL CLIENTE È TENUTO A PROVVEDERE DIRETTAMENTE ALLA TRADUZIONE.
- SI PROCEDA ALLA MANUTENZIONE DELL'APPARECCHIATURA SOLO DOPO AVER CONSULTATO IL PRESENTE MANUALE ED AVERNE COMPRESO IL CONTENUTO.
- NON TENERE CONTO DELLA PRESENTE AVVERTENZA POTREBBE FAR COMPIERE OPERAZIONI DA CUI DERIVINO LESIONI ALL'ADDETTO ALLA MANUTENZIONE, ALL'UTILIZZATORE ED AL PAZIENTE PER FOLGORAZIONE ELETTRICA, PER URTI MECCANICI OD ALTRI RISCHI.

警告

- ・このサービスマニュアルには英語版しかありません。
- ・GEMS以外でサービスを担当される業者が英語以外の言語を要求される場合、翻訳作業はその業者の責任で行うものとさせていただきます。
- ・このサービスマニュアルを熟読し理解せずに、装置のサービスを行わないで下さい。
- ・この警告に従わない場合、サービスを担当される方、操作員あるいは患者さんが、感電や機械的又はその他の危険により負傷する可能性があります。

注意:

- 本维修手册仅存有英文本。
- 非 GEMS 公司的维修员要求非英文本的维修手册时，客户需自行负责翻译。
- 未详细阅读和完全了解本手册之前，不得进行维修。
- 忽略本注意事项会对维修员，操作员或病人造成触电，机械伤害或其他伤害。

REVISION HISTORY

REV	DATE	PRIMARY REASON FOR CHANGE
0	August 24, 2001	Initial release.
1	October 31, 2001	Second release: 1. SYSTEM Tab Page 1-3 and 1-4; Many circuit corrections, is redrawn. 2. GRADIENT Tab Page 1-2; Many circuit corrections, is redrawn. Corrected header mistake also.
2	March 24, 2003	Third release, various corrections to match current product, specifically: 1. SYSTEM Tab Page 1-3 and 1-4; Removed the Vacuum Sensor and connections. 2. Many changes to Tab TPS/ISE (1.5T).

LIST OF EFFECTIVE PAGES

<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Title Page	2	2-1	0	1-1 to 1-2	0	PDU		1-3 to 1-4	0
GE Logo Page	-	3-1 to 3-6	0			Tab MR406	-	1-5 to 1-14	2
2128126	0*	SYSTEM		PATIENT HANDLING		i	0	2-1 to 2-6	2
A	2	Tab MR364	-	Tab MR405	-				
i	0	1-1 to 1-2	0	1-1 to 1-4	0				
		1-3 to 1-4	2	GRADIENT		TPS/ISE (1.5T)		COOLING SYSTEM	
OVERVIEW		OPERATOR		Tab MR418	-	Tab 2163921	-	Tab 2308173	-
Tab MR367	-	WORKSPACE		1-1	0	i	2	i	0
i	0	Tab 2158195	-	1-2	1	1-1 to 1-2	2	1-1	0
1-1	0								

* This revision number/letter corresponds to the indicated document's revision control system.

CONTENTS

OVERVIEW

Purpose
Standard Conventions
Abbreviations



SYSTEM

Overall System Block Diagrams



OPERATOR WORKSPACE

Operator Workspace Block Diagrams



PATIENT HANDLING

Patient Handling Block Diagrams



GRADIENT

Gradient Block Diagram (ACGD)



PDU

PDU Block Diagrams

Note:

Phoenix PDU found in Vendor Manuals only. (SEE Service Methods CDROM)



TPS/ISE (1.5T)

TPS/ISE Block Diagrams (1.5T)
RF Supplemental Schematics (1.5T)



COOLING SYSTEM

Water Chillers TGWC & TSCC Block Diagrams

Note:

For all Chillers and Cooling System variations (SEE Vendor Manuals on MR Service Methods CDROM)



OVERVIEW

TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
SECTION 1 – PURPOSE	1-1
SECTION 2 – STANDARD CONVENTIONS	2-1
SECTION 3 – ABBREVIATIONS	3-1

SECTION 1 – PURPOSE


Note

This revision of Direction 2304882 is compatible with 1.5T TWINSPEED hardware.

The purpose of this document is to provide an overview of the hardware which makes up the MR Signa System, and a quick means to troubleshoot major signals within the system.

It is recommended that you page through this manual to become familiar with its content and organization, beginning with System Tab, OVERALL SYSTEM. This section shows how the major subsystems are grouped by cabinet, and the main communication lines between them.

The remaining block diagrams are organized into sections which correspond to major subsystems. Some diagrams show hardware from another subsystem in order to complete an entire signal flow route. For such exceptions, the board/module title of the circuitry shown will generally help identify the applicable subsystem. See Illustration 1–1 for tab organization of this manual.




TAB ORGANIZATION FOR THIS DIRECTION
ILLUSTRATION 1

These block diagrams, along with provided diagnostics, are intended to be the primary tools for troubleshooting the MR System from a system or subsystem level on down to the specific module or board level. Since some of the RF subsystem block diagrams are presently insufficient for troubleshooting to the FRU (Field Replaceable Unit) level, schematics have been provided to supplement those block diagrams. When using the diagnostics, the block diagrams can provide a better understanding of what circuitry the diagnostics are testing, and how to further isolate a problem if required.


SECTION 2 – STANDARD CONVENTIONS

The following standard conventions are used within the block diagrams:

- The type of line used to enclose a major box on a block diagram indicates the following:



- The component designator for cabinets, modules, circuit boards, etc. is enclosed in “()” after or below their name on the block diagram.
- Logic levels: A “0” indicates a logic “low” level (typ. <0.8VDC), a “1” indicates a logic “high” level (typ. > 2.0VDC).
- An “*” is used to identify an “active low” signal (e.g. “INT*” signal when at a logic “0” indicates that an interrupt is present; a logic “1” level indicates that no interrupt is present).
- A “#” indicates the number of wires in a cable (e.g. $\frac{9}{\text{---}}$ indicates 9 wires).
- Signal flow is generally from left to right on a block diagram. Arrows are used to identify exceptions to this rule and to aid in identifying input/output signal lines on the diagram.
- The “x” symbol (e.g. $\text{---}x$) shown on boxes in a block diagram are used to indicate that the function within the box is enabled by the “active” signal on this input line.



For this example, when signal TEMPEXTREME goes active (“1”), relay K42 will close (become enabled) and the signal SETIND now becomes active to set a temperature indicator.

SECTION 3 – ABBREVIATIONS

Note

This list has been compiled from a wide variety of sources. Not all of these abbreviations will apply to your system. For an electronic copy of over 1000 acronyms, see *Direction 2124201-2*.

A

AC = Alternating Current
 ACK = Acknowledge
 ACQ = Acquisition
 A/D = Analog-to-Digital
 ADC = Analog to Digital Converter
 ADD or ADDR = Address
 ADJ = Adjust
 ADV = Advance
 AGC = Automatic Gain Control
 AI = Applications Interface, Artificial Intelligence
 ALGN = Align
 ALU = Arithmetic Logic Unit
 AM = Amplitude Modulation
 AMP = ampere, amplifier
 A/N = Alpha-numeric
 ANSI = American National Standards Institute
 AP = Array Processor
 APM = Analog Power Monitor
 APS = Auto Pre-Scan
 APSGT = Auto Pre-Scan Gain Control
 AQT = Auto QTUNE
 ARCH = Archive/Remove Process
 ASC = Automated Service Center
 ASCII = American Standard Code for Information Interchange
 ASM or ASSY = Assembly or Analog Service Module
 ASM = Analog Service Module
 ATTEN = Attenuator
 AUTO = Automatic
 AUX = Auxiliary

B

B = magnetic field
 BATT = Battery
 BB = Buffered Data or Broadband
 BBT = BroadBand Tool
 BC = Body Coil

BD = Board
 BDS = Buffered Data Strobe
 BLK = Black
 BNC = BNC Connector
 BRT = Brightness
 BUF = Buffer
 B/W = BandWidth

C

C = centigrade, Celsius
 CAB = Cabinet
 CAL = Calibration
 CALC = Calculate
 CATTEN = Course Attenuation
 CB = Circuit Breaker
 CCW = Counter Clockwise
 CD-ROM = Compact Disk Read-Only Memory
 CE = Chip Enable
 CERD = Combined Exciter, Receiver, DAB
 CF = Center Frequency
 CFA = Center Frequency Adjust
 CHAN or CHNL = Channel
 CHK = Check
 CKT = Circuit
 CLR = Clear
 CLI = Command Line Interpreter
 CMD = Command
 CMMR = Common Mode Rejection Ratio
 CMOS = Complimentary Symmetry Metal Oxide Semiconductor
 CMP = Compare or Compressor
 CM/PM = Communications Manager / Power Monitor Board
 CNT = Contrast
 COAX = Coaxial
 COF = Cut Off Frequency
 CON or CONS = Console
 CONFIG = Configuration
 CONT = Continuous

CONV = Converter
 CPD = Communication Pin Driver
 CPU = Central Processing Unit
 CRT = Cathode-Ray Tube
 CT = Computed Tomography
 CTL = Clock Time Control
 CTRL = Control
 CV = Configuration Variable
 CW = Clockwise
 CYL = Cylinder

D

D/A = Digital-to-Analog
 DAC = Digital to Analog Converter
 dB = Decibel
 dBm = Decibels above or below one Milliwatt
 DC = Direct Current
 DCI = Digital Control Interface
 DCKSTP = Dock Stop
 DD = Dynamic Disable
 DDS = Direct Digital Synthesizer
 DDSB = Dynamic Disable Switch Box
 DEC = Decimal
 DECC = Digital Eddy Current Compensation
 DEG = Degree(s)
 DEMUX = Demultiplexer
 DEV = device
 DIA = Diagnostic Imaging Accessories
 DIAG = Diagnostic
 DIFF = Differential
 DIG = Digit
 Dig B0 = Digital B0 Adjustment
 DIP = Dual In-line Package
 DIR = Direction
 DIS = Disable
 DMA = Direct Memory Access
 DMAC = Direct Memory Access Controller
 DMM = Digital Multi-Meter
 DPR or DPRAM = Dual Port Random Access Memory
 DQA = Daily Quality Assurance
 DRAM = Dynamic Random Access Memory
 DRC = Dynamic Ram Controller
 DRVR = Driver
 DSP = Digital Signa Processor
 DTS = Display touch Screen
 DVM = Digital Volt Meter

DWG = Drawing

E

ECC = Extended Cursor Control or Error Checking and Correction
 ECG = Electro-Cardiogram
 EDM = Equipment Diagnostic Monitor
 EEPROM = Electrically Erasable Programmable Read Only Memory
 EFB = Envelop Feedback
 EPI = Echo Planer Imaging
 e.g. = for example
 EM or EMERG = Emergency
 EMI = Electro-Magnetic Interference
 EN or ENA = Enable
 ENV = Envelope
 EOS = End Of Sequence
 EOT = End Of Travel
 EOW = End Of Waveform
 EPO = Emergency Power Off
 EPROM = Erasable Programmable Read Only Memory
 ERU = Emergency Rundown Unit
 EXER = Exerciser

F

f = farad
 F = Fahrenheit, Fuse (designator)
 FBD = Functional Block Diagram
 FBK = Feedback
 FC = footcandle
 Fc = Center Frequency
 FDBK = Feedback
 FDO = Future Delivery Order
 FE = Field Engineer
 FET = Field-Effect Transistor
 FF or FIFO = First In – First Out
 FFD = Full Field Distortion
 FFT = Fast Fourier Transform
 FID = Free Induction Decay
 FIDAT = FIFO In Data
 FIFO or FF = First In – First Out
 FILT = Filter
 FM = Frequency Modulated
 FMI = Field Modification Instruction
 FOV = Field Of View
 FPU = Floating Point Unit
 FREQ = Frequency

FRESBECC = Frequency Shift B0 Eddy Current Compensation
FRNT = Front
FRU = Field Replaceable Unit
ft = foot
FT = Feet
FTP = File Transfer Protocol

G

G = gauss
G/A = Gradient Amplifier
GAP = Gradient Amplifier Processor
GASM = Gram Analog Service Module
G(ASM) = GASM and ASM
GEN = Generator
GEPI = Gradient Echo Planer Imaging
GIF = Gradient Interface
GIP = Gradient Interface Processor
GND = Ground
GPIB = General Purpose Interface Board
GPM = Gradient Power Module (aka 8651)
GRAD = Gradient
GRAM = Gradient Ramp Accelerator Module
GRN = Green
GRY = Grey
GP = Gradient Processor

H

H = henry
H = Hexadecimal
HC = Head Coil
HD = Head
He = Helium
He-ALARM = Helium Alarm
HEX = Hexadecimal
He-WARNING = Helium Warning
HLFC = Host Load From Cold
HORIZ or HORZ = Horizontal
HPC = Hardware and Protocol Control
HPFS = Host Platform Software
hr = hour
HSN = High Speed Network

HSS (D) = High Speed Serial Data Link
HTR = Heater
HV = High Voltage
HV PSU = High Voltage Power Supply Unit
HVAC = Heating, Ventilation, & Air Conditioning
Hz = hertz

I

I = Current or Phase or Real
IAMP = Amplifier Current (Gradient)
IC = Independent Console, Integrated Circuit
ICD = Installation Certification Document
ICL = Coil Overcurrent (Gradient)
ID = inner diameter
ID = Identity, Identifier
i.e. = that is to say
IEC = International Electrical Code
IEEE = Institute of Electrical & Electronics Engineers
I/F = Interface
IGRAD = GRAM HSSD D/A converter voltage output
in. = inch
INFO = Information
INS, INST, or INSTR = Instruction
INT or INTR = Interrupt
I/O = Input/Output
IP = Image Processor or Plate Current (RF amp) or Instruction Processor or Internet Protocol
IPA = Intermediate Power Amplifier
IPC = Interprocess Communication
IPG = Integrated Pulse Generator
IQ = Image Quality
IRQ = Interrupt Request
ISA = Interrupt Start Address
ISE = Integrated Systems Electronics
ISI = Inter-Sequence Interrupt
ISO = Isolator or Isolation
IVI = Interactive Vascular Imaging
IVME = Internal VME

J

J = Jack Connector (designator)
JP = Jumper

K

K = Kelvin (unit of temperature) or Kilo (Thousand)
 KB = one thousand bytes of memory
 KCal = kiloCalorie
 kg = kilogram
 kHz = kilohertz
 KVA = kilo Volt Amperes
 kW = Kilowatt

L

L or LFT = Left
 LAN = Local Area Network
 lb = pound
 LBK = Loopback
 LCA = Logic Cell Arrays
 LED = Light Emitting Diode
 LFC = Load From Cold
 LHe =Liquid Helium
 LIN = Linear
 LM or LMK = Landmark
 LMP = Lamp
 LNK = Link
 LO = Local Oscillator
 LOEP = List of Effective Pages
 LONG = Longitudinal
 LVLE = Low Voltage, Low Energy

M

m = meter
 M = mega – one million
 mA = millamps
 MAG = Magnet
 MB = megabyte
 MCD = Multi-Coil switch Driver
 MDS = Multidrop Serial interface
 MEM = Memory
 Mfg. = Manufacturing
 MFU = Multi-Format Unit
 MGMNT = Management
 mGRAM = mini Gradient Ramp Accelerator Module
 MHz = megahertz

MIF = Master Interface
 mm = millimeter
 MOD = Modifier
 MON = Monitor
 MOSFET = Metal Oxide Semiconductor Field Effect Transistor
 MOV = Metal Oxide Varistor
 MPS = Manual Pre-Scan
 MR = Magnetic Resonance
 MSG = Message
 MTR = Motor
 MUX = Multiplexer
 mV = millivolts
 MW = megawatt
 mW = milliwatt

N

N = Negative
 N2 = Nitrogen
 n.c. = normally closed
 NEC = National Electrical Code
 Ni-Cad = Nickel Cadmium
 n.o. = normally open
 NO. = Number
 NOPROC = No Processing
 NOREC = No Reconstruction
 nS = Nanosecond

O

O2 = oxygen
 OC = Operator's Console
 OD = outside diameter
 OPI = Oblique Plane Imaging
 OPTO = Optical
 ORG = Orange
 OT = Overtemp
 OV = Over voltage
 OVR = Over
 OVRD or OL = Overload
 oz = ounce
 OW = Operator Workspace

P

P = Positive
 PA = Patient or Power Amplifier
 PAC = Physiological Acquisition Controller
 PAL = Programmable Array Logic or Patient Alignment Lights
 PAN = Primary Archive Node
 PAT = Patient
 PCI = Programmable Communications Interface
 PCT = Pulse Controller Task
 PDI = Product Delivery Instruction
 PDU = Power Distribution Unit
 PEN = Penetration
 PIN = Positive–Intrinsic–Negative
 PIO = Programmed Input/Output
 PIXEL = Picture Element
 PM = Phase Modulation, Periodic Maintenance
 P/N = Positive/Negative
 PNL = Panel
 POS = Position
 POT = Potentiometer
 P–P = Peak to Peak
 ppm = parts per million
 PREAMP = Preamplifier
 PROC = Process, Processor
 PROG = Program, –ed, –able
 PROM = Programmable Read Only Memory or Program
 PS = Power Supply
 PSD = Pulse Sequence Database
 PSG = Pulse Sequence Generator
 PSU = Power Supply Unit
 PSS = Power Supply System Board
 PUR = Purple
 PVC = Polyvinylchloride
 PW = Pulse Width
 PWA = Printed Wiring Assembly
 PWB = printed wiring board
 PWR = Power

Q

Q = Quadrature, Quality Factor or Imaginary
 Qty = Quantity
 QUAD = Quadrature or Quadrant

R

R or RT = Right

RAM = Random Access Memory
 RCV = Receive
 RD = Read
 RDY = Ready
 RECON = Reconstruction
 REF = Reference or Reflected
 REG = Register or Regulator
 RET or RTN = Return
 REV = Revision or Reverse
 RF = Radio Frequency
 RFAMP = Radio Frequency Amplifier
 RFI = Radio Frequency Interference
 RF/PEN = Combined RF/Penetration Cabinet
 RFSC = RF System Controller
 RGB = Red Green Blue
 RLA or RLY = Relay
 rms = root mean square
 ROI = Region of Interest
 ROM = Read Only Memory
 RSR = Register
 RST = Reset
 R/W = Read/Write
 RX = Receive
 RXD = Receive Data

S

S = Switch (designator)
 SACK = Slave Acknowledge
 SAG = Sagittal
 SAR = Specific Absorption Rate
 SARR = Stand–Alone Retrospective Reconstruction
 S/C = Superconducting
 SCA = Scan Control Assembly
 SCB = Signal Collection Board
 SCR = Silicon Controlled Rectifier
 SCSI = Small Computer System Interface
 SEC = second
 SEG = Segment
 SEL = Select
 SEPI = Spin Echo Planer Imaging
 SGI = Silicon Graphics Inc.
 SIG = Signal
 SII = "S–2" (type of magnet)
 SIII = "S–3" (type of magnet)
 SIMMS = Single In–line Memory Module
 SIP = Single In–line Package

REV 0

DIRECTION 2304882

SL = Slice
 SLE = System Level Exerciser
 SMPTE = Society of Motion Picture and Television Engineers
 SN = Service Note
 SNR = Signal-to-Noise Ratio
 SPI = Serial Peripheral Interface
 SPT = System Performance Test
 SPU = Signal Processing Unit
 SRAM = Static Random Access Memory
 SRF = Solid State RF Amplifier SGA
 SRI = Scan Room Interface
 SRV = Service
 SSM = System Support Module
 SSP = Sequence Synchronous Protocol
 SST = Small Sample Test
 STAB = Stability
 STAT = Status
 STD = Standard
 STD DEV = Standard Deviation
 SUB-D = Subminiature D (connector)
 SW = Switch or Software
 S/W = Software
 SYNC = Synchronous
 SYS = System
 SYSID = System Identification

T

T = tesla
 TBD = To Be Determined
 TBL = Table
 TEMP = Temperature
 TFU = Thyristor Firing Unit
 TG = Transmit Gain
 THD = TPS Host Diagnostic
 TNF = Transient Noise Filter
 TP = Test Point (designator)
 TPS = Transceiver Processing and Storage
 TR or T/R = Transmit-Receive
 TRMAP = Transmit/Receive Map
 TRM = Twin Coil Resonator Module
 TS = Terminal Strip
 TST = Test
 TTL = Transistor to Transistor Logic
 TX = Transmit
 TXD = Transmit Data

U

UART = Universal Asynchronous Receiver-Transmitter

uf = microfarad
 UFI = Ultra Fast Imaging
 UNBLK = Unblank
 uP = Microprocessor
 UPLMT = Up Limit
 USART = Universal Synchronous/Asynchronous Receiver-Transmitter
 uS or usec = microsecond
 UV = Under voltage

V

VAC = Volts Alternating Current
 VAR = Variable
 VBUS = Bus Voltage
 VCR = Video Cassette Recorder
 VCNTL = Control Voltage for 8645 Techrons
 VDC = Volts Direct Current
 VERT = Vertical
 VFILT = GRAM Output Voltage
 VIN = DAC Idi/dt to command the GRAM PWM voltage
 VME = Bus Standard Veras Module Europe
 VMIF = MIF Control Voltage
 Vout-H/L = Power Module High/Low Voltage Output
 Vp-p = Volts Peak-to-Peak
 VRMS = Volts, Root Mean Square
 VSA = VSB Address
 VSB = VME Subsystem Bus
 VTAC = Vacuum Tube Amplifier Cavity
 VTR = Video Tape Recorder

W


WARP = Waveform And Rotation Processor
 WR = Write
 WHT = White
 WND = Window
 WORM = Write Once Read Many
 WRT = Write

X

XCVR = Transceiver
 XFMR = Transformer
 XMT or XMIT = Transmit


Y

YEL = Yellow
 YMS = Yokogawa Medical Systems
 YP = Yellow Pages (Sun Computer)





OVERALL SYSTEM – OCTANE
(OPERATOR WORKSPACE CABINET/TABLE – OCTANE)

OPTIONAL CABINET COMMUNICATION LINES INDICATED ON NEXT SHEET




OVERALL SYSTEM
(ONLY MAIN COMMUNICATION LINES SHOWN)
SHEET 2 OF 4



NOTES:
 1 PHOENIX 45KVA SEE VENDOR MANUAL.
 2 OPTIONAL SHIM SUPPLY

OVERALL SYSTEM
 (ONLY MAIN COMMUNICATION LINES SHOWN)
 SHEET 3 OF 4




- 1 ALL POWER AND SIGNAL LINES ENTERING AND LEAVING THE RF ENCLOSURE PASS THROUGH THE PENETRATION PANEL. FOR EASE OF USE, ON THIS DIAGRAM, THEY ARE NOT SHOWN GOING TO THE PENETRATION PANEL.
- 2 S/C SHIM POWER SUPPLY AND S/C SHIM COILS FOR GE MAGNETS ONLY.
- 3 FOR MOBILE SYSTEMS, SENSOR IS LOCATED IN CEILING BOX.
- 4 OPTIONAL RESISTIVE SHIM SUPPLY AND COIL.
- 5 DESIGN DIFFERENCES BETWEEN 1.5T & 1.0T SYSTEMS.
- 6 OPTIONS AVAILABLE FOR 1.5 T SYSTEMS ONLY.
- 7 MULTI-COIL HARDWARE IS ONLY PRESENT WHEN MULTI-COIL COMPATIBILITY KIT HAS BEEN INSTALLED.
- 8 HORIZON SYSTEMS HAVE LITRON HEAT EXCHANGER MR/; AND CV/; FIXED SITE SYSTEMS HAVE NESLAB WATER CHILLER. MOBILES HAVE E & W HEAT EXCHANGER.
- 9 SIGNA MR/; AND CV/; SYSTEMS MANUFACTURED AFTER JUNE 15, 2000 HAVE MICROPHONE IN THE HEAD CARRIAGE INSTEAD OF THE REAR END BELL.
- 10 FOR TSSC AND TGWC SEE WATER CHILLER.

EMERGENCY RUNDOWN UNIT (ER1 FOR OXFORD MAGNET) (MS4 FOR GE MAGNET)
(SEE NOTE 4)

REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAM FOR ADDITIONAL INFORMATION:

- NOTE:
- 1) PATIENT HANDLING
 - 2) TPS/ISE
 - 3) GRADIENT
 - 4) MAGNET AND CRYOGENS
 - 5) REFER TO SPECTROSCOPY MANUALS
 - 6) OXYGEN MONITOR SUBSYSTEM (DIR. 15336).

OVERALL SYSTEM
(ONLY MAIN COMMUNICATION LINES SHOWN)




OPERATOR WORKSPACE
(OPERATOR WORKSPACE CABINET/TABLE – OCTANE)

SHEET 1 OF 2

- NOTES:
- 1 NETWORK CONNECTION AVAILABLE FOR EXTERNAL SUITE CONNECTION
 - 2 NETWORK CONNECTION AVAILABLE AT FRONT OF DUPLEX ASSEMBLY FOR THE SERVICE LAPTOP
 - 3 FOR VIDEO CONFERENCING
 - 4 MV618 PCI CARD HAS ON-BOARD FIBER OPTICS SO THERE IS NO BIT 3 FIBRE OPTIC INTERFACE.

REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAM FOR ADDITIONAL INFORMATION:

- NOTE:
- 1) PATIENT HANDLING
 - 2) TPS/ISE





REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAM FOR ADDITIONAL INFORMATION:



NOTE:
1) PATIENT HANDLING
2) TPS/ISE


OPERATOR WORKSPACE
(OPERATOR WORKSPACE TABLE)

SHEET 2 OF 2




REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAMS FOR ADDITIONAL INFORMATION.
 NOTES:
 1) TPS/ISE
 2) SYSTEM






REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAM FOR ADDITIONAL INFORMATION:

NOTE:
1) GRADIENT



NOTES:
1) REFER TO THE TPS/SE BLOCK DIAGRAMS FOR ADDITIONAL INFORMATION.

1 FIBER OPTIC CABLES ONLY FEEDTHROUGH THE CABINET I/F PANEL; THERE ARE NO TERMINATIONS THERE.



PDU

TABLE OF CONTENTS

Note:


Phoenix PDU found in Vendor Manuals only. (SEE SERVICE METHODS CDROM)

TPS/ISE (1.5T)

TABLE OF CONTENTS


<u>SECTION</u>	<u>PAGE</u>
SECTION 1 – ISE/RF BLOCK DIAGRAMS (1.5T)	1-1
SECTION 2 – RF SUPPLEMENTAL SCHEMATICS (1.5T)	
DUAL DIRECTIONAL COUPLER SCHEMATIC (46-264288S, sheets 1 & 2)	2-1 & 2-2
HEAD & BODY PREAMP SCHEMATIC (46-264442S)	2-3
DYNAMIC DISABLE FILTER MODULE SCHEMATIC (46-264776S)	2-4
HEAD T/R SWITCH SCHEMATIC (46-264988S)	2-5
BODY T/R SWITCH/SPLITTER SCHEMATIC (46-321314S)	2-6



SECTION 1 - ISE/RF BLOCK DIAGRAMS (1.5T)



REFER TO FOLLOWING BLOCK DIAGRAMS FOR ADDITIONAL INFORMATION.


- NOTE:
- 1) OPERATOR WORKSPACE
 - 2) MAGNET AND CRYOGENS
 - 3) POWER AND ON/OFF CONTROL
 - 4) PATIENT HANDLING







REFER TO FOLLOWING BLOCK DIAGRAMS FOR ADDITIONAL INFORMATION:

- NOTE:
- 1) GRADIENT
 - 2) PATIENT HANDLING
 - 3) OPERATOR WORKSPACE







TPS/ISE/RF SUBSYSTEM
UNIVERSAL COMBINED EXCITER RECEIVER DAB (UCERD)
 SHEET 5 OF 13



TPS/ISE/RF SUBSYSTEM (IPG BOARD) SHEET 6 OF 13




TPS/ISE/RF SUBSYSTEM
SRFD II CABINET – SRFD II RF AMPLIFIER
SHEET 7A OF 13



TPS/ISE/RF SUBSYSTEM
 (RF POWER CABINET – SYSTEM SUPPORT MODULE)
 SHEET 7B OF 13

REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAM FOR ADDITIONAL INFORMATION:


NOTE:
 1) PATIENT HANDLING



TPS/ISE/RF SUBSYSTEM
(RF POWER CABINET - SYSTEM SUPPORT MODULE)
SHEET 8 OF 13

REFER TO THE FOLLOWING FUNCTIONAL BLOCK DIAGRAM FOR ADDITIONAL INFORMATION:


NOTE:
1) GRADIENT




TPS/ISE/RF SUBSYSTEM
 (RF/PDU CABINET – RFI SOLID STATE AMPS)
 SHEET 9 OF 13


**ANALOGIC
RF
AMPLIFIER
#1**

**ANALOGIC
RF
AMPLIFIER
#2**




TPS/ISE/RF SUBSYSTEM
(MAGNET ROOM & RELATED COMPONENTS)





1 MULTIPLE WIRES/PINS WITH SAME SIGNAL ARE NOT SHOWN

SECTION 2 - RF SUPPLEMENTAL SCHEMATICS (1.5T)



46-264288-S 1 2 3 4 5

CONT ON SHEET 2 SHT NO. 1
46-264288-S
DRAWING NO.

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-

APPLIED PRACTICES				
46-208600				

REV	AH
46-264288-S	
CONT ON SHEET 2	SHT NO. 1

TITLE	46-264288-S
SCHMATIC DIAGRAM	CONT ON SHEET 2 SHT NO. 1
DUAL DIRECTIONAL COUPLER	
FIRST MADE FOR MRDDAMP	

THE REVISION NUMBER OF THIS SHEET IS THE CONTROL NUMBER OF THIS ENTIRE DOCUMENT. THIS REVISION NUMBER MUST INCREASE WHENEVER ANY SHEET IS CHANGED. THE CHANGE MUST BE RECORDED IN THE FORM BELOW.

PCN#	158051	158445	171514	171706	171687	
DATE	12MAY89	13MAR90	21FEB91	02MAY91	08MAY91	
SHT	2	AA	AB	AC	AD	AE
						AF

NOTES:
1. UNLESS OTHERWISE SPECIFIED:
ALL RESISTORS ARE 0.25W, 5%
AND ARE IN OHMS.

REVISIONS		PRINTS TO	
AD	FEP 16MAY89	AG	FEP 2MAY91
PCN	158051	PCN	171706
AE	DLW 13MAR90	AH	FEP 8MAY91
PCN	158445	PCN	171687
AF	DPB 21FEB91	NR	DLW 30JAN89
PCN	171514	ADDED	G2

MADE BY	C. LUEBBE 4AUG87
ISSUED	R. BECERRA 14OCT87

APPROVALS	R. BECERRA 14OCT87	MEDICAL SYSTEMS	DIV OR DEPT	46-264288-S
		MILWAUKEE, WI	LOCATION	CONT ON SHEET 2 SHT NO. 1


46-264288-S
 2 - ON SHEET - SHT NO. 2
 DRAWING NO.

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-
 APPLIED PRACTICES
 46-208600

REV AF
 46-264288-S
 CONT ON SHEET - SHT NO. 2

TITLE
 SCHEMATIC DIAGRAM
 DUAL DIRECTIONAL COUPLER
 FIRST MADE FOR MRDDAMP

46-264288-S
 CONT ON SHEET - SHT NO. 2



46-264288G2-B

REVISIONS		PRINTS TO	
AE	FEP 2MAY91	AB	FEP 16MAY89
	PCN 171706		PCN 158051
AF	FEP 8MAY91	AC	DLW 13MAR90
	PCN 171687		PCN 158445
		AD	DPB 21FEB91
			PCN 171514

MADE BY
 D. WALK 30JAN89
 ISSUED

APPROVALS
 MEDICAL SYSTEMS DIV OR DEPT
 MILWAUKEE, WI. LOCATION
 46-264288-S
 CONT ON SHEET - SHT NO. 2

46-264288-S

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-

APPLIED PRACTICES			
46-208600			

REV ~~AE~~ **AG AH**

46-264442-S

CONT ON SHEET - SHT NO. 1


TITLE

SCHEMATIC DIAGRAM
 PREAMP


FIRST MADE FOR MR TPS

MG2A16A3 A1 (HEAD)
 MG3A3A5 A1 (BODY)




46-264442G2-A
 46-264442G3-A



- NOTES:
- UNLESS OTHERWISE SPECIFIED:
 ALL RESISTORS ARE IN OHMS AND ARE 0.125W, 5%.
 ALL CAPACITORS ARE IN MFD AND ARE 50V, 10%.
 - "*" INDICATES SURFACE MOUNTED COMPONENT.
 - FOR COMPONENT VALUES SEE 46-264442 SECTION 9.1.




REVISIONS		PRINTS TO	
AE DLW 10DEC91	PCN 173213	740	
AG DRL 2 FEB 97	PCN 197402		
AF WJK 8 JAN 96	PCN 186807		
AC DLW 13JAN90	PCN 158383		
AD DLW 5FEB90	PCN 158676		
MADE BY R. LISOWSKI 14JAN88 ISSUED R. BECERRA 24JUN88		APPROVALS R.B. 20MAR89 MEDICAL SYSTEMS DEPT MILWAUKEE, WI. LOCATION DIV OR DEPT 46-264442-S CONT ON SHEET - SHT NO. 1	



NOTES:
 1. UNLESS OTHERWISE SPECIFIED:
 ALL RESISTORS ARE 3.0W, 1% AND ARE IN OHMS.
 ALL CAPACITORS ARE IN pF.
 ALL INDUCTORS ARE IN uH.
 ALL DIODES ARE UM9415.
 2. DARKER CONNECTIONS INDICATE RF SIGNAL PATHS.

REVISIONS		PRINTS TO
AD DRL 16 JUN 98	B DLW 5DEC89	740
PCN 197927	GEN. CHNGS.	
	AB DLW 20MAR90	
	PCN 158759	
	AC WJK 8 JUN 93	
	PCN 185690	

MADE BY D. WALK 27JUL89	APPROVALS RCR 12MAR90	MEDICAL SYSTEMS DEPT MILWAUKEE, WI	DIV OR 46-264988-S
ISSUED J. WALTER 23JAN90			LOCATION CONT ON SHEET - SHT NO. 1



46-321312
ASSEMBLIES

	R11	R12
G1, G3	100 ohm	100 ohm
G2, G6	50 ohm	NOT USED

1. THE FOLLOWING CONNECTIONS ARE MADE ON THE NEXT HIGHER ASSEMBLY.
 3. THE FOLLOWING COMPONENTS ARE NOT ASSEMBLED FOR GROUP 1: C28, C29, CR6, L12 AND R17.

BLOCK PATHNAME		/user/body_hyb/tr_switch_1.5t SHEET 1 OF 1			
REV 4	1.5T T/R BOARD		LOCATION CODE	APPROVALS	GE MEDICAL SYSTEMS
DRAWING NO.	FIRST MADE FOR	DATE	ISSUED	DATE	PRINTS TO
46-321314-S	MRPWB	31/MAR/93	DANIEL S. SMALL	3/MAY/93	740
REV. 1	MADE BY	DATE		DATE	
	DAVID LYLE	31/MAR/93		3/MAY/93	
		MILWAUKEE WI		R. JACOBS 30/ APR/ 93	
		PCN 185887		PCN 186271	
		PCN 186057		PCN 186408	

COOLING SYSTEM

TABLE OF CONTENTS

Water Chillers TGWC & TSCC Block Diagrams

Note:

For all Chillers and Cooling System variations (SEE Vendor Manuals on MR Service Methods CDROM)

