g

GE Medical Systems

Technical Publications

**Direction 2327645
Revision 0**

Signa Lx Patch for FIESTA2D Option Patch Software Load Instructions

Copyright© 2001 by General Electric Co.

Operating Documentation

DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

Immediately complete a "Damage Loss Claim Form" (available via Public Folders on the GEMS intranet) after the damage is found.

GEMS Home Page / Public Folders:

All Public Folders/Medical Systems/!Global Initiatives/Information Management/Forms /Common Forms/ DAMAGE LOSS CLAIM FORM.

Send the completed form to the email address listed in the form.

For more information about the Transportation Claim Procedure, access the GE Medical Systems Intranet and enter the following URL address (case sensitive):

<ftp://3.87.40.2/globepro/qualsys/Docs/190016MF.PDF>

Rev. 12/17/2001

Direction 2128126 - Language Policy For Service Documentation, Rev 0

WARNING

- THIS SERVICE MANUAL IS AVAILABLE IN ENGLISH ONLY.
- IF A CUSTOMER'S SERVICE PROVIDER REQUIRES A LANGUAGE OTHER THAN ENGLISH, IT IS THE CUSTOMER'S RESPONSIBILITY TO PROVIDE TRANSLATION SERVICES.
- DO NOT ATTEMPT TO SERVICE THE EQUIPMENT UNLESS THIS SERVICE MANUAL HAS BEEN CONSULTED AND IS UNDERSTOOD.
- FAILURE TO HEED THIS WARNING MAY RESULT IN INJURY TO THE SERVICE PROVIDER, OPERATOR OR PATIENT FROM ELECTRIC SHOCK, MECHANICAL OR OTHER HAZARDS.

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REVISION HISTORY

<u>REV</u>	<u>DATE</u>	<u>PRIMARY REASON FOR CHANGE</u>
A.....	Nov. 29, 2001	Initial version of procedure.
0.....	Dec. 17, 2001	Revised for Fiesta 2D specifically

<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>	<u>PAGE</u>	<u>REV</u>
Title Page	0						
Damage in Trans.....	-						
Direction 2128126 ...	0						
5 to 10	0*						

* This revision/letter corresponds to the indicated document's revision control system.

TABLE OF CONTENTS

1- OVERVIEW	7
2- LOADING SIGNA LX RELEASE 9.0, CNV3, CNV4 PATCHES FROM CD-ROM.....	7

1- OVERVIEW

The Signa LX Release 9.0, CNV3, CNV4 Patch03 CD contains critical updates for the Fiesta 2D option key.

- Enables the 2D fiesta feature when the 2D fiesta key is installed

Note

Signa software must be at the M4 version of Release 9.0, CNV3, or CNV4 to install this patch CD.

2- LOADING SIGNA LX RELEASE 9.0, CNV3, CNV4 PATCHES FROM CD-ROM

1. With the system idle (not scanning, filming or archiving – reboot if unsure)
2. Inside a C Shell window, log in as the root user by typing the following:


su <Enter>

Password: **operator <Enter>** ("operator" is the system default root password. Use whatever password is appropriate for the system you are working on if it has been changed.)

3. Place the Signa Patch CD in the system (not PC) CD drive. Wait 20 seconds for the system to see the CD until the light on the front of the CD drive stops blinking.
4. Inside the C Shell, type the following:

mediad <Enter> (This command mounts the CDROM.)

swmgr <Enter> (This command starts the software manager which is used to load the patch.)


SOFTWARE MANAGER INITIALIZATION
ILLUSTRATION 2-1

5. Once the Software Manager starts, begin typing the following (case sensitive) in the **[Software Location Edit Window]**. See Illustration 2-1.

/CDROM (The system will automatically enter **/dist** after you type /CDROM.)

6. Click the **[Lookup]** button found to the right of the **[Software Location Edit Window]**. See Illustration 2-1.
7. Click the **[Customize Installation]** button at the top of the screen. See Illustration 2-1.

SOFTWARE MANAGER CUSTOM INSTALLATION (EXAMPLE)
ILLUSTRATION 2-2


8. All available patches should be visible at this time in the selection window. Any patches that apply to your system and are not loaded will be selected (check marked). Verify the appropriate box is checked corresponding the system software version. Note: only the correct software version should be checked. See Illustration 2-2.

9. If the **[Start]** button activates, proceed to 9.a. If the **[Conflicts]** button activates proceed to section 9.b

9.a The **[Start]** button will activate (The letters in "Start" will change colors or shades). If the **[Start]** button activates, push it and the patch load will begin. A progress bar will open and count up from 0 to 100%. The load takes less than 1 minute. When the progress bar goes away, a confirmation window will be displayed. Select **[Exit]** to quit the Software manager.

9.b After the patch is selected, if there is a problem, the **[Conflicts]** button will activate, (the letters in "Conflicts" change colors or shades). If this happens, this indicates that your system does not need this patch and the indicated patch will not be loaded. Push the **[Conflicts]** to view the error if desired.

10. After installation, select **Exit, Installations completed** checkbox, and select **[OK]**. See Illustration 2-3.

INSTALLATION DECISION WINDOW
ILLUSTRATION 2-3

11. In the C Shell type the following (case sensitive) to eject the CD-ROM from the drive:

umount /CDROM <Enter> (Un-mounts the CD from the system)

eject /CDROM <Enter> (Opens the CD drawer)

12. If the patches were required and loaded, verify they are installed by typing in the C Shell:

getver <Enter> You should see an entry for the LX patch as follows (example only – 9.0 may be replaced with CNV3, CNV4):

Build number for MrpApps is 90.29F_M4_0139.z

Build number for PostSdC is 90.29F_M4_0139.z

Build number for cclass is 90.29F_M4_0139.z

Build number for driverSupport is 90.29F_M4_0139.z

Build number for insite is 2.3-mips

Build number for install is 90.29F_M4_0139.z

Build number for os_cd is 6_11.1.a

Build number for os_u_cd is 6u11.6.a

Build number for patch_90_M4_03 is 90.29F_M4_0139.z.PC

Build number for pc is 90.29F_M4_0139.z

Build number for z is Upgrade

MR Software release: 9.0.0139z

Note: CNV3 and CNV4 will have entries like this from the getver command:

Build number for patch_CNV3_M4_03 is . . .

Build number for patch_CNV4_M4_01 is . . .

12. Reboot Signa for the patch to take effect.