


### MR2 A15 A20 A3 BROADBAND EXCITER

46-264230G1-D, G2-E




#### Description

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:		REV. REF.		TITLE	
APPLIED PRACTICES	SURFACES	DECIMALS (2 PLACES)	DECIMALS (3 PLACES)	ANGLES	46-264230-S
46-28898B	✓				CONT ON SHEET 2 SH NO. 1

P/L ISSUED		REV. REF.		TITLE	
46-264230G1-D		46-264230G2-E		SCHEMATIC DIAGRAM BROADBAND EXCITER	
FIRST MADE FOR: MRSPEC-1.5					

THE REVISION NUMBER OF THIS SHEET IS THE CONTROL NUMBER OF THIS ENTIRE DOCUMENT. THIS REVISION NUMBER MUST INCREASE WHENEVER ANY SHEET IS CHANGED.

PCN #	DATE	BY	REASON
171751	1JUN91	AE	RA
158733	9FEB98	AR	AB
158952	21JUN90	AC	AC
171127	12SEP90	AD	AD
171496	6MAR91		





NOTES:  
 1. UNLESS OTHERWISE SPECIFIED:  
 ALL RESISTORS SURFACE MOUNT  
 1/8W, 5% AND ARE IN OHMS  
 CAPS-ALL FRACTIONAL VALUES - UF,  
 ALL OTHERS-PF  
 ALL INDUCTORS ARE IN NH.  
 ALL DIODES ARE 1N4148  
 2. \*INDICATES AN ACTIVE LOW.

REVISIONS		PRINTS TO
AE	GT_6MAR91	AB
PCN_171496		DLW_5MAR90
AF	DPB_1JUN91	PCN_158733
PCN_171751		AC
		FEP_21JUN90
		PCN_158952
		AD
		FEP_12SEP90
		PCN_171127

APPROVALS	ISSUED	MEDICAL SYSTEM	46-264230-S
RCA	SFEB98	MILWAUKEE, WISCONSIN	CONT ON SHEET 2 SH NO. 1

MR2 A15 A20 A3

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING :-		REV AE	TITLE
APPLIED PRACTICES	SURFACES	46-264230-S	SCHEMATIC DIAGRAM
46-208600	✓	CONT ON SHEET - SH NO. 2	BROADBAND EXCITER
TOLERANCES ON DIMENSIONS		FIRST MADE FOR MASPEC-4.5	
DECIMALS (2 PLACES)	DECIMALS (3 PLACES)	P/L ISSUED	




FOR G1: 46-301880P1  
 FOR G2: 46-301880P2

REVISIONS		PRINTS TO
AE DPB_1JUN91	AB DLW_3MAR90	740
PCN_171751	PCN_158733	
	AC FEP_21JUN90	
	PCN_158952	
	AD GT_6MAR91	
	PCN_171496	


DESIGNED BY E.P. FERRE_18OCT89	APPROVED BY RCR SFEB90	MEDICAL SYSTEMS MILLIS, MASSACHUSETTS	46-264230-S CONT ON SHEET - SH NO. 2
-----------------------------------	---------------------------	--	---

### MR2 A15 A18 A3 BROADBAND RECEIVER

46-264232G1-B, G2-B




#### Description



GAIN ADJUSTMENT CHART

ATTENUATION -dB	R20	R21	R22
0	JUMPER	OPEN	JUMPER
1	JUMPER	470 OHMS	6.2 OHMS
2	6.2 OHMS	220	6.2
3	9.1	150	9.1
4	12	100	12
5	12	82	15
6	15	68	18
7	18	56	22
8	22	47	22
9	24	39	24
10	27	36	27

DEVICES SHOWN TOP VIEW:




NOTES:  
1. UNLESS OTHERWISE SPECIFIED:  
ALL RESISTORS ARE 0.25W, 1%  
AND ARE IN OHMS.  
ALL CAPACITORS ARE IN MFD.  
ALL INDUCTORS ARE IN nH.

REVISIONS	PRINTS TO
NR/DLW 19JAN90	740
ADDED G2	
AB/DLW 30MAR90	
PCN 158784	

**MG2 A16 A6**  
**MG2 A16 A7 A3**  
**SPECTROSCOPY ADAPTOR**

46-264272G1-A



**Description**


46-264272-S  
CONT ON SHEET - SHT NO. 1  
DRAWING NO.

46-264272-S  
CONT ON SHEET - SHT NO. 1

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-			REV AB	TITLE
APPLIED PRACTICES			46-264272-S	SCHEMATIC DIAGRAM SPECTROSCOPY ADAPTOR
46-208600			CONT ON SHEET - SHT NO. 1	FIRST MADE FOR SPECTROSCOPY COIL QUICK DISCONNECT

MG2A16 A6  
MG2A16 A7 A3

46-264272G1-A



REVISIONS		PRINTS TO	
AB	DPB 24OCT90	740	
	PCN 171248		

MADE BY R. LISOWSKI 10JUN87	APPROVALS RCR 27JUL87	MEDICAL SYSTEMS MILWAUKEE, WI	DIV OR DEPT LOCATION	46-264272-S
ISSUED B. AHONEN 14JUL87				CONT ON SHEET - SHT NO. 1

46-264272-S<sup>1</sup>

2


3

4

5

### MG2 A16 A7 A2 SPECTRO T/R SWITCH

46-264762G1-B




#### Description

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING--		REV AC
APPLIED PRACTICES	46-208600	46-264762-S
		CONT. ON SHEET - SHT NO. 1

TITLE		46-264762-S
SCHEMATIC DIAGRAM		
SPECTRO T/R SWITCH		
FIRST MADE FOR MR SPECT 4.5		

46-264762G1-B




NOTES:  
1. UNLESS OTHERWISE SPECIFIED:  
ALL CAPACITORS ARE IN MFD.  
ALL DIODES ARE UM7001E.  
ALL RESISTORS ARE 0.4W, 5%  
AND ARE IN OHMS.

REVISIONS	PRINTS TO
B   DLW 12JUN89	740
GEN. CHNGS.	
AB   DLW 20FEB90	
PCN 158677	
AC   FEP 12FEB91	
PCN 171522	

MADE BY G. TESKE 17MAY89	APPROVALS RCH 20OCT89	MEDICAL SYSTEMS DIV DEPT MILWAUKEE, WI	46-264762-S CONT. ON SHEET - SHT NO. 1
-----------------------------	-----------------------------	---	---

**PP1 A13**  
**SPECTRO BROADBAND 20dB GAIN BLOCK**

46-264814G1-A




**Description**

The Spectroscopy BroadBand 20 dB Gain Block was designed to augment receive channel gain from the RF Preamp to the RF Receiver. Located at the Penetration Panel it provides roughly 20 dB of signal amplification at all frequencies of multinuclear interest while ensuring a total noise figure of less than 1 dB for the composite receive channel. In general, the +15 VDC bias applied to port J2 (Output) enables the active device for proper RF operation over the frequency range. Capacitors C6 and C10 accomplish DC blocking at the input and output of the device while inductors L5, L8, and L11 act as RF chokes. Capacitors C7 and C12 route any residual RF in the DC bypass path to ground. This configuration bridges the DC supply from ports J2 (Output) to J1 (Input) so as to supply DC current to the preamp.

As a continuity check with no bias applied, the DC resistance at ports J1 and J2 to Ground (common) should be observed infinite (1 to 10 Megaohms). The resistance measured through the device, J1 center conductor to J2 center conductor, should be observed at 4 to 6 ohms.

EXPLANATION OF +15 VDC PATH FROM RECEIVER BD. TO SPECTRO PREAMP



CONT ON SHEET - SHT NO 1  
 46-264814-S  
 DRAWING NO.


PP1 A13

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-  
 APPLIED PRACTICES

REV A EAA  
 46-264814-S  
 CONT ON SHEET - SHT NO. 1


GENERAL ELECTRIC  
 46-264814-S  
 CONT ON SHEET - SHT NO 1  
 TITLE  
 SCHEMATIC DIAGRAM  
 SPECTRO BROAD BAND 20dB GAIN BLOCK  
 FIRST MADE FOR MR  
 P/L ISSUED

46-264814G1-A



NOTES:  
 1. UNLESS OTHERWISE SPECIFIED:  
 ALL CAPACITORS ARE IN MFD.  
 2. "\*" INDICATES SURFACE MOUNTED  
 COMPONENT.

ALL DEVICES SHOWN TOP VIEW




REVISIONS		PRINTS TO	
		740	

MADE BY G. TESKE 10MAR89  
 ISSUED  
 APPROVALS  
 MEDICAL SYSTEMS DEPT MILWAUKEE, WI LOCATION  
 46-264814-S  
 CONT ON SHEET - SHT NO. 1

**MR6 A2 A4  
RF INTERFACE**

46-264818G1-C




**Description**

The RF Interface Board connects to a bi-directional communication (Bicycle) link via connector J6. The RF Interface Board serves as an interface that controls/monitors the Spectroscopy amplifier when commanded by the host computer and it links up with the power monitor to inhibit the Spectroscopy amplifier output in the event an unsafe condition is detected. This board communicates with the power monitor by sending an echo of the Unblank signal. If the unsafe condition is detected an Rflock signal is transmitted to the RF Interface Board which in turn generates a shutdown signal to inhibit the amplifier output. As long as the shutdown signal is not cleared, the Unblank signal will remain inhibited from being applied to the amplifier.

- J1 — connects to the band selectable filter (Fist1rtn, Fist2rtn, Fist3rtn, +15 VDC)
- J2 — connects to the Power Monitor (Ublks\*-p/n, Rflck-p/n, Rflck\*-p/n)
- J3 — connects to the dynamic disable board (UBLNS\*-P/N)
- J4 — connects to the power supply (+5 VDC, +15 VDC, Common)
- J5 — connects to the BroadBand RF Amplifier
- J6 — connects to the bi-directional Bicycle link

APPLIED PRACTICES 46-206600	SURFACES ✓	TOLERANCES ON DIMENSIONS DECIMALS (2 PLACES)	DECIMALS (3 PLACES)	ANGLES
UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING :				

REV. NO. 46-264818-S CONT. ON SHEET 2 SH. NO. 1	TITLE SCHEMATIC DIAGRAM RF_INTERFACE_BOARD FIRST MADE FOR SPECTROSCOPY_AMPLIFIER P/L ISSUED 46-264818G1-C
---	---





- NOTES:
- UNLESS OTHERWISE SPECIFIED:  
ALL RESISTORS ARE 0.25W, 5%  
AND ARE IN OHMS  
ALL CAPACITORS ARE IN MFD  
ALL INVERTERS ARE 1N7404  
ALL DIODES ARE 1N4148
  - # INDICATES AN ACTIVE LOW.

THE REVISION NUMBER OF THIS SHEET IS THE CONTROL NUMBER OF THIS REVISION DOCUMENT. THIS REVISION NUMBER MUST INCREASE WHENEVER ANY SHEET IS CHANGED.

PCN #	DATE	BY	REV	DESCRIPTION
1	24JUL89	AA	AB	158492
2	21SEP89	AA	AB	158687
3	16JAN90	AB	AC	171340
4	27FEB90	AC	AD	31DEC90
5	31DEC90	AD	-	-


REV	DATE	BY	DESCRIPTION	PRINTS TO
AB	GJO	16JAN90	740	
AC	FEP	27FEB90		
AD	FEP	31DEC90		



REVISIONS		PRINTS TO
RB	GJ0...16JAN90	740
PCN	...158492	
AC	FEP...27FEB90	
PCN	...158687	
AD	FEP...31DEC90	
PCN	...171340	

### MR1 A3 A7 POWER MONITOR I/F

46-264820G1-B



#### Description

#### Note

FMI 60420 REMOVES THIS BOARD (M1040FF) ONLY.

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING :-		REV. AB	46-264820-S	TITLE SCHEMATIC DIAGRAM
APPLIED PRACTICES 46-286500	SURFACES ✓	DECIMALS (2 PLACES)	CONT. ON SHEET 2 SH. NO. 1	POWER-MONITOR-1/F FIRST MADE FOR MR-SPECT-4.5
		TOLERANCES ON DIMENSIONS DIMENSIONS (3 PLACES)		P/L ISS 46-264820

FROM BROAD BAND RF CABINET

FROM NARROW BAND RF CABINET

FROM BROAD BAND RF CABINET


FROM NARROW BAND RF CABINET

TO ANALOG SENSE BD (POWER MONITOR)

TO ANALOG SENSE BD (POWER MONITOR)

THE REVISION NUMBER OF THIS SHEET IS THE CONTROL NUMBER OF THIS ENTIRE DOCUMENT. THIS REVISION NUMBER MUST INCREASE WHENEVER ANY SHEET IS CHANGED.


PCN #	DATE	REV	BY
-	11-10-89	11-10-89	1MAY90
2	A	AA	AB



MR1 A3 A7

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING :-		REV. A	46-264820-S	TITLE SCHEMATIC DIAGRAM
APPLIED PRACTICES 46-206600	SURFACES ✓	DECIMALS (2 PLACES)	DECIMALS (3 PLACES)	ANGLES
		CONT. ON SHEET - SH. NO. 2		FIRST MADE FOR MR-SPECT-4.5 P/L ISS

A  
B  
C  
D  
E  
F



TO BROAD  
BAND  
RF  
CABINET


TO  
NARROW  
BAND  
RF  
CABINET

TO CONTROL BOARD  
(POWER MONITOR)

SCHEMATICS/INTERCONNECTS

**MG2 A16 A7 A2**  
**SPECTROSCOPY PREAMP (Phosphorus)**

46-264862G1-B



**Description**

The preamp provides the first stage of amplification for the MR signal in the receive path. The preamp is part of the Spectroscopy TR Module which is placed at the front of the head carriage slide trolley when Multinuclear BroadBand scanning is selected. The J2 RF output port of the preamp requires +15 VDC, +/- 10% @ 100 mAmps. The +15 VDC preamp bias originates in the Systems Cabinet TPS Power Supply and is supplied to the Receiver Board and BB Receiver Module via the backplane. The Spectroscopy AUX port preamp bias will measure +15 VDC after a Broad-Band Spectroscopy protocol is selected and pulsed once. The spectroscopy preamp bias is sent down the selected AUX receive coaxial cable to power the 20 dB gain block and the preamp. The total nominal gain of the preamp is +35 dB, +/- 3 dB, this gain can be measured by injecting an RF signal at the specified frequency with a maximum level of -10 dBm (nominally this is -30 dBm). The noise figure should be less than 0.5 dB.

CONT ON SHEET - SHT NO. 1  
 S-298492-97  
 ON 9NIMV8D

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-  
 APPLIED PRACTICES  
 46-208600


REV AB  
 46-264862-S  
 CONT ON SHEET - SHT NO. 1

TITLE  
 SCHEMATIC DIAGRAM  
 SPECTROSCOPY PREAMP (PHOSPHORUS)  
 FIRST MADE FOR MR

46-264862-S  
 CONT ON SHEET - SHT NO 1

MG2A16A7 A2

46-264862G1-B



OPERATING FREQUENCY  
 (25.8 MHZ)

MAXIMUM  
 RESISTANCE  
 DC =  
 8.1 OHMS



MAXIMUM  
 RESISTANCE  
 DC =  
 8.1 OHMS

MAXIMUM  
 RESISTANCE  
 DC =  
 8.1 OHMS

MAXIMUM  
 RESISTANCE  
 DC =  
 8.1 OHMS

- NOTES:
- UNLESS OTHERWISE SPECIFIED:  
 ALL RESISTORS ARE IN OHMS AND  
 ARE 0.125W, 5%.  
 ALL CAPACITORS ARE IN MFD AND  
 ARE 50V, 10%.
  - "\*" INDICATES SURFACE MOUNTED  
 COMPONENT.

ALL DEVICES SHOWN TOP VIEW




REVISIONS		PRINTS TO	
AB	DLW 16FEB90	740	
PCN	158690		

MADE BY G.TESKE 7MAR89  
 ISSUED R.BECERRA 2MAY89  
 APPROVALS RCR 4MAY89  
 MEDICAL SYSTEMS MILWAUKEE, WI  
 DIV OR DEPT LOCATION  
 46-264862-S  
 CONT ON SHEET - SHT NO. 1

46-264862-S

**MR6 A2 A2**  
**THREE BAND LOW PASS FILTER**

46-264866G1-B



**Description**

The Three Band Low Pass Filter is used to remove the harmonics at the output of the Spectroscopy RF Amplifier. The Three Band Low Pass Filter has three independent low pass filters. A set of relay switches select which of the three frequency bands the RF will pass through.

The first band (25 MHz) has input J5 and output J6. The loss within the pass-band (DC to 24 MHz) is  $\leq$  0.25 dB. The stop-band has a minimum attenuation of 40 dB for frequencies greater than 50 MHz.

The second band (40 MHz) has input J3 and output J4. The loss within the pass-band (DC to 38 MHz) is  $\leq$  0.25 dB. The stop-band has a minimum attenuation of 40 dB for frequencies greater than 80 MHz.

The third band (64 MHz) has input J1 and output J2. The loss within the pass-band (DC to 64 MHz) is  $\leq$  0.25 dB. The stop-band has a minimum attenuation of 40 dB for frequencies greater than 130 MHz.

The remaining connector J7 is used to introduce the DC bias for the Spectroscopy TR Module and coil. This DC bias originates in the RF Cabinet as the Spectro TR Bias signal and is controlled by the unblank signal.


MR6 A2 A2

GENERAL ELECTRIC

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING :-		REV AC	
APPLIED PRACTICES 46-264866	SURFACES ✓	DECIMALS (2 PLACES)	INCHES (16)
		46-264866-S	
		CONT. ON SHEET - SH. NO. 1	

TITLE SCHEMATIC DIAGRAM  
THREE\_BAND\_LOW\_PASS\_FILTER  
FIRST MADE FOR RF\_CABINET

P/L ISSUED  
46-264866G1-B




REVISIONS	PRINTS TO
AB   DPB_22OCT90	740
PCN_171257	
AC   DPB_24MAY91	
PCN_171761	

MADE BY F. PIERCE_13JUL89	APPROVALS RCH 30OCT89	REGIONAL SYSTEMS MILWAUKEE, WISCONSIN	46-264866-S CONT. ON SHEET - SH. NO. 1
------------------------------	-----------------------------	--	---

**MR6 A2 A7  
RF INPUT**

46-264920G1-A



**Description**



The RF Input board's purpose is to accept +10 dBm of RF at J1 from the exciter and to provide an adjustable RF Input to the BroadBand RF Amplifier or the Narrowband RF Amplifier. The J3 output port sends the RF signal to the Narrowband RF Amplifier, AT38 adjusts this RF signal level range between +14 to -1 dBm. The J2 output port sends the RF signal to the Multi-Nuclear BroadBand RF Amplifier, AT25 adjusts this RF signal level range between +6 to -9 dBm.

CONT ON SHEET - SHT NO. 1  
S-026492-97  
MR6A2 A7

46-264920-S  
CONT ON SHEET - SHT NO. 1

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-				REV AB	TITLE
APPLIED PRACTICES				46-264920-S	SCHEMATIC DIAGRAM
46-208600				CONT ON SHEET - SHT NO. 1	RF INPUT
					FIRST MADE FOR MR

P/L ISSUED  
46-264920G1-A




NOTES:  
1. UNLESS OTHERWISE SPECIFIED:  
ALL RESISTORS ARE IN OHMS  
ALL CAPACITORS ARE IN MFD

REVISIONS		PRINTS TO	
	BDLW 14 JUN 89	740	
	GEN. CHNGS.		
	ABDLW 20 JAN 90		
	PCN 158634		

MADE BY G. TESKE 19 MAY 89	APPROVALS RCR 13 SEP 89	MEDICAL SYSTEMS MILWAUKEE, WI	DIV OR DEPT LOCATION	46-264920-S
ISSUED R. MAROLIA				CONT ON SHEET - SHT NO. 1

**MR6 A2 A5  
SENSOR DRIVER**

46-288004G1-A



**Description**

**Note**

**FOR M1040FF:** The Sensor Driver Board was removed by FMI 60420, SPECTRO POWER MONITOR UPGRADE (FOR PRE 5.4 PRODUCTION SYSTEMS, M1040FF). If Sensor Driver Board is present in the system then FMI 60420 needs to be performed on the system (4.X – 5.3 HARDWARE / M1040FF ONLY).




**FOR M1040FK:** The Sensor Driver Board may be present but is not connected or used.

46-288004-S

MR6 A2 A5

UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-		REV	AA	TITLE	46-288004-S
APPLIED PRACTICES	46-288000			SCHEMATIC DIAGRAM	
				SENSOR DRIVER	
				FIRST MADE FOR MR SPECT 4.5	

46-288004G1-A



NOTES:  
 1. UNLESS OTHERWISE SPECIFIED:  
 ALL RESISTORS ARE 0.25W, 1%  
 AND ARE IN OHMS  
 ALL CAPACITORS ARE IN MFD

REVISIONS	PRINTS TO
	740

MADE BY	D. WALK 17AUG89	APPROVALS	MEDICAL SYSTEMS DEPT MILWAUKEE, WI	DIV OF	46-288004-S
ISSUED				LOCATION	CONT ON SHEET - SHT NO. 1

46-288004-S


**POWER I/F PANEL**

46-301698S

**Description**

46-301698S  
DRAWING NO.

THIRD ANGLE PROJECTION



UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-

APPLIED PRACTICES	SURFACES	TOLERANCES ON DIMENSIONS		
	✓	DECIMALS (2 PLACES)	DECIMALS (3 PLACES)	ANGLES
		+/-	+/-	+/-

REV. ~~A~~ B 0  
46-301698S  
CONT ON SHEET - SH NO. 1

GE Medical Systems  
46-301698S  
CONT ON SHEET - SH NO. 1  
TITLE  
SCHEMATIC  
FIRST MADE FOR MRSPECT4.5


//STA6/USER/MRSPECT4.5/S301698D01

P/R

U/O

TO  
ENI  
AMPLIFIER  
AC LINE

TO  
PDU (PD1)  
A2TS1-4,5,6



APPROVALS		
	NAME	DATE
TECH		
DES		
ENG		
MFG		

MADE BY  
MITCH. MARAGOS 90MAR19  
ISSUED

REVISIONS			PRINTS TO	
B	MITCH. MARAGOS 90AUG15	NCN	742	
208 & 110 VAC 1 PHASE WERE 3 PHASE. ADDED Z,Y, X & N TO J1. J14N ADDED. CB2 & CB3 LINE & LOADS ADDED. J12 L1 & L2 ADDED. MR6 A3 CONNECTION POINTS TABLE DELETED.				
APPROVALS			GE MEDICAL SYSTEMS MILWAUKEE, WISCONSIN	
			DIV. OR DEPT. LOCATION	
			46-301698S CONT ON SHEET - SH NO. 1	

CAD

**FILTER RELAY**

46-301704S

**Description**

The Filter Relay Assembly is controlled through the software selection of the frequency and/or nuclei (the NUC Control Variable ultimately overrides any frequency selection). The RF enters into the Assembly at a common point and exits the Assembly at a common point. Three sets of relays are available which directly correspond to three frequency sensitive Low-Pass Filter networks. The Filter Relay Assembly essentially steers the BroadBand RF frequency to the proper frequency sensitive circuitry (3 Band Low-Pass Filter).

46-301704S  
DRAWING NO. 46-301704S  
CONT. ON SHEET - SH. NO. 1

2


3

4

5

GE Medical Systems


46-301704S  
CONT. ON SHEET - SH. NO. 1



UNLESS OTHERWISE SPECIFIED USE THE FOLLOWING:-  
APPLIED PRACTICES SURFACES TOLERANCES ON DIMENSIONS  
DECIMALS 2 PLACES DECIMALS 3 PLACES ANGLES

TITLE  
**FILTER, RELAY**  
FIRST MADE FOR MRSPECT4.5

A  
B  
C  
D  
P/R  
E  
U/O



DRAWING NO. 46-301704S  
CONT. ON SHEET - SH. NO. 1  
REV. A

APPROVALS		
	NAME	DATE
TECH		
DES		
ENG		
MFG		

REVISIONS		PRINTS TO	
		740	

MADE BY A. GRULKE 89NOV30  
ISSUED APPROVALS GE MEDICAL SYSTEMS MILWAUKEE, WISCONSIN DIV OR DEPT. LOCATION 46-301704S  
-CONT. ON SHEET - SH. NO. 1

CAD

2

3


4

5

6

**MG2 A16 A7 A2  
SPECTRO TR SWITCH**


46-321316G1-B



**Description**

The Spectro TR Switch is used to switch between the transmit and receive modes. The UNBLANK signal drives the TR Bias voltage mode. A bias voltage that originates from the TR Driver circuitry (under UNBLANK control) in the RF Amplifier Cabinet is used to either forward or reverse bias PIN diodes. When transmitting, a positive voltage is supplied via the transmit heliach. This forward biases the PIN diodes which in turn connect the transmitter to the coil. It also shorts the input to the preamplifier, which protects it from being damaged by the high level of RF present during the transmit cycle. When receiving, a negative voltage is supplied via the transmit heliach. This reverse biases the PIN diodes which in turn disconnects the transmit heliach from the coil and connects the preamplifier to the coil.

REVISION	PCN NUMBER	DATE	SIGN OFF
0	-----	27 / OCT / 93	



**NOTES :**

1. THE GROUP NUMBER AND GROUP REVISION FOR THIS PWA IS 46-321316G1-A.
2. THIS SCHEMATIC MEETS THE SPECIFICATIONS OF GE SCHEMATIC STANDARDS 46-017486.
3. UNLESS OTHERWISE SPECIFIED:  
RESISTORS ARE IN OHMS  
CAPACITORS ARE IN FARADS.  
INDUCTORS ARE IN HENRIES.

BLOCK PATHNAME		/user/body_hyb/spectr_sw SHEET 1 OF 1			
REV 0	SPECTRO T/R SWITCH	LOCATION CODE	APPROVALS	GE MEDICAL SYSTEMS	REVISIONS
DRAWING NO. 46-321316-S	FIRST MADE FOR MRSPECT4.5 (31P)	MG2-A16-A7-A1		MILWAUKEE WI	
BY: 1	MADE BY: Bill Kostolni	DATE: 27-OCT-93	ISSUED	DATE	PRINTS TO 740

Body Gain—>  $7.5K/5.11K * 20K/2K * (1+ 14300/1580) = 1.4677 * 10 * 10.051 = 147.5$   
 If input signal = 40 mVDC \* 147.5 =about 5.900 VDC

Head Gain—>  $21K/10K * 20K/2K * (1+ 14300/1580) = 2.1 * 10 * 10.051 = 211.1$   
 If input signal = 40 mVDC \* 211.1 =about 8.443 VDC

