

TABLE OF CONTENTS

TABLE OF CONTENTS	1
1- SST 1.5T PHOSPHORUS OVERVIEW	2
1-1 SST Test Hardware	2
2- SST 1.5T PHOSPHORUS TEST SET-UP	3
3- ³¹ P SST SCAN PRESCRIPTION.....	4
4- HELP	5
5- SYSTEM RESTORATION.....	5
REVISION HISTORY.....	6

1- SST 1.5T PHOSPHORUS OVERVIEW

The SST 1.5T Phosphorus Coil (~25.85 MHz) may be run to check Multi-Nuclear Spectroscopy system functions, excluding the 1.5T Phosphorus TR Flex Coil and associated MNS Quick Disconnect Adapter Box.

This procedure will attempt to test the receiver (via the AUX1 input), and the MNSpectroscopy TR Switch / preamplifier (housed in the ³¹P Spectroscopy Module Assembly). This represents the receive circuitry of the MNS system. This procedure will not test the 1.5T Phosphorus TR Flex Coils or associated MNS Quick Disconnect Adapter Box.

This procedure does not attempt to test the transmit portion of the MNSpectroscopy system. As the transmit portion of the system can be evaluated using the RF Power Measurement Kit (use the Set-up and Calibration section RF Out protocol (change the RF pulse on the User CV's Screen to "hard" or square pulse when Troubleshooting).

The test principle is to pulse the RF into a small sample (14.7 M H₃ PO₄ solution / coil), and then determine if the received echo is present or not. If a receive signal is not present, it is advisable to determine if the transmit signal is present. If the Transmit signal is present and meets acceptable parameters the receiver should be swapped (swap with another receiver module) as a next step. See Analysis for more help.

1-1 SST Test Hardware

The Universal SST Kit (46-287357G2) is a collection of all items needed to set up and perform the various SST test modes. The major items required for testing the 1.5T Phosphorus are listed in Tables 1-1 and 1-2.

TABLE 1-1
UNIVERSAL SST KIT 46-320383G2 HARDWARE

ITEM	DESCRIPTION	PART NUMBER	QUANTITY
1	Universal Small Sample 1.5T Phosphorus (~25.85 MHz) 14.7M H ₃ PO ₄	46-320153G3	1
2	Quick Disconnect Adapter for 1.5T systems	46-282468G3	1

TABLE 1-2
ADDITIONAL HARDWARE

ITEM	DESCRIPTION	PART NUMBER	QUANTITY
1	Grafidy Holder Base Plate	46-271410G1	1

Note

The G1 version of the SST Kit does not contain the additional coil for MNSpectroscopy (49-320153G3, 1.5T Phosphorus).

2- SST 1.5T PHOSPHORUS TEST SET-UP



PROPERTY DAMAGE! TO PREVENT COIL AND ASSOCIATED SWITCH DAMAGE, REMOVE ALL PHANTOMS AND HARDWARE (I.E., HEAD COIL, SURFACE COIL...) FROM THE MAGNET BORE.

1. Click on **[New Pt]**.
2. Move the cradle to the "HOME" position.
3. Remove all phantoms and coils (and associated hardware) from the bore.
4. Remove the head holder from the cradle.
5. Place the Grafidy Holder Base Plate on patient table near head end of cradle.
6. Place the Universal Small Sample "1.5T Phosphorus" (14.7M H₃ PO, 46-320153G3) Coil into the "center position" on the Grafidy Holder Base Plate.
7. Unwrap entire length of cable from around coil.
8. With the "1.5T Phosphorus" Coil column on base plate, turn brass lever to lock in place.



POISON HAZARD! SAMPLE CONTAINS PHOSPHORIC ACID, A SUSPECTED CARCINOGEN. DO NOT INGEST. DISPOSE OF AS A HAZARDOUS WASTE ACCORDING TO STATE AND FEDERAL REGULATIONS.

9. Connect cable from coil to the 1.5T adapter (use Extremity Coil/Linear Head Coil Adapter or 1.5T Service Tool Interface).
10. Insert Quick Disconnect Adapter into the "1.5T Spectroscopy Phosphorus Module" (white box used for Multi-Nuclear Spectroscopy).
11. Verify the MNSpectroscopy Module (white box) is connected properly to the transmit and receive cable-take-up lines (connections vary for depending on carriage cover type...)
12. Landmark at center of base plate/coil.
13. Press MOVE TO SCAN to send the coil to isocenter.

3- ³¹P SST SCAN PRESCRIPTION

- Set up Phosphorus (³¹P) SST Test scan protocol as follows:

TABLE 3-1
SCAN PROTOCOL: ³¹P SST TEST

[New Series]		<u>ACQUISITION TIMING</u>	
<u>PATIENT INFORMATION</u>		NEX	[1.00]
Patient Id	geservice	Freq DIR	[R/L]
Patient Name	31p sst	Auto CF	[Water]
Weight (Lb)	300 ~ <i>IMPORTANT</i>	Autoshim	[OFF]
<u>PATIENT POSITION</u>		<u>SCANNING RANGE</u>	
Patient Position	[Supine]	FOV	[20]
Patient Entry	[Head First]	Slice Thickness	[20]
Coil [Coil Type] [Surface] [31P_FLEX] [Accept]		Spacing	[1.5]
<u>IMAGING PARAMETERS</u>		Start	0
Plane	[Axial]	End	0
Mode	[2D]	# Slices	1 (default)
Pulse Seq	[Fid CSI (MRS)]	L/R Center	0 (default)
	[Accept]	P/A Center	0 (default)
Imaging Options	[Extended Dynamic Range]	Table Delta	0.00 (default)
	[Accept]	<u>(lower portion of monitor)</u>	
PSD Name	none (default)	[Save Series]	
<u>SCAN TIMING</u>		[Prepare To Scan]	
# of Echoes	1 (default)	[Spectro Prescan]	
TR	[1500]	Entry Point	[single1]
<u>USER CVs SCREEN</u>		Nucleus	[31]
CV0 spectral width	2500	R1	[9]
CV1 number of points	1024	R2	[28]
CV2 nucleus	31	TG	[50]
CV3 scan mode	1.00	AX	[Default] (~25.85 MHz)
CV4 total # of scans	32.00	Top Screen	[Magnitude] [Hz]
CV5 rl resolution - csi scans	1.00	Bottom Screen	[Q Chan Raw] [Pts]
CV6 ap resolution - csi scans	1.00	[Start] (we saw 3 x10e8 Magnitude signal)	
CV7 si resolution - csi scans	1.00	[Stop]	
CV14 rfpulse (soft)	1.00	[Done]	
CV16 CSI grid (acq)	2.00		
	[Accept]		

4- HELP

The following is a list of suggestions for troubleshooting if a receive signal echo is not present.

Note

SST is not used to measure SNR.

1. If a receive signal is present the 1.5T Phosphorus TR Flex Coil and associated MNS Quick Disconnect Adapter Box are potentially good.
2. If a receive signal is not present then troubleshoot:
 - ❖ Verify +15 VDC (receive bias) is present at the input to the MNS TR Module on the receive side lemo cable.
 - ❖ Determine if the Transmit signal out to the MNS Q. D. Adapter Box is good. Measure the RF Out of the MNS RF Amplifier at TG=200 using the MNS Set-up and Cal protocol, however, set the “rfpulse” in the User CV’s page to hard (from sinc to square). Next, measure the MNS signal at the rear pedestal cable-take-up bracket (there will be approximately -1.0 dB to -1.7 dB of cable loss). Next, procure a Q.D. Box with 2 BNC connections. The Quad Normal is acceptable (this adds in some other load circuitry not normally used located in the MNS TR Module). Connect test equipment to the transmit side of the Q.D. Box and measure. RF value should be very close to the value measured at the rear pedestal. The RF pulse should be square.
 - ❖ If transmit circuitry (see overview section) is good then swap receivers.

5- SYSTEM RESTORATION

1. When completed click on **[End Exam]**.
2. Verify scan desktop icon is displaying the “Idle” message.
3. Remove all test equipment.
4. Successfully complete one body scan.
5. Successfully complete one head scan.

