GE Healthcare

Direction 24239893TPH
1.5T HNS Coil Installation

Class A

Direction 24239893TPH
Rev 4
© 2008, General Electric Company,
All Rights Reserved

Revision History

Date	Rev	Change Description
5-Sept-08	0	Updated to meet compliance
27-Oct-08	1	Initial Release
19-March-09	2	<ul style="list-style-type: none"> • Updated Table of Contents with the “Scope” section • Updated “Damage in Transportation” statement • Extended the scope of installation to HDx and Discovery MR450 & MR750 • Updated the “Coil Config” section with a note related to pop-up box specific to Discovery MR systems • Removed the HDx designator from Section 5 • Update Table 5.3.1 with phantom part number • Added Section 6 - HNS Coil Setup for Reverse Ramped System
9-Sept-09	3	<ul style="list-style-type: none"> • Updated Language Policy per DOC0620154 • Added information per Phantom Elimination Program
8-June-10	4	<ul style="list-style-type: none"> • Removed MCQA section and referenced service docs per SPR MRIhc46403 • Updated Language Policy to DOC0620154 rev 2 • Updated Damage in Transportation • Removed requirement to install coil configs

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ATTENTION!

Installation of this coil MUST be performed by a GE field engineer. Please contact your field engineer to arrange for installation.

Attempts to perform installation of this coil without proper system configuration will result in system errors.

Table of Contents

1- DAMAGE IN TRANSPORTATION 11

2- Scope 12

3- COIL HARDWARE CONFIGURATIONS 13

4- COIL CONFIG..... 14

 4-1 Personnel Requirements 14

 4-2 Overview 14

 4-3 Preliminary Requirements 14

 4-4 Procedure 14

5- HNS Coil MCQA Setup..... 18

6- HNS Coil Setup for Reverse Ramped System..... 19

 6-1 Personnel Requirements 19

 6-2 Overview..... 19

 6-3 Tools and Test Equipment 19

 6-4 Procedure 19

1- DAMAGE IN TRANSPORTATION

All packages should be closely examined at time of delivery. If damage is apparent, have notation "**damage in shipment**" written on **all** copies of the freight or express bill **before** delivery is accepted or "signed for" by a General Electric representative or a hospital receiving agent. Whether noted or concealed, damage **MUST** be reported to the carrier **immediately** upon discovery, or in any event, within **14** days after receipt, and the contents and containers held for inspection by the carrier. A transportation company will not pay a claim for damage if an inspection is not requested within this **14** day period.

File a report with

- Call 1-800-548-3366 and use option 8.
- Fill out a report on <http://egems.med.ge.com/edq/home.jsp>
- Contact your local service coordinator for more information on this process.

Rev. October 28, 2009

2- SCOPE

This coil installation service manual is applicable to the HDx/HDxt Upgrade System, HDx/HDxt Forward Production System, Discovery MR450, and Optima MR450w.

The HNS Coil has 5 rating plates so there are 4 sets of Product Locator cards included with the shipment.

The Global Install Base Database (also known as Product Locator System) tracks shipment, trans-shipment, and field location of the serialized models. There are now two methods for submitting Product Locator information.

1. At this time, for U.S. ONLY, the preferred method of submitting information is the *FE Site Verification Web Site* at:

<http://qib.gehealthcare.com>

The FE Site Verification consists of three components that are available on the web from the main menu. They are:

- _ Install/deinstall product locator model and serial numbers
- _ Add/modify ship to address information
- _ Update CARES FE data for primary/secondary FEs

To obtain a copy of the FE training tutorial for using this website, a downloadable copy is available at the following:

Product Locator Support Central Page – <http://supportcentral.ge.com/15563>

Note: One “Shipping Card” is filled out and submitted when shipped (extra cards are supplied for trans-shipments between storage and distribution points), and the “Installation Card” and extra shipment cards are attached.

2. Verify the serial and model number on each rating plate matches the installation card numbers before removing the installation cards.

Note: There may be one or more shipment cards and bar code labels with the installation card. These shipment cards are used to trace the transfer of serialized units between various inventory storage and distribution points until the product reaches its final installation destination. Process just the installation card and discard any extra shipment cards and labels.

3. Enter the information from the installation cards into the GIB database for this customer.

3- COIL HARDWARE CONFIGURATIONS

The HNS coil has 5 components (or called units). The components are listed below:

- Posterior Head Component (or called Posterior Head Neck Unit)
- Face Component (or called Anterior Head Neck Unit)
- TL Component (or called TL unit)
- Horseshoe Component (or called adaptor block)
- Chest Component (or called Neck Chest unit)

The HNS coil may be configured in 9 different ways using its modular design of its components. The configurations are listed below:

- Posterior Head and Face
- Posterior Head, Face and TL
- Posterior Head and Horseshoe
- Posterior Head, Horseshoe and TL
- Posterior Head, Face and Chest
- Posterior Head, Face, Chest and TL
- Posterior Head, Horseshoe and Chest
- Posterior Head, Horseshoe, Chest and TL
- TL only

4- COIL CONFIG

4-1 Personnel Requirements

Required Persons	Preliminary Reqs	Procedure	Finalization
1	Timing Not Applicable	30 Min	Timing Not Applicable

4-2 Overview

Installation is not required for the 1.5T HNS Coil. All coil modes are present in the cfg file location for HDxt, and in the coil database for MR450/MR450w on the host software. The steps below are required only if software prior to 15M4 is on the system.

4-3 Preliminary Requirements

4-4 Procedure

1. Go to the scan desktop and select the [New Patient] button.
2. Assemble the coil in one of the 5 valid hardware setup configurations below and make sure the coil components have a snug fit with each other.

Note: Each configuration of the HNS coil must be installed for the system to recognize it. However, only 5 configurations are needed to setup the coil to get all coil configurations installed. Connecting all the coil components together may not result in the coil getting completely configured.

Illustration 1 – Setup 1; Posterior Head Component + Face Component


Illustration 2 – Setup 2; Posterior Head Component + Face Component + TL Component


Illustration 3 – Setup 3; Posterior Head Component + Horseshoe Component + TL Component


Illustration 4 – Setup 4; Chest Component Added To Setup 3



Illustration 5 – Setup 5; Posterior Head Component + Face Component + TL Component + Neck Component

- 3. Plug the connector, from the Head Neck Unit (HNU) of the coil, into the appropriate Port connector on the LPCA Port. If the coil is not configured on the system a pop-up box will appear. See Illustration 6. If the box does not appear, then proceed to step 5.

Illustration 6 – New Coil Detected Pop-up

- 4. Select **[Yes]** to configure the new coil. If the process is successful a pop-up message will be displayed as shown in Illustration 7. Continue the process of plugging in each of the five coil set-ups per illustrations 1-5.

Illustration 7 – Configuration Successful Pop-up


5. Prescribe a new exam and enter patient information.
6. Click on the coil selection button and verify that the HNS coil is present under currently connected on the scanner user interface. When the five coil Setups have been configured, you will see all the coil entries in the Coil Names window as shown in Illustrations 8a / 8b.

Illustration 8a - HNS coil modes present under currently connected on Discovery MR450 and Optima MR450w Systems

Illustration 8b - HNS coil modes present under currently connected on a HDxt system

Note: The presence of coil entries like those shown is a confirmation that the coil installed properly. If you do not see the HNS coil modes, **do not attempt to install the coil manually**, call the OnLine Center for coil configuration support.

7. After successfully configuring the coil proceed to Section 5, HNS Coil MCQA Setup.

5- HNS COIL MCQA SETUP

Refer to Service Methods for MCQA set-up instructions found in the following paths:

- Signa HDxt: Installation > Coils > HDx & HDxt Surface Coils > *1.5T Coil Vendor Manuals* > Table 6 > [1.5T Head Neck Spine \(HNS\) Coil MCQA Setup](#)
- MR450/MR450w: Coils > MR450/MR450w 1.5T Coil Vendor Manuals > Table 3 > [1.5T Head Neck Spine \(HNS\) Coil MCQA Setup](#)

6- HNS COIL SETUP FOR REVERSE RAMPED SYSTEM

6-1 Personnel Requirements

Required Persons	Preliminary Requirements	Procedure	Finalization
1	N/A	15 mins	N/A

6-2 Overview

Procedure to make the coil compatible with reverse ramped magnets for HNS coil is provided.

6-3 Tools and Test Equipment


Item	Qty	Effectivity	Part#	Manufacturer
Standard Tool Kit-Flathead screwdriver	1	N/A	N/A	N/A

6-4 Procedure


NOTICE

Observe anti-static precautions: Use wrist strap and ESD workstation (ESD mat connected to ground and place to connect ESD wrist strap).

Illustration 10 - Posterior Unit Interface board cable swap (HNS Rev 3)

The Rev 3 HNS coil can be made compatible with reverse ramped magnets for any system. To make the coil compatible the following steps are to be performed;

1. Put an ESD strap on and ground.
2. Unscrew the posterior cover and remove from coil. Retain the screws.
3. Swap cables E1-E2. Disconnect cables labeled E1 and E2. Connect the cable labeled E2 to location E1 on Interface sub-assembly. Similarly, connect the cable labeled E1 to location E2 on Interface sub-assembly.
4. Swap cables E3-E4. Disconnect cables labeled E3 and E4. Connect the cable labeled E3 to location E4 on Interface sub-assembly. Similarly, connect the cable labeled E4 to location E3 on Interface sub-assembly.
5. Swap cables E5-E6. Disconnect cables labeled E5 and E6. Connect the cable labeled E5 to location E6 on Interface sub-assembly. Similarly, connect the cable labeled E6 to location E5 on Interface sub-assembly.

Illustration 6-2 Posterior Unit- Interface board setup post-cable swap (HNS Rev 2)

The Rev 2 HNS coil can be made compatible with reverse ramped magnets for any system except DVMR system. To make the coil compatible the following steps are to be performed;

1. Put an ESD strap on and ground.
2. Unscrew the posterior cover and remove from coil. Retain the screws.
3. Swap cables E1-E2. Disconnect cables labeled E1 and E2. Connect the cable labeled E2 to location E2 on Interface sub-assembly. Similarly, connect the cable labeled E1 to location E1 on Interface sub-assembly.
4. Swap cables E3-E4. Disconnect cables labeled E3 and E4. Connect the cable labeled E3 to location E3 on Interface sub-assembly. Similarly, connect the cable labeled E4 to location E4 on Interface sub-assembly.
5. Swap cables E5-E6. Disconnect cables labeled E5 and E6. Connect the cable labeled E5 to location E5 on Interface sub-assembly. Similarly, connect the cable labeled E6 to location E6 on Interface sub-assembly.

www.gehealthcare.com

Imagination at work