GE Healthcare

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This page left intentionally blank.

Revision History

Rev	Date	Description
1.0	06/29/2012	Initial Release
2.0	09/12/2012	Release for M3

This page left intentionally blank.

Table of Contents

Chapter 1 INTRODUCTION.....15

- 1 GE 8CH Foot Ankle Coil Introduction.....15
 - 1.1 Product Identification and Shipping List.....15
 - 1.2 Warranty Statement.....17
 - 1.3 Safety.....18
 - 1.4 Compatibility.....18
 - 1.5 Related Documentation.....18
 - 1.6 Environmental Requirements.....18

Chapter 2 COIL SETUP.....23

- 1 GE 8CH Foot Ankle Coil MCQA Test Setup.....23
- 2 GE 8CH Foot Ankle Coil Shock Sensor Installation.....28

Chapter 3 TROUBLESHOOTING.....31

- 1 GE 8CH Foot Ankle Coil Troubleshooting Tips.....31
 - 1.4.1 Receiving No Signal.....31
 - 1.4.1.1 Problem:.....31
 - 1.4.1.2 Possible Solution:.....32
 - 1.4.2 Image Quality.....32
 - 1.4.2.1 Problem:.....32
 - 1.4.2.2 Possible Solution:.....32
 - 1.4.3 Artifacts.....33
 - 1.4.3.1 Problem 1:.....33
 - 1.4.3.2 Possible Solution for problem1:.....33
 - 1.4.3.3 Problem2:.....33
 - 1.4.3.4 Possible Solution for problem2:.....33
 - 1.4.4 External Cable Wear.....33
 - 1.4.4.1 Problem:.....33
 - 1.4.4.2 Possible Solution:.....33
 - 1.4.5 Mechanical Parts Wear.....35

Chapter 4 FRU AND REPLACE PROCEDURES.....37

- 1 GE 8CH Foot Ankle Coil FRUs.....37
 - 1.1 Field Replaceable Units.....37
 - 1.237
- 2 GE 8CH Foot Ankle Coil Replacements.....38
 - 2.3.1 Cable Replacement.....38
 - 2.3.2 Base Plate Replacement.....40


Chapter 1 Introduction

1 GE 8CH Foot Ankle Coil Introduction

1.1 Product Identification and Shipping List

Identify the GE 8CH Foot Ankle Coil by the coil label.

Illustration 1-1: Coil Label

The GE 8CH Foot Ankle Coil (M50002GF) is shown below.

Illustration 1-2: GE 8CH Foot Ankle Coil

This table provides a the list of items shipped with the GE 8CH Foot Ankle Coil.

Table 1-1: Contents

Description	GE Part Number	Quantity
GE 8CH Foot Ankle Coil (Base Plate, Foot Pad, Ankle Pad Knee Pad and operator manual include)	5459510	1
GE 8CH Foot Ankle Coil Service Manual	5443918-1EN	1

NOTE: Coils do not ship with phantoms. Phantoms come in a unified phantom set with the MR system.

1.2 Warranty Statement

Proper performance of this coil is guaranteed only while the coil is being used on the MR system (hardware/software level) specified at the time of purchase. Any modification of this coil may affect compatibility. Prior to modifying the coil, please contact your GE Healthcare representative to discuss coil compatibility issues. Failure to do so may void your warranty.

1.3 Safety

	Warning: Hand crushing
	Warning: Hot surface
	Warning: General sign
	Prohibited: do not loop cable
	Mandatory: refer to instructions

1.4 Compatibility

Refer to Coil Compatibility Matrix.

1.5 Related Documentation

GE 8CH Foot Ankle Coil Operator Manual, Direction 5445498.

1.6 Environmental Requirements

This equipment must be stored under the following conditions:

- Ambient temperature range of -40° to +70°C.
- Relative humidity range of 5% to 95%, noncondensing.
- Atmospheric pressure range of 525 to 1,013 hPa.

Storage Requirements: To store the coil and baseplate, a storage space greater than 38.1cm x 40.6cm x 37.5cm (Width X Length X Height), 15 in. x 16.0 in. x 14.8 in. is required.

Table 1-2: Dimensions (Width x Length x Height)

Coil:	179.5mm x 332.7mm x 314.0mm	7.0 in. x 13.0 in. x 12.2 in.
Baseplate	336mm x 406mm x 335mm	13.1 in. x 15.8 in. x 13.1 in.

Table 1-3: Weight

Coil with Cable:	2.70 kg	5.94 lb
Baseplate:	2.65 kg	5.83 lb

The images below show more detailed dimensions of the GE 8CH Foot Ankle Coil.

Illustration 1-3: Coil Dimensions


Illustration 1-4: Baseplate Dimensions



Illustration 1-5: Dimensions

This page left intentionally blank.

Chapter 2 Coil Setup

1 GE 8CH Foot Ankle Coil MCQA Test Setup

1.1 Personnel Requirements

Personnel Requirements	Preliminary Reqs	Procedure	Finalization
1	0 min	15 mins	0 min

1.2 Overview

NOTE: Coils do not ship with phantoms. Phantoms come in a unified phantom set with the MR system.

1.3 Preliminary Requirements

1.3.1 Tools and Test Equipment

Item	Qty	Effectivity	Part#	Manufacturer
Small Cylindrical Unified Phantom	1	-	5342680	-
Cubical Unified Phantom	1	-	5342681	-

1.3.2 Replacement Parts

Item	Qty	Effectivity	Part#	Manufacturer
See GE 8CH Foot Ankle Coil FRUs	-	-	-	-

1.3.3 Safety

WARNING

ALL RF COIL-RELATED TESTS MUST BE RUN ON A SYSTEM THAT IS WELL CALIBRATED AND PASSES ALL SYSTEM TESTS. THE SYSTEM SHOULD HAVE PASSED INSTALL IN SPEC (IIS), ESPECIALLY WHITE PIXEL, CORRELATED NOISE, AND MULTI-COIL-RECEIVE (MCR) TOOL.

1.3.4 Required Conditions

Condition	Reference	Effectivity
The following coil configuration name must be installed to run this tool: 8CH Foot Ankle . Illustration 2-1: Coil DB Explorer 	-	-

1.4 Procedure

1. Place the Foot Ankle Coil baseplate on the table as shown below.

Illustration 2-2: Foot Ankle Coil Baseplate on Table

2. Place the Small Cylindrical Unified Phantom horizontally on the baseplate, place the Cubical Unified Phantom vertically, and fix the phantom to the baseplate using the Velcro as shown below.

Illustration 2-3: Cylindrical Unified Phantom on Baseplate

3. Cover the baseplate with the coil as shown below.

Illustration 2-4: Coil Over Baseplate

4. Connect the Foot Ankle Coil to Port A of the system LPCA, and landmark the coil at the crossmark on the coil handle.
 - a. Advance to scan.
 - b. Perform the MCQA Test.

Illustration 2-5: Crossmark on Coil Handle

5. After the test is completed, carefully remove the Cubical Phantom out of the baseplate, then the Cylindrical Phantom to avoid damage to either phantom.

1.5 Finalization

No finalization steps.

2 GE 8CH Foot Ankle Coil Shock Sensor Installation

2.1 Personnel Requirements

Personnel Requirements	Preliminary Reqs	Procedure	Finalization
1	Not Applicable	5 mins	5 mins

2.2 Preliminary Requirements


2.2.1 Tools and Test Equipment

Item	Qty	Effectivity	Part#	Manufacturer
Shock Sensor Sticker	1 box of 50	-	5337009(100G)	ShockWatch
Isopropyl Alcohol	1	-	46-183000P164	-

2.3 Procedure

1. Place the sensor per [Illustration 2-6](#)

Illustration 2-6: GE 8CH Foot Ankle Coil shock sensor placement

2. Check if the shock sensor has been triggered (center tube turns red – see right side of [Illustration 2-7](#))
 - a. If the shock sensor has NOT been triggered (center tube is white – see left side of [Illustration 2-7](#)), continue to Finalization.
 - b. If the shock sensor has been triggered (center tube is red – see left side of [Illustration 2-7](#)) or contaminated in any way, continue to [Step 3](#).

NOTE: A shock sensor may be triggered and the tube is not completely red. Any red or pink color in the center tube is considered triggered. See [Illustration 2-8](#) for an example.

Illustration 2-7: Shock Sensor

Illustration 2-8: Triggered Shock Sensor

3. Remove the triggered shock sensor from the coil and dispose of it using local procedures.
4. Verbally educate the customer on safe use of coils and implications to warranty on mishandling of coils.
5. Re-qualify the coil by running the Multi-Coil Quality Assurance Tool test.
 - a. If the coil passes the MCQA test, apply a new shock sensor sticker (follow [Step 1](#)) and continue to Finalization.
 - b. If the coil fails the MCQA test, continue to [Step 6](#)

6. Replace the coil following replacement procedures. Charge the customer for repair/ replacement of the damaged coil per local procedures. Verbally educate the customer on safe use of coils and implications to warranty on mishandling of coils.

NOTE: REMOVE THE SHOCK SENSOR BEFORE RETURNING THE DAMAGED COIL!

2.4 Finalization

1. Make sure that you have completed Section 2.3 for the coil at the customer site. Return the coil for use to the customer.
2. Re-use the original packaging of the coil shock sensors for transportation. Order a new box of shock sensor stickers through SSP purchasing if you are running low. Order a couple weeks ahead of your next PM to allow for delivery.

Chapter 3 Troubleshooting

1 GE 8CH Foot Ankle Coil Troubleshooting Tips

1.1 Personnel Requirements

Personnel Requirements	Preliminary Reqs	Procedure	Finalization
1	0 min	60 mins	0 min

1.2 Overview

The following tips can be used to troubleshoot common problems with the GE 8CH Foot Ankle Coil, catalog M50002GF.

1.3 Preliminary Requirements

1.3.1 Tools and Test Equipment

Item	Qty	Effectivity	Part#	Manufacturer
Digital Volt Meter	1	-	-	-

1.3.2 Required Conditions

Condition	Reference	Effectivity
<p>The following coil configuration name must be installed to run this tool: 8CH Foot Ankle</p> <p>Illustration 3-1: Coil DB Explorer</p> 	-	-

1.4 Procedure

1.4.1 Receiving No Signal

1.4.1.1 Problem:

You are unable to pre-scan or are scanning and yet receiving no signal.

1.4.1.2 Possible Solution:

- Verify that you have made the appropriate system coil selection. Refer to the Operators Manual for additional information.
- Verify the green light above Port A is illuminated. This indicates the coil is properly plugged into the system.
- Verify that the landmark is correct and that the cradle has not unlatched.
- Verify that the scan locations and any FOV offsets are correct.
- Perform a continuity check on the external cable (below). The continuity check must be performed by a GE authorized Service Engineer.
- Verify that the coil is properly positioned.
- If you still cannot get a signal, try to scan (transmit and receive) with the body coil. For this test, be sure to remove the imaging coil from the magnet bore before you scan with the body coil. If you still receive no signal, the problem probably lies with your MR system. If the scan completes successfully, there is probably a problem with the Invivo coil. Contact GE for further assistance. If you are unable to scan with the substitute coil, there may be a system problem related to this particular coil type.

1.4.2 Image Quality

1.4.2.1 Problem:

Image quality is not what you expect.

1.4.2.2 Possible Solution:

- Review the selected protocol. If you are performing the Periodic Quality Assurance, be sure your protocol is identical to the protocol provided in your Systems Recommended Protocols. If you are performing diagnostic images, you may need to increase NEX or FOV.
- Perform a continuity check on the external cable (below). The continuity check must be performed by a GE authorized Service Engineer.
- Verify that there are no loops in the cables.
- Verify that there are no metal or ferromagnetic objects close to the coil, patient or magnet (i.e., safety pin, hair pin).
- Verify that the coil is properly positioned.
- Verify that your center frequency is within the frequency adjustment range for your system.
- Verify that the R1, R2 and TG values from the pre-scan are within normally expected ranges.
- Cleaning the electrical contacts that connect the upper and lower halves with Isopropyl Alcohol and a cotton swab has been demonstrated to improve the coil's SNR on occasion. If you have not done so already, perform a MCQA (Multi Coil Quality Assurance) Test. If the values you

obtain do not fall within normal operating parameters, investigate this further by performing a phantom scan with the body coil. For this test, be sure to remove the imaging coil from the magnet bore before you scan with the body coil. If you still have the same problems, there is probably an MR system problem. If the body coil scan is satisfactory, acquire a scan using both another coil of the same type plugged into the same system port (port A). If the image quality is visibly improved, there may be a problem with the coil. Contact GE for further assistance. If the image quality still suffers, there may be a system problem related to imaging with this type of coil.

1.4.3 Artifacts

1.4.3.1 Problem 1:

There is a black line or signal void on the image.

1.4.3.2 Possible Solution for problem1:

Verify that there is no metal present in the area being scanned in or on the patient.

1.4.3.3 Problem2:

Some or all of the images appear shaded or exhibit uneven signal or banding.

1.4.3.4 Possible Solution for problem2:

Confirm that no metallic objects are located nearby, outside the FOV. This is especially important on images utilizing Fat Saturation. If Fat Saturation is being used, verify that the Center Frequency fine adjustment has been optimized.

1.4.4 External Cable Wear

1.4.4.1 Problem:

The system cannot recognize the coil or scan with the coil attached.

1.4.4.2 Possible Solution:

Perform the cable continuity check. As follows:

- RF Signal Channels Check:** The center pins of the System Connector receive coaxial connectors are extremely fragile. Do not attempt to make any measurements at these center pins.


NOTE: With the coil cable disconnected, use a multimeter to assure an open circuit condition exists between the designated pin and GND for each signal channel shown in Table 3–1.

Table 3-1: Coil Cable Connector Short Circuit and Open Circuit Check

Signal Channel	Coil Cable Connector Pin (Illustration 3–1)	System Connector (Illustration 3–2)
RF 1	CH1	C1
RF 2	CH2	C2

RF 3	CH3	C3
RF 4	CH4	C4
RF 5	CH5	D1
RF 6	CH6	D2
RF 7	CH7	D3
RF 8	CH8	D4
GND		K1-K8, M6, MC Shield

Illustration 3-2: Coil Cable Connector

If one or more of the readings indicate a short, replace the cable. If all of the above readings indicate an open circuit condition, proceed to the PIN Diode Control Channels Check.


2. **MC bias line (DC line) Check:** With the coil cable disconnected, use a multimeter to determine continuity between the designated pins for each signal channel below. Also verify that the control pins are not shorted to GND.

Table 3-2: Coil Cable Connector Continuity Check

Signal Channel	Coil Cable Connector (Illustration 1)	System Connector (Illustration 2)
MC bias 1	YEL	J1
MC bias 2	GRN	J2
MC bias 3	PURPLE	J3
MC bias 4	BRN	J4

MC bias 5	YEL/WHT	J5
MC bias 6	GRY	J6
MC bias 7	GRN/WHT	J7
MC bias 8	BLU	J8
GND		K1-K8, M6, MC Shield

Illustration 3-3: System Connector

1.4.5 Mechanical Parts Wear

Inspect the coil for visible signs of wear. If the coil housings cannot seat with each other or if the latch mechanism cannot work properly, contact Service.

1.5 Finalization

No finalization steps.

This page left intentionally blank.

Chapter 4 FRU and Replace procedures

1 GE 8CH Foot Ankle Coil FRUs

The following tables list the Field Replaceable Units (FRUs) and other replaceable parts for the GE 8CH Foot Ankle Coil.

1.1 Field Replaceable Units

Table 4-1: Field Replaceable Units

Description	GE Part Number
GE 8ch Foot Ankle Coil Assembly	5436920
GE 8ch Foot Ankle Coil Cable	5436148
GE 8ch Foot Ankle Coil Base Plate	5458752
Small Cylindrical Unified Phantom	5342680
Cubical Unified Phantom	5342681
Strap for GE 8CH Foot Ankle Coil	5439452
Ankle Pad for GE 8CH Foot Ankle coil	5448720
Foot Pad for GE 8CH Foot Ankle Coil	5448720-2
Knee Pad for GE 8CH Foot Ankle Coil	5448720-3

1.2

2 GE 8CH Foot Ankle Coil Replacements

2.1 Personnel Requirements

Personnel Requirements	Preliminary Reqs	Procedure	Finalization
1	0 min	30 mins	0 min

2.2 Preliminary Requirements

2.2.1 Replacement Parts

Item	Qty	Effectivity	Part#	Manufacturer
Refer to GE 8CH Foot Ankle Coil Field Replaceable Units		-	-	-

2.2.2 Safety

NOTICE


This procedure must be performed by an authorized GE Service Engineer only.

2.3 Procedure

2.3.1 Cable Replacement

1. Proceed as follows for the Cable Replacement procedure:
 - a. Remove the six Coil Connector mounting screws as detailed below.

Illustration 4-1: Removing Screws on Coil Connector

- b. Disconnect the cable by grasping the coil connector housing, and pulling the cable away from the coil connector in a downward direction.

Illustration 4-2: Disconnecting Cable from Coil Connector



Illustration 4-3: Coil Cable Connector

- c. Perform MCQA after replacement

2.3.2 Base Plate Replacement

Do the following steps to check the new base plate.

1. Exercise the latches on the baseplate by moving the blue latch handle back and forth. (The foot support and coil is capable of a 15 degree plantar flexion.)
 - Latches must move smoothly and lock positively into place.
 - Latches must not loosen or open during the course of an examination
2. Move the coil from side to side on the baseplate, making sure the coil moves smoothly.
 - a. Swivel the coil on the baseplate. The coil must move smoothly within a 10 to 15 degree rotation.
 - b. Lock the coil into position on the baseplate. The coil must remain in place while locked down, and resume smooth movement when unlocked.

2.4 Finalization

No finalization steps.

© 2012 General Electric Company.
General Electric Company, doing business as GE Healthcare.

www.gehealthcare.com