

LightSpeed 4.X (LS 16) GOC3 (GRE) Console Information Sheet

- Site ID #: _____ Console Serial Number _____
- DAS Type: MDAS16
- Tube Type (15): Performix
- PDU Type: CT COMPACT
- KV Values: 80 100 120 140
- Collimator/Filter Type: Helios G1
- Scan Recon Hardware Configuration: GLOBAL RECON ENGINE
- GANTRY TYPE: H2
- Target Noise Index Table: Table 1 Table 2 Table 3
- Dose Info Display: On On Without Total DLP Off
- Modified in Room Start: Off
- HIPAA Present: On Off
- Default Archive Device: DICOM
- Suite Name: _____
- Host Name: _____
- IP Address: ____ . ____ . ____ . ____
- Netmask: ____ . ____ . ____ . ____
- Broadcast Address: ____ . ____ . ____ . ____
- Default Gateway: ____ . ____ . ____ . ____
- Printer: _____
- Printer IP: ____ . ____ . ____ . ____
- Internal Option – Internal Subnet: "this item may not appear" do not modify contents
- Advanced Option: this item is checked
- Use NIS? Enter Domain Name: _____
IP Address of Internal Server: ____ . ____ . ____ . ____
- ConnectPRO Option Info:
CTTITLE: _____
PORTNUM: _____
AETITLE: _____
IPADDR: ____ . ____ . ____ . ____