

## LightSpeed 5.X (LS Pro<sup>16</sup> (80kW)) GOC3 (GRE) Console Information Sheet

- Site ID #: \_\_\_\_\_ Console Serial Number \_\_\_\_\_
- DAS Type: GDAS16 or MDAS16
- Tube Type (16): OEM
- PDU Type: NG PDU
- KV Values: 80 100 120 140
- Collimator/Filter Type: Helios G1
- Scan Recon Hardware Configuration: GLOBAL RECON ENGINE
- Gantry Type: HP
- Target Noise Index Table: Table 1 Table 2 Table 3
- Dose Info Display: On On Without Total DLP Off
- Modified in Room Start: Off
- HIPAA Present: On Off
- Default Archive Device: DICOM
- Suite Name: \_\_\_\_\_
- Host Name: \_\_\_\_\_
- IP Address: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_
- Netmask: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_
- Broadcast Address: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_
- Default Gateway: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_
- Printer: \_\_\_\_\_
- Printer IP: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_
- Internal Option – Internal Subnet: "this item may not appear" do not modify contents
- Advanced Option: this item is checked
- Use NIS? Enter Domain Name: \_\_\_\_\_  
IP Address of Internal Server: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_
- ConnectPRO Option Info:  
CTTITLE: \_\_\_\_\_  
PORTNUM: \_\_\_\_\_  
AETITLE: \_\_\_\_\_  
IPADDR: \_\_\_\_ . \_\_\_\_ . \_\_\_\_ . \_\_\_\_