


State: RELEASE - Document is released and under formal Change Control. Changes are subject to the ECR/ECO Process. See the GEHC Myworkshop System to determine the status of this document.

GOC6 CONSOLE OUTLET BOX SCHEMATIC

PART NUMBER 5272833

REVISION	ECR/ECO	DATE	SIGN	OFF
4	N/A	4/20/2008	M. ROSE	
5	2066102	5/02/2008	M. ROSE	
6	N/A	5/29/2010	M. ROSE	



120/208Y VAC OR 120/240 VAC

NEMA L14-20R

TO GANTRY INTERFACE CABLE (J20)

SYSTEM SIGNAL

SATELLITE SIGNAL

3-POLE CONTACTOR CONT_1

POWER SWITCH PASS-THROUGH SW1

PERFORMANCE EMC FILTER FL_1

PERFORMANCE EMC FILTER FL_2

GROUND BAR GB1

CIRCUIT BREAKER CB1

CIRCUIT BREAKER CB2

CIRCUIT BREAKER CB3

ACC 1 GROUND

ACC 2 GROUND

ESSENTIAL CLASSIFICATION: MINOR

GE Healthcare Company Restricted

Rev 6

Dwg No 5272833SCH

Block 1 of 1

Sheet 1 of 1

Mon May 03 17:16:46 2010