

ULTRA 16 APPLICATIONS FET MODE SELECTIONS - 16 SLICE MODE

A12	A11	A10	A9	A8	A7	A6	A5	A4	A3	A2	A1	B1	B2	B3	B4	B5	B6	B7	B8	B9	B10	B11	B12	Slice mode	FET mode	FET 1 LED	FET 2 LED	FET 3 LED	FET 4 LED	FET 5 LED	FET 6 LED
2	1									4	3	3	4									1	2	Z-Chan	49	OFF	ON	ON	ON	OFF	OFF
8	7	6	5	4	4	3	3	2	2	1	1	1	1	2	2	3	3	4	4	5	6	7	8	16x1.25	54	ON	OFF	OFF	ON	OFF	OFF
2	1									4	3	3	4									1	2	Z-Chan	49	OFF	ON	ON	ON	OFF	OFF
4	4	3	3	2	2	2	2	1	1	1	1	1	1	1	1	2	2	2	2	3	3	4	4	8x2.5	1	OFF	ON	ON	ON	ON	ON
3	2	1								4	4										1	2	3	Z-Chan	48	ON	ON	ON	ON	OFF	OFF
		2	2	2	2	1	1	1	1	1	1	1	1	1	1	2	2	2	2					4x3.75	7	OFF	OFF	OFF	ON	ON	ON
			2	1	1													1	1	2				Z-Chan	39	OFF	OFF	OFF	ON	ON	OFF
				8	7	6	5	4	3	2	1	1	2	3	4	5	6	7	8					16x0.625	53	OFF	ON	OFF	ON	OFF	OFF
			2	1	1													1	1	2				Z-Chan	39	OFF	OFF	OFF	ON	ON	OFF
				4	4	3	3	2	2	1	1	1	1	2	2	3	3	4	4					8x1.25	54	ON	OFF	OFF	ON	OFF	OFF
4	4	3	3	2	2	2	2	1	1	1	1	1	1	1	1	2	2	2	2	3	3	4	4	Z-Chan	1	OFF	ON	ON	ON	ON	ON
								4	3	2	1	1	2	3	4									8x0.625	3	OFF	OFF	ON	ON	ON	ON