


REV 0

## SECTION COIL – 9INCH GP COIL



**P9424LA 9 INCH GP ASSY      2269507-2      REV 000**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2269510-2	1	9' GP Coil Dom	1	P9424LC 9 INCH GP COIL FOR JAPAN
6	2210268	1	9inch ITEM2	1	P9364RR 9INCH GP SV ITEM2
7	2198505	2	Coil Holder	1	P9369QA COIL HOLDER
8	2197185	DIA	PAD	1	P9369PN W-PAD-2
9	2198911	DIA	HEAD STRAP	2	P9369QW HEAD STRAP
10	2198912	DIA	CUSHION	1	P9369QX HEAD CUP CUSHION

**P9424LB 9 INCH GP ASSY      2269508-2      REV 000**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2269512-2	1	9' GP Coil Exp	1	P9424LD 9 INCH GP COIL FOR OUTSIDE JAPAN
6	2210268	1	9inch ITEM2	1	P9364RR 9INCH GP SV ITEM2
7	2198505	2	Coil Holder	1	P9369QA COIL HOLDER
8	2197185	DIA	PAD	1	P9369PN W-PAD-2
9	2198911	DIA	HEAD STRAP	2	P9369QW HEAD STRAP
10	2198912	DIA	CUSHION	1	P9369QX HEAD CUP CUSHION