


REV 1

SECTION COIL – BODY COIL

P9424BA BODY COIL ASSY 2276942 REV 0

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2283072	1	Body Coil	1	P9424DA COIL ASSY
2	2283063	1	SV Kit	1	P9424DM TUNING TOOL
3	2285783	1	JIG	1	P9424BG B COIL JIG

P9424DA COIL ASSY 2283072 REV 001



P9424DA COIL ASSY 2283072 REV 001


Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2283073	N	Coil Top	1	P9424DB TOP ASSY
2	2283075	N	Coil Bottom	1	P9424DC BOTTOM ASSY

REV 1

P9424DB TOP ASSY 2283073 REV 002

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2283102	N	P9424DK FTOP ASSY	1	P9424DK FTOP ASSY
4	2285779	1	AIR FILTER	0	P9424BC AIR FILTER
5	2285780	1	SEAL	0	P9424BD SEAL A
6	2285781	1	SEAL	0	P9424BE SEAL B
7	2285782	1	SEAL	0	P9424BF SEAL C
8	2285865	2	COVER	0	P9424BH CONDENSER COVER
9	2285866	1	COVER	0	P9424BJ FO BOARD COVER

P9424DM TUNING TOOL 2283063 REV 002



F0 Tune Board

2283173 consists of 4 series of
F0 Tune Boards

P9424DM TUNING TOOL 2283063 REV 003

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	Z9817BK	2	SHORT CONNECTER	3	SHORT CONNECTER
3	2283173	N	FO BRD ASSY	8	P9424CU FO BRD ASSY
4	Z9807DB	2	BNC-003	4	BNC-003
5	Z9807QA	2	PD003	2	PD003

REV 1

Revision History

Rev	Date	Author	Primary reasons for change
0	May 23, 2001	M.Ikeda	Initial release
1	Jan. 31, 2002	M.Ikeda	Corrected quantity of FO BRD ASSY