


REV 0

## SECTION COIL – HEAD COIL



**P9424EA HEAD COIL ASSY      2276966      REV 001**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2280068	1	HEAD COIL	1	P9424EC HEAD COIL
2	2284805	1	MIRROR2	1	P9424FD MIRROR 2
3	2284794	1	HEAD REST	1	P9424ED HEAD REST ASSY
4	P9140TE	DIA	CHIN BAND	2	CHIN BAND
5	P9140TD	DIA	HEAD BAND	2	HEAD BAND
6	2198920	DIA	COIL PAD	1	P9300UU COIL PAD
7	2198833	DIA	HEAD CUSHION2	1	P9360JP HEAD CUSHION 2

**P9424EC HEAD COIL      2280068      REV 001**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
12	U0432AE	2	CATCH	1	CATCH
16	2284524	N	HEAD COIL IN ASSY	1	P9424EU HEAD COIL IN ASSY

**P9424EU HEAD COIL IN ASSY      2284524      REV 003**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
107	U0201AC	2	CLAMP	4	CLAMP
202	2284523	1	BOX ASSY	1	P9424ET HEAD BOX ASSY

**P9424ED HEAD REST ASSY      2284794      REV 002**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
33	U0092AA	1	SCREW	2	SCREW