

REV 0

## SECTION MENC – MAGNET ENCLOSURE

**P9423CA ENCLOSURE ASSY      2282871      REV 008**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2284788	N	FLOOR ASM	1	P9423BA FLOOR ASSY
5	2284792	N	POST L ASM	1	P9423BE POST L ASSY
6	2284793	N	POST R ASM	1	P9423BF POST R ASSY
33	P9123PC	1	CAP	18	CAP /
34	2284919	1	CAP	6	U0255AE CAP
35	2284920	1	CAP	4	U0256AE CAP
36	2151637	1	CAP	8	RUBBER CAP
38	2290637	1	CAP	4	P9423JW CAP

**P9423BA FLOOR ASSY      2284788      REV 002**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
2	2283179	N	FRONT TS ASM	1	P9423VA FRONT TS ASSY

**P9423VA FRONT TS ASSY      2283179      REV 003**

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2284625	1	FR TS	4	P9423VC F TOUCH SENSOR
2	2284627	1	FR TS COVER	1	P9423VG F TS COVER

**P9423BF POST R ASSY      2284793      REV 001**


Item	Part Number	FRU	Name	Quantity	Description (Remarks)
2	2283180	1	SIDE TS ASM	1	P9423VB SIDE TS ASSY

REV 0

P9423EA FRAME ASSY

2282872

REV 007



P9423EA FRAME ASSY

2282872

REV 007

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
2	2284778	1	LAMP 12V20W	2	P9423SY LAMP 12V20W
158	2291196	N	LIGHT BRKT	2	P9423BU LIGHT BRKT ASSY

P9423BU LIGHT BRKT ASSY

2291196

REV 000

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2284777	1	CABLE	1	P9423SX SOCKET CBL ASSY