

## SECTION SCBL – SYSTEM CABLES

### SYSTEM CABLES

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
1	2283545	N	CABLE ASM	1	P9420SE MFO FIXED SITE DATA CABLES
2	2198993	1	CABLE (#818)	1	(818) MR2-A11-J31 TO OW1-
3	46-301125G807	1	CABLE (#701)	1	(701T) MR2-A11-J15 TO RF
4	46-258431G810	1	CABLE (#702)	1	(702M) MR2-A11-J18 TO MR1
5	46-258431G809	1	CABLE (#703)	1	(703M) MR2-A11-J23 TO MR1
6	46-301125G98	1	CABLE (#229)	1	(229T) MR2-A11-J1 TO MR1-
7	46-301125G100	1	CABLE (#231)	1	(231T) MR2-A11-J2 TO PP1
8	2263200	1	CABLE (#837)	1	(837) OW1-A17-J4 TO PP1-A
9	2263200-2	1	CABLE (#838)	1	(838) OW1-A17-J3 TO PP1-J
10	2263200-3	1	CABLE (#841)	1	(841) PT1-J5 TO PP1-J51
11	2263200-14	1	CABLE (#906)	1	(906) PT1-J11 TO PP1-J47
12	2263200-15	1	CABLE (#910)	1	(910)OW1-A17-J2 TO OW1-A
13	2263200-16	1	CABLE (#911)	1	(911)OW1-A17-J2 TO OW1-A
14	2263200-18	1	CABLE (#921)	1	(921)PT1-J8 TO PP1-J11
15	2263200-26	1	CABLE (#929)	1	(929)MR2-A11-J5 TO PP1-J
16	2263200-23	1	CABLE (#926)	1	(926) PP1-A15-J78 TO PT1-
17	2263200-24	1	CABLE (#927)	1	(927) PP1-A16-J77 TO PT1-
18	2263200-27	1	CABLE (#934)	1	(934) PP1-J8 TO MG3-A3-A5
19	46-328000G902	1	CABLE (#768)	1	(768TR) PP1-J20 TO MR1-A7
20	46-328000G905	1	CABLE (#771)	1	(771TR) PP1-J47 TO MR1-A7
21	46-328000G904	1	CABLE (#770)	1	(770) PP1-J11 TO MR1-A7-J
22	46-328000G903	1	CABLE (#769)	1	(769) PP1-J15 TO MR1-A7-J
23	46-258431G908	1	CABLE (#774)	1	(774M) MR2-A11-J25 TO MR1
24	46-258431G808	1	CABLE (#705)	1	(705M) MR2-a11-J14 TO PD1
25	46-328000G906	1	CABLE (#772)	1	(772) PP1-J14 TO MR1-A7-J
26	46-328000G907	1	CABLE (#773)	1	(773) PP1-A17-J1 TO MR1-A
27	46-317359G931	1	CABLE (#716)	1	(716TR) MG2-A11-A1-J2 TO
28	46-328000G915	1	CABLE (#775)	1	(775) MR1-A7-J42 TO PP1-A
29	46-301125G764	1	CABLE (#487)	1	(487T) MR1-A7-J20 TO TP PP1-
30	46-301125G765	1	CABLE (#488)	1	(488T) MR2-A11-J8 TO PP1-
31	46-301125G766	1	CABLE (#489)	1	(489T) MR2-A11-J9 TO PP1-
32	46-301125G767	1	CABLE (#490)	1	(490T) MR2-A11-J10 TO PP1
33	46-301125G831	1	CABLE (#726)	1	(726T) MR1-A7-J18 TO PP1-A
34	46-328000G937	1	CABLE (#788)	1	(788) PP1-J12 TO OW1-A1-A2
35	46-328000G938	1	CABLE (#789)	1	(789) MR2-A11-J26 TO OW1-A
36	46-328000G939	1	CABLE (#791)	1	(791) MR2-A11-J30 TO OW1-
37	46-271692G601	1	CABLE (#296)	1	(296T) PP1-J18 TO FACILITY DISCONN
38	46-271693G1	1	CABLE (#297)	1	RUN 297, EMERGENCY OFF SW TO PP1-J18
39	2283548	N	CABLE ASM	1	P9420SC MFO PDU CABLES
40	46-271364G517	1	CABLE (#030)	1	(030M) MR2-A12 TO PD1-A2

REV 0

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
41	46-255947G607	1	CABLE (#037)	1	(037) MR2-A12-GND TO PD1-
42	2103427	1	CABLE (#040)	1	(040) MS1 TO RF COM GROUN
43	46-255947G9	1	CABLE (#044)	1	(044) MR1-A7-GND TO PP1-G
44	46-307832P4	1	CABLE (#047)	1	(047) OW1A2A5J1 TO PD1-A2
45	46-255947G11	1	CABLE (#048)	1	(048) OW1-GND TO PD1-GND
46	46-317249G39	N	CABLES F.O.	1	HFO FIXED SITE FIBEROPTIC
47	46-317076G602	2	FO CABLE ASM	1	(708T) MR2-A11-J16 TO MR1-A7-TX
48	46-307245P2	2	1/2" CONDUIT	0.2	100 FOOT ROLL, HEYCO #8431
49	46-301125G923	2	FO CABLE ASM	1	(708T-1) MR2-A23-A2-J5 TO MR1-A15-A1-
50	46-307564P12	N	FO CABLE ASM	1	41'0" +/- 6"
51	46-208781P2	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
52	46-221856P13	2	TAPE	AR	TAPE, ELECTRICAL, PVC, CO
53	46-221856P14	2	TAPE	AR	TAPE, ELECTRICAL, PVC, CO
54	46-221856P15	2	TAPE	AR	TAPE, ELECTRICAL, PVC, CO
55	46-230653P2	2	TUBING OLD	20	1/4" DIA. SPIRALLY CUT POLYETHYLENE
56	46-221856P13	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, COLOR-BLUE
57	46-221856P14	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, COLOR-GREEN
58	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, COLOR-ORANGE
59	46-317085G603	2	CABLE	1	(710T)MR2-A11-J17 TO MR3-
60	46-307245P3	2	3/4" CONDUIT	0.2	100 FOOT ROLL, HEYCO #8431
61	46-301125G952	2	CABLE	1	(710T-1)MR2-A15-A23-J13 O
62	46-307584P5	N	FO CABLE ASM	1	36'0" +/- 6"
63	46-208781P2	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
64	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
65	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
66	46-301125G953	2	CABLE	1	(710T-2)MR2-A15-A23-J15 O
67	46-307584P5	N	FO CABLE ASM	1	36'0" +/- 6"
68	46-208781P2	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
69	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
70	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
71	46-301125G954	2	CABLE	1	(710T-3)MR2-A15-A23-J17 O
72	46-307584P5	N	FO CABLE ASM	1	36'0" +/- 6"
73	46-208781P2	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
74	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
75	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
76	46-301125G955	2	CABLE	1	(710T-4)MR2-A15-A23-J14 O
77	46-307584P5	N	LABEL OLD	1	36'0" +/- 6"
78	46-208781P2	2	TAPE OLD	2	WRITE-ON LABEL FOR MARKIN
79	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
80	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
81	46-301125G956	2	CABLE	1	(710T-5)MR3-GIP-MDS-JS TO
82	46-307584P5	N	LABEL OLD	1	36'0" +/- 6"
83	46-208781P2	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
84	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
85	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
86	46-208781P4	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
87	46-230653P3	2	TUBING OLD	15	SPIRALLY CUT POLYETHYLENE
88	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
89	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
90	2126180	2	CONNECTOR KT	1	FIBER-OPTIC CONDUIT PANEL

REV 0

Item	Part Number	FRU	Name	Quantity	Description (Remarks)
91	46-307244P3	2	CONNECTOR	4	PANEL MOUNT FOR 1/2" F.O.
92	46-307244P4	2	CONNECTOR	3	PANEL MOUNT FOR 3/4" F.O.
93	508A865P5	N	PACKING	1	PACKING
94	2123149	N	DIRECTION	1	FIBER OPTIC CONNECTOR INS
95	46-220191P1	2	CLAMP OLD	2	CLAMP FOR FLAT CABLE
96	46-230653P2	2	TUBING OLD	3	1/4" DIA. SPIRALLY CUT PO
97	46-208758P2	2	TY-RAP OLD	9	5.5" X .140" SELF-LOCKING
98	46-328078G5	2	FO CABLE ASM	1	(711V/712V) MR2A11-J13 TO
99	46-328073P9	2	FO CABLE ASM	1	(712V-1) MR2-A15-A23-J9,1
100	46-328073P10	2	FO CABLE ASM	1	(711V-1) MR2-A15-A23-J6,7
101	46-328074P5	2	FO CABLE ASM	1	(711V-2) MR2-A15-A23-J8 T
102	46-230653P2	2	TUBING OLD	5.375	1/4" DIA. SPIRALLY CUT PO
103	46-307245P3	2	3/4" CONDUIT	0.4	100 FOOT ROLL, HEYCO #843
104	46-208781P4	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
105	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
106	46-221856P15	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
107	46-328079G9	2	FO CABLE ASM	1	(711M/712M) PP1-J91 TO MG
108	46-328076P5	2	FO CABLE ASM	1	(711M-1) MG2-A33-J12 TO P
109	46-328077P5	2	FO CABLE ASM	1	(711M-2)MG2-A33-J13 TO PP
110	46-307245P3	2	3/4" CONDUIT	0.77	100 FOOT ROLL, HEYCO #843
111	46-208781P4	2	LABEL OLD	2	WRITE-ON LABEL FOR MARKIN
112	46-328075P9	2	FO CABLE ASM	1	PP1-A17-A1-J1,2 TO MG2-A1
113	46-230653P2	2	TUBING OLD	19	1/4" DIA. SPIRALLY CUT PO
114	46-221856P11	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
115	46-221856P13	2	TAPE OLD	AR	TAPE, ELECTRICAL, PVC, CO
116	2241052	2	FO CABLE ASM	1	(836) OW1A15-BIT3 TO MR2A
117	46-317185G1	2	SPEAKER BOX	1	PATIENT COMMUNICATION BOX
118	46-317177G1	2	HARNESS	1	MG3 A2 J8 TO MG2 A2 A1 J3
119	46-317184G1	2	HARNESS	1	MG3 A2 J7 TO MG2 A2 A1 J1
120	46-317184G2	2	HARNESS	1	MG3 A2 J9 TO MG2 A2 A1 J4
121	2288234	1	CABLE		P9420TL CABLE ASSY
122	2298512	1	CORE		
123	U0028MC	1	CORE		