

DATA SHEETS

TABLE OF CONTENTS

<u>SECTION</u>	<u>PAGE</u>
SECTION 1 - PAC-II LEAKAGE CURRENT	DS-2
1-1 PAC-II LEAKAGE CURRENT MEASUREMENTS	DS-2
SECTION 2 - SHIM	DS-3
2-1 MAGNET SHIM VALUES.....	DS-3
SECTION 3 - GRADCAL DATA	DS-4
3-1 GRADIENT CALIBRATION VALUES.....	DS-4

SECTION 1 - PAC-II LEAKAGE CURRENT

DATA TABLE 1-1
 PAC-II LEAKAGE CURRENT MEASUREMENTS

LEAD SELECT POSITION	METER READING NORMAL POSITION (microamps)	METER READING REVERSE POSITION (microamps)
Patient Leakage Current (Lead to Ground)		
Right Leg (RL)		
Right Arm (RA)		
Left Arm (LA)		
Left Leg (LL)		
Patient Auxillary Current (Lead to Lead)		
RL		
RA		
LA		
LL		
M.A.P. (ISOLATION Tests)		
ALL LEADS M.A.P. (ISO) Test		

SECTION 2 - SHIM

DATA TABLE 2-1
 MAGNET SHIM DATA

SITE _____							
SPHERICAL HARMONICS	SIGNA COORDINATE SYSTEM - SPHERICAL HARMONIC ERROR IN HERTZ						
Z1							
Z2							
Z3							
Z4							
Z5							
Z6							
X							
Y							
ZX							
ZY							
X2-Y2							
XY							
Z2X							
Z2Y							
Z(X2-Y2)							
ZXY							
X3							
Y3							
HOMOGENEITY STD DEV (HZ)							

SECTION 3 - GRADCAL DATA

DATA TABLE 3-1
 GRADIENT CALIBRATION VALUES

DATE _____ NAME _____

		INITIAL VALUE	NEXT VALUE	NEXT VALUE	FINAL VALUE	FINAL DISTANCE

DATE _____ NAME _____

		INITIAL VALUE	NEXT VALUE	NEXT VALUE	FINAL VALUE	FINAL DISTANCE

DATE _____ NAME _____

		INITIAL VALUE	NEXT VALUE	NEXT VALUE	FINAL VALUE	FINAL DISTANCE