

Measurements for Ground Leakage Current Test

Record the measurements taken in Ground Leakage Current Test. Mark any fields that are not applicable to the current location as “not applicable” or “N/A”.

Safety Analyzer Information		
Brand of Safety Analyzer		
Safety Analyzer Model No.		
Safety Analyzer Serial No.		
Calibration Date		
Date of Tests		
Name of Tester (Printed)		
Signature of Tester		
DV Product Unit	Max. Leakage Current	If < 5.0 mA, PASS
PEN Cabinet		
Global Operation Console		
Brainwave Cabinet (Optional Equipment)		
MNS Cabinet (Optional Equipment)		