

Run#	Part#	From		To	
M4502	6611001	System Cabinet	GND	Blower	GND
M0503	6611002	System Cabinet	F5-F6	Blower	J1
M0504	6611003	System Cabinet	F7-F8	Blower	J2
M3529	6611015	Gradient Coil-1,2	GFILTER-1/2	System Cabinet-1,2	X-/X+
M3530	6611016	Gradient Coil-3,4	GFILTER-3/4	System Cabinet-3,4	Y-/Y+
M4505	6613100	System Cabinet	J1000	Magnet side kit	J1000
	6613100-1	System Cabinet	J1000	Magnet side kit	J1001
M4503	6111117	Magnet	GND	Patient Table	GND
M3531	6611017	Gradient Coil-5,6	GFILTER-5/6	System Cabinet-5,6	Z-/Z+
M4500	6611050	MS1		RF-COM-GRND	
E3503	6611085	PP1	J12	PCIAA	J7
E3047	5308574	PCIAA	J6	SCIM	
				Cable to Cabinet	J26
E3500	6611096	FACILITY DISCONNECT		PP1	J18
M3527	6611097	EMERGENCY OFF SW		PP1	J18
E3506	6611010	System Cabinet	J15	RF-Door-switch	RF DOOR SWITCH
E0501	6611041	WC2		System Cabinet	T2
E0503	6611043	WC2		System Cabinet	T2
E0502	6611042	WC2		LCS1 (for Cabinet)	
E3502	6611036	Host-Eth1		System Cabinet	J63
E3504	6611044	WC2		System Cabinet	J301
E3505	6611505	WC1 and WC2		Cable #2044	J302
E3501	6611034	System Cabinet	J26	Run #3047	n.a
E3507	6611045	To Run #826		System Cabinet	J303
				To Run #827	Water Flow
E0507	6611020				
E0500	6611033	GOC	J1	System Cabinet	T4
M1500	6611054	System Cabinet	J44	MG3-A3-A1	J1
M1504	6611056	System Cabinet	J4	MG3-PEDIF	J1
M1501	6611055	MG3-A3	J1	RUN M1500	N/A
M1505	6611057	MG3-PEDIF	J1	RUN M1504	N/A
M4501	5333378	Magnet	GND	Rear Pedestal	Frame GND
M3535	5333388	System Cabinet	J11	Motor Drive	MG3-A7
M3534	5333386	System Cabinet	J74	Body QHS	J6
M3532	5333384	System Cabinet	J75	RF Body coil	J5
M3533	5333385	System Cabinet	J76	RF Body coil	J6
M3536	5333376	Cable #M3535	N/A	Motor Drive	N/A
M3526	5333379	Patient Comm Box	J7	MIC	N/A
M1516	5333381	LPCA Port A+	N/A	Dummy LOAD	
M1506	5333380	LPCA MIF	J18	PED IO	N/A
M3512	5333360	Mega Switch	J15	LPCA Port A&B LED	LEDS
M3510	5333361	Cable #M3537	N/A	MC Box	J21
M3537	5371203	Mega Switch	J18	Cable #M3510	N/A
M1508	5333362	Mega Switch	J21	MC Box	J9-J16
M3511	5333363	Mega Switch	J17	MAC Port Interface Board	P5
M1509	5442420	Mega Switch	J22	PA coil I/F	P4
M1503	5333382	Body QHS	J3	RF Body coil	J1
M1502	5333383	Body QHS	J4	RF Body coil	J2
M3528	5333387	Patient Comm Box	J6	PP1	J12
M1515	5333353	Mega Switch	J103	PP	J32
M1514	5333354	Mega Switch	J105	PP	J30

M1513	5333355	Mega Switch	J106	PP	J29
M3519	5333356	Mega Switch	J5	PP	J77
M3520	5333358	Mega Switch	J11	PP	J78
M1507	5333366	Mega Switch	J100	Body Preamp	
M1512	5333367	Mega Switch	J104	TNS Antenna	
M1510	5333368	Mega Switch	J1	RRX	J1
M1511	5333369	Mega Switch	J2	RRX	J2
M2502	5333370	Mega Switch	J107	RRX	J8
M3517	5333371	Mega Switch	J108	RRX	J7
M3515	5333372	Mega Switch	J102-1	Reg Box	J2-1
M3516	5333373	Mega Switch	J102-2	Reg Box	J2-2
M3518	6611518	Mega Switch	J20	SRI3	J20
M3513	5333374	Reg Box	J1-1	SCIF on PP	J120
M3514	5333375	Reg Box	J1-2	SCIF on PP	J121
P2500	6611500	RRX	J5	ICE	J32
	6611500-1	RRX	J5	ICE	J32
M3523	5338128	LED PWR Box	J3	SRI3	J1
M3524	6611524	LED PWR Box	J4	PAC3	J2
M3525	5338129	LED PWR Box	J5	PAL	PAL
414	46-317204G1	Front Bell MIC		Patient Comm Box	J9
421	46-287908G1	SRI3	J3	Table Control Box	J1
422	46-243775G668	SRI3	J2	Encoder Box-J4	J4
601	46-260723G1	Magnet P.S.		Magnet	
602	46-260723G2	Magnet P.S.		Magnet	
603	46-260726G1	S/C Shim P.S.		Terminal Box	
604	46-260724G1	Magnet P.S.		Terminal Box	
606	2128699	MRU		Terminal Box	
624	2218292	Cold Head		PP1	F1-F4
823	2263200-41	Magnet Monitor	J10	System Cabinet	J24
824TR	46-317359G975	PP1	J10	FJ1	
825TR	46-317359G977	PP1	J48	Magnet Monito	J8
826	46-328000G978	Magnet Monitor	J9	FJ3	
827	46-328578P1	FJ3		Water Flow-J1 / FJ4	MS5 A5 A1 J1 / FJ4
828	46-328000G974	PP1	J10	Terminal Box P403	
829	46-328000G976	PP1	J48	Magnet-FJ2	
830	2222797	MS1 FJ2		Magnet	J2
				Magnet	J1
				Magnet	TB1
831	2204488-5	FJ1		Magnet Monitor	J7
832	2271499	Magnet Pre Amplifier		Magnet SLEEVE/ COLDHEAD	
848	2221945	Operator Display	J5	PAC3	J4
862/863	2221946	Operator Display	J1, J2	Right/Left Operator Control	J1,J2
869	46-317252G4	Sperker		Patient Comm Box	J8
897	2221948	Operator Display	J3	SRI3	J4
1260	5123812	SRI3	J10	Bore Temp Sensor	
1261	5123813	Operator Display	J4	SRI3	J5
2009	5145249	System Cabinet	J47	Dock ASM	J1

2011/2012	5493045	System Cabinet	J92	SRI3-J11-J13/PAC3-J9,J11	CAM
2018	5145356	System Cabinet	J72	RF Body Coil	J3
2019	5145358	System Cabinet	J73	RF Body Coil	J4
457	46-317992G1	OM1		PP1	J17
458	46-243775G707	PP1	J17	OM3	
623	2155316-2	MS5		PP1	F1-F4
833	5508127-Rohs version (old version 46-328000G981)	FJ4		MS5 A1 A6 JR	
942	2271497	Mag Mon		Site Ethernet	
M3522	5441863-1	LED PWR Box	J1	RUN# M3544	
				PP	J122
M3544	5448314	RUN# 3544		RUN# 3545	
M3540	5393392	Mega Switch	J10	PP	J79
M1520	5413723	Mega Switch	J23	Port B	J23
M1521	5413724	Mega Switch	J24	Port B	J21
M3541	5333365	Mega Switch	J110	Port B	J22
1383	5143792-3	USB 1 wire		I Button	
1384	5143794-2	I Button		Cabinet Sensor	
E3037	5725077	ISC	J25	GOC	USB